Ecosite Bond

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English

Product description

Ecosite Bond is a light-curing single-component bonding agent. The material can be used as an adhesive on enamel and dentin with the self-etch, total-etch or selective-etch technique in combination with light-curing filling materials.

Intended purpose

Bonding agent for direct restorations of all cavity classes with light-curing composites, compomers, and ormocers and for adhesive luting of indirect restorations with light-curing veneer cements.

Limitations of use

Do not use the material if a sufficiently dry working area or the recommended application technique is not possible.

Indication

- Loss of tooth substance due to caries, trauma, tooth abrasion or development disorders: for direct restorations made of composites, compomers and ormocers or in adhesive luting of indirect restorations with light-curing veneer cements, e.g. with Vitique
- Tooth or periodontal injuries which require composite splints
- Esthetic requirements (e.g. Skyce tooth jewelry)

Contraindication

- Do not use directly on the exposed pulp.
- Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended user

Dentist

Notes for use

- Light-curing units should have an output of 450 nm and should be checked regularly.
- The light intensity should be a minimum of 400 mW/cm². Place the light-curing unit as close as possible to the material. Observe the manufacturer's instructions.
- In order to prevent premature polymerization due to ambient light, the material must be added to the mixing palette immediately prior to application.
- Ensure that the etched or bonded surface does not become contaminated prior to each subsequent working step.
- Do not remove the oxygen inhibition layer caused by the application because this is important for the bonding process with the material applied on top.
- Observe the manufacturer's instructions for all materials that are used with the product.
- Do not subject the dentin to excessive drying, particularly after etching.
- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

Recommended use

Prepare tooth structure and etch if necessary

1. Prepare the cavity in accordance with the general rules for adhesive technique. It is advisable to use a rubber dam. In deep cavities, protect areas close to pulp with a suitable liner.
2. Where necessary, etch as follows:
 - Selective-etch technique: Etch affected enamel with 37% phosphoric acid gel. The contact time of the etching gel on the enamel should be 20 to 60 seconds. The etched enamel should have a chalky white appearance. If this is not the case, etching must be repeated. Accidental etching of dentin has no effect on the bonding values, but the reaction time of the etching gel on dentin should not exceed 15 seconds.
 - Total-etch technique: Etch dentin and enamel with 37% phosphoric acid gel. Start the etching process with the enamel and then expand this to include the dentin. The reaction time on the enamel should be 20 to 60 seconds. The etched enamel should have a chalky white appearance. If this is not the case, etching must be repeated. The contact time on the dentin must not exceed 15 seconds.
3. Rinse with water for at least 15 seconds. Remove large areas of excess water.
4. After drying, an even, moist layer must remain on the etched dentin surface requiring treatment. If necessary, wet the dentin surface using a lightly moistened micro brush.

Apply adhesive to the tooth structure

1. Place 2–3 drops of Ecosite Bond on a mixing palette. Using a micro brush, rub the bond in for 10 seconds before leaving to work for 20 seconds.
2. Dry the surface with compressed air for 5 seconds to evaporate the solvent.

Note: Ensure that enough material remains on the tooth.

3. Light-cure the bond with a suitable LED polymerization device or a halogen polymerization device for 10 seconds.
4. Apply light-curing composite material to the cured Ecosite Bond and light-cure in accordance with manufacturer's specifications.

Interactions

Moisture, materials containing eugenol and oily air may inhibit polymerization at the contact area.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.
- Avoid contamination of the bottle content from unclean disposable brushes.

Product performance characteristics

SBS bonding values on etched and non-etched tooth substance with light-curing composite materials measured in accordance with ISO 29022:2013 are ≥ 15 MPa

Storage/disposal

- Store in a dry place at 2 to 25 °C (36 to 77 °F)!
- Disposal must comply with national regulations.

Composition

HEMA, Bis-GMA, ethanol, MDP, water, polymerized aminobenzoate derivative, additives.

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