DMG Etching Gel

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English

Product description

DMG Etching Gel is a universal etching gel (37 % phosphoric acid gel) with optimum flowability and thixotropy. DMG Etching Gel delivers reliable conditioning results and, due to its blue coloring, contrasts clearly with the tooth structure.

Intended purpose

Etching of enamel and dentine

Limitations of use

DMG Etching Gel is not suitable for etching ceramics.

Indication

Loss of tooth substance due to caries, fractures, tooth abrasion or developmental disorders.

Contraindication

Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons treated in the course of a dental procedure.

Intended user

Dentist

Notes for use

- Prior to the treatment continuing - e.g. through the application of a bonding agent - no contact should take place with the etched surface and this should not be contaminated with saliva. Should contamination occur after drying, the surface should be re-etched for approximately 10 seconds.
- If the material is applied into the patient's mouth with a single use applicator, the single use applicator must only be used on one patient for hygiene reasons.

Recommended use

1. Clean the tooth prior to preparation. Remove all cleaning residue with water spray. Prepare the cavity in accordance with the general rules for an adhesive technique. Bevel the enamel edges. It is advisable to use a rubber dam. In deep cavities in areas close to pulp provide protection in the form of a suitable liner, e.g. by using a small quantity of calcium hydroxide

and a thin layer of glass ionomer cement.

2. Apply DMG Etching Gel to the enamel edges of the cavity and allow to act for 20 – 60 seconds. When using the Total-Etch technique start the etching process with the enamel and then expand this to the dentin, so that the application time on the enamel is between 20 and 60 seconds, whilst the application time on the dentin does not exceed 15 seconds. The etched enamel should have a chalky white appearance. If this is not the case, etching must be repeated. Please observe the instructions for use applicable to the materials used with the DMG Etching Gel.
3. Rinse with water for a minimum of 15 seconds. Then dry using oil-free and water-free air.

Interactions

DMG Etching Gel must not be applied on acid-soluble linings.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

Etches dentin and enamel optimally because of its acid content and thixotropic behavior.

Storage/disposal

- Store in a dry place at room temperature 15-25 °C (59-77 °F).
- Do not use after the expiration date.
- Disposal must comply with national regulations.

Composition

Water, phosphoric acid, SiO₂, additives.

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