# EcuSphere

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English

### Product description

The EcuSphere System consists of four light-curing composites that have been specially developed for different purposes:

The radiopaque, high viscosity EcuSphere Carat is preferred for permanent fillings, especially in the posterior area.

The radiopaque EcuSphere Shape is best suited for modeling, as dentine material for high quality esthetic anterior tooth care and as a universal composite.

Radiopaque EcuSphere Flow is suitable for all indications that require an easy-flowing composite and it perfectly matches the other EcuSphere colors.

The color coding on the syringe shaft or Safetip plunger makes it easy to tell the difference between the materials (EcuSphere Carat= black, EcuSphere Shape = grey, EcuSphere Flow in the Luer-Lock Syringe).

### Intended purpose

EcuSphere Carat	Restorations of all cavity classes I – V, production of direct inlays and onlays
EcuSphere Shape	Restorations of all cavity classes I – V
EcuSphere Flow	Small restorations of cavity classes I, III, IV, V (cervical caries, root erosions, wedge-shaped defects), deciduous tooth cavities, minimally invasive fillings, base material, extended fissure sealing, Blocking and filling of cavity undercuts, Repairs of composite restorations

### Limitations of use

Applying EcuSphere is contraindicated if a dry working area or the recommended application technique is not possible.

### Indications

Loss of dental enamel due to caries, trauma, tooth abrasion or development disorders

### Contraindication

- Do not use directly on the exposed pulp.
- Do not use the material where there are known allergies to any of the components or contact allergies.

### Patient target group

Persons receiving treatment as part of a dental procedure.

### Intended user

Dentist

### Notes for use

- Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 600 mW/cm<sup>2</sup>. Position the light-curing unit as close to the material as possible.
- Do not remove the oxygen inhibited layer because this is required for the bond with the next layer.
- Observe the manufacturer's instructions for any other materials that are used with EcuSphere.
- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

### Recommended use

#### Preparation

1. Before preparation, clean the tooth and remove all cleaning residue with water spray. Determine the VITA tooth color.

**Note:** DMG recommends the use of a rubber dam.

2. Prepare the cavity in accordance with the general rules for adhesive technique. Bevel the enamel edges. Protect the areas close to the pulp with a suitable material.
3. Use matrices to protect the neighboring teeth.
4. Condition the dental enamel and apply the bonding agent according to the manufacturer's instructions for the material used. DMG advises using the bonding agent Ecosite Bond.

#### Application

**Note:** Apply all three EcuSphere materials in layers measuring max. 2 mm thick each and light-cure each layer for 20 seconds at a light intensity of  $\geq 1,000$  mW/cm<sup>2</sup>.

1. If necessary, apply a liner of thin-flowing composite or compomer (e.g. EcuSphere Flow or Ionosit Baseline).
2. Fill the cavity with EcuSphere Carat or EcuSphere Shape in the desired shade with layers measuring max. 2 mm thick. Light-cure each layer for 20 seconds at a light intensity of  $\geq 1,000$  mW/cm<sup>2</sup> or model the dentin core made of opaque material (such as EcuSphere Shape or EcuSphere Carat) in the desired shade.

#### Finishing

- Remove excesses, shape and finish, e.g. with diamond finishing burs and flexible discs with decreasing grit sizes.
- Carry out high-gloss polishing of EcuSphere Carat and EcuSphere Shape with polishing pastes using polishing cups/discs.

#### Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid this.
- The use of mouthwashes, plaque indicators and chlorhexidine may lead to discoloration.

#### Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

#### Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Reuse: Standard surface disinfection as required. No specific maintenance required.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

#### Product performance characteristics

EcuSphere Carat complies with DIN EN ISO 4049:2019 (Type 1, Class 2, Group 1). Polymerization depth  $\geq 2$  mm with a 20 second light-curing time when the lamp's light intensity is  $\geq 1,000$  mW/cm<sup>2</sup>. Polymerization depth  $\geq 2$  mm with a 40 second light-curing time when the lamp's light intensity is  $\geq 600$  mW/cm<sup>2</sup>.

EcuSphere Shape complies with DIN EN ISO 4049:2019 (Type 2, Class 2, Group 1). Polymerization depth  $\geq 2$  mm with a 20 second light-curing time when the light intensity of the lamp is  $\geq 1000$  mW/cm<sup>2</sup>. Polymerization depth  $\geq 2$  mm with a 40 second light-curing time when the light intensity of the lamp is  $\geq 600$  mW/cm<sup>2</sup>.

EcuSphere Flow complies with DIN EN ISO 4049:2019 (Type 2, Class 2, Group 1). Polymerization depth  $\geq 2$  mm with a 20 second light-curing time when the light intensity of the lamp is  $\geq 1000$  mW/cm<sup>2</sup>. Polymerization depth  $\geq 2$  mm with a 40 second light-curing time when the light intensity of the lamp is  $\geq 600$  mW/cm<sup>2</sup>.

#### Storage/disposal

- Store in a dry place at room temperature 15-25 °C (59-77 °F).
- Disposal must comply with national regulations.

#### Composition

##### EcuSphere Carat

Dental glass, bis-GMA, EBPADMA, UDMA, TEDMA, SiO<sub>2</sub>, additives. Inorganic filling materials: approx. 57 vol%, 0.02 – 0.7 µm.

## EcuSphere Shape

Dental glass, bis-GMA, TEDMA, UDMA, EBPADMA, SiO<sub>2</sub>, additives. Inorganic filling materials: approx. 57 vol%, 0.02 – 1.5 µm.

## EcuSphere Flow

Dental glass, bis-GMA, UDMA, HDDMA, DDDMA, SiO<sub>2</sub>, additives. Contains TPO. Inorganic filling materials: approx. 41 vol%, 0.02 – 3 µm.

Contains: 2-ethylhexyl 4-(dimethylamino)benzoate

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