# LuxaBond Total Etch

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English

### Product description

LuxaBond Total Etch is a dual-curing bonding agent for enamel and dentin, suitable above all for use with dual and self-curing composites, such as LuxaCore Z Dual. It is necessary to etch enamel and dentin surfaces (total etch technique).

Due to the dual-curing nature of the material, LuxaBond Total Etch is particularly useful in root canals, where light-curing is not possible.

### Intended purpose

Dual-curing bonding agent for use in conjunction with the total etch technique and a composite for core build-ups, for root post cementation and for the cementation of crowns/bridges, inlays and onlays.

### Limitations of use

Do not use the material if a sufficiently dry working area or the recommended application technique is not possible.

### Indication

Loss of tooth substance/tooth due to caries, trauma, tooth abrasion or development disorders

### Contraindication

- Do not use directly on the exposed pulp.
- Do not use the material if pulpitis is present.
- Do not apply the material in cases of allergies to any of the components, in particular benzoyl peroxide, or in the event of contact allergies.

### Patient target group

People treated in the course of a dental procedure.

### Intended user

Dentist

### Notes for use

- Light-curing units should have an output of 450 nm and should be checked regularly.
- The light intensity should be at least 600 mW/cm<sup>2</sup>. Place the light-curing unit as close as possible to the material.
- Always use etching gel on enamel and dentin surfaces.
- In order to prevent premature polymerization due to ambient light, the material should be placed onto the mixing pad immediately prior to the application.
- Ensure that the etched or bonded surface does not become contaminated prior to each subsequent working step.
- Do not remove the oxygen inhibition layer caused by the application because this is important for the bonding process with the material applied on top.
- Observe the manufacturer's instructions for all materials that are used with LuxaBond Total Etch.
- If the material is applied into the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

### Recommended use

#### For use in root canals, for cementing root posts, for cementing inlays

- Prepare the root canal or cavity in accordance with the general rules for adhesive procedures. It is advisable to use a rubber dam.
- Etch dentin and enamel - insofar as is present - with 37% phosphoric acid gel. Start the etching process with the enamel and then expand this to include the dentin. The reaction time on the enamel should be 20 to 60 seconds. The reaction time on the dentin must not exceed 15 seconds.

**Note:** When using LuxaBond Total Etch it is essential to etch the affected enamel and dentin surfaces.

- Rinse with water for a minimum of 15 seconds.
- Remove excess water. After the rinsing process an even moist layer must remain on the dentin surface to be treated. If necessary, wet the dentin surface using a lightly moistened micro brush.
- Use a brush to work 1 to 2 drops of Pre-bond into the etched tooth structure for 15 seconds. Remove any excess, e.g. with a paper tip. Blow gently using oil-free and water-free air. The surface should appear evenly moist.
- Mix 1 to 2 drops of Bond A and Bond B in a trough in the mixing palette at a ratio of 1:1 for approx. 5 seconds and work into the tooth structure for 20 seconds using a micro-brush.
- Carefully remove bond excess from the preparation, e.g. with a paper tip, and then gently blow with oil-free and water-free air until the surface appears evenly moist.
- Apply the luting cement in accordance with the manufacturer's instructions.

**Note:** The working time of the luting cement may be shortened due to the contact with LuxaBond Total Etch.

**Note:** The use of LuxaCore Z Dual is recommended for the cementation of root posts.

#### Core build-ups, cementation of crowns and bridges

**Note:** Since excess material can form at the border between tooth and matrix (pooling effect), use a matrix only after application if possible.

- Prepare the cavity in accordance with the general rules for the adhesive technique. Bevel the enamel edges. It is advisable to use a rubber dam. In deep cavities in areas close to pulp provide protection in the form of a suitable liner, e.g. by using a small quantity of calcium hydroxide and a thin layer of glass ionomer cement.
- Etch dentin and enamel with 37% phosphoric acid gel. Start the etching process with the enamel and then expand this to include the dentin. The contact time on the enamel should be 20 to 60 seconds. The contact time on the dentin must not exceed 15 seconds.

**Note:** When using LuxaBond Total Etch it is essential to etch the affected enamel and dentin surfaces.

- Rinse with water for at least 15 seconds.
- Remove excess water. After the rinsing process an even moist layer must remain on the dentin surface to be treated. If necessary, wet the dentin surface using a lightly moistened micro brush.
- Use a brush to work 1 to 2 drops of Pre-bond into the etched tooth structure for 15 seconds. Remove any excess, e.g. with a paper tip. Blow gently using oil-free and water-free air. The surface should appear evenly moist.
- Mix 1 to 2 drops of Bond A and Bond B in a trough in the mixing palette at a ratio of 1:1 for approx. 5 seconds and work into the tooth structure for 20 seconds using a micro-brush.
- Carefully remove any bond excess and then gently blow with oil-free and water-free air until the surface appears evenly moist.
- If required, light-cure LuxaBond Total Etch for a minimum of 10 seconds using a suitable light-curing unit.
- Apply core build-up material (e.g. LuxaCore Z Dual) or luting cement in accordance with the

manufacturer's instructions.

**Note:** The working time of the luting cement or the core build-up material may be shortened due to the contact with LuxaBond Total Etch.

## Timing

Working time of the mixture of Bond A + Bond B when protected from light	5 min
Contact time of the etching gel on enamel	20-60 s
Contact time of the etching gel on dentin	max. 15 s
Pre-bond application time	15 s
Mixing time of Bond A + Bond B	5 s
Application time (mixture of bond A and bond B)	20 s
Light-curing time (optional)	min. 10 s

## Interactions

- Products containing eugenol, for example temporary cements or temporary root canal filling materials, may inhibit the polymerization.
- Products containing eugenol and other substances that inhibit polymerization must not be used together with the material.

## Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

## Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

## Product performance characteristics

SBS bonding properties on etched tooth substance measured according to ISO 29022:2013 are  $\geq 10$  MPa

## Storage/disposal

- Store LuxaBond Total Etch in a refrigerator at 2 - 8 °C.
- Close bottles immediately after use.
- Do not use the material past the stipulated shelf life or if the fluid assumes a gelatinous consistency.
- Ensure that the lids of bond bottles are not switched in order to avoid any cross contamination.
- Dispose of in accordance with national regulations.

## Composition

Pre-bond: Ethanol, sodium benzenesulfinate.

Bond A: EBPADMA, HEMA, polymeric aminobenzoate derivative, CQ, additives.

Bond B: HEMA, bis-GMA, GDMA, glycidyl methacrylate-modified polyacrylic acid, BPO, additives.

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