LuxaBond Universal

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Bond A: HEMA, Bis-GMA, MDP, BPO, Additive.

Bond B: Ethanol, Wasser, Na-toluolsulfonat, DABE, CQ, Additive.

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English

Product description

LuxaBond Universal is a dual-curing 2-bottle bonding agent. The material can be used as an adhesive on enamel and dentine with the self-etch, total-etch or selective-etch technique in combination with dual-curing, self-curing or light-curing composites (e.g. PermaCem Universal, LuxaCore Z Dual, Ecosite Elements). Thanks to its dual-curing property, the material can be used in the root canal area where it is difficult for light to penetrate.

The material can also be used as a primer for indirect restorations made of metal (base metal, precious metal), oxide and silicate ceramics or composites.

Intended purpose

As an adhesive for:

- Light-curing fissure-sealing materials
- Direct restorations made of light-curing composites and compomers
- Light-curing, self-curing and dual-curing core build-up materials
- Cementing of indirect restorations (inlays, onlays, crowns and bridges) and root posts with light-curing, self-curing and dual-curing composite-based luting cements
- Cementing of veneers with light-curing luting composite cements in conjunction with the enamel etching technique (e.g. with Vitique Cement)

As a primer for indirect restorations made from:

- Metals (base metal, precious metal)
- Zirconium oxide and aluminum oxide
- Silicate ceramic
- Composite

Limitations of use

Do not use the material if a sufficiently dry working area or the recommended application technique is not possible.

Indication

- Loss of tooth substance/tooth due to caries, trauma, tooth abrasion or development disorders

Contraindications

- Do not use directly on the exposed pulp.
- Do not use the material if pulpitis is present.
- Do not use the material in cases of known allergies to any of the components, in particular benzoyl peroxide, or in the event of contact allergies.

Patient target group

People treated in the course of a dental procedure.

Intended user

Dentist

Notes for use

- Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be at least 600 mW/cm². Place the light-curing unit as close as possible to the material. Observe the manufacturer's instructions.
- Bring the material to room temperature (15 – 25 °C/59 – 77 °F) before use.
- In order to prevent premature polymerization due to ambient light, the material should be placed onto the mixing pad immediately prior to the application.
- Ensure that the etched or bonded surface does not become contaminated prior to each subsequent working step.
- Do not remove the oxygen inhibition layer caused by the application because this is important for the bonding process with the material applied on top.
- In deep cavities, material excesses can form at the margin between tooth and matrix (»pooling effect«). Therefore, place a matrix band, if possible, only after the bond is applied.
- It is advisable to use a rubber dam.
- If the material is applied into the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

Recommended use

Prepare tooth structure and etch if necessary

- Carry out cavity or endodontic treatment and prepare the root canal for the root post insertion in accordance with the general rules of adhesive technology.
- In deep cavities, protect areas close to pulp with a suitable liner.

Note: Clinically sufficient bonding values are achieved (self-etch) by applying LuxaBond Universal. Additional etching is no longer needed. The bonding values can be increased, if necessary, by etching the enamel and dentin.

3. Where necessary, etch as follows:

- Selective-etch technique: Etch affected enamel with 37% phosphoric acid gel. The contact time of the etching gel on the enamel should be 20 to 60 seconds. The etched enamel should have a chalky white appearance. If this is not the case, etching must be repeated. Accidental etching of dentin has no effect on the bonding values, but the contact time of the etching gel on dentin should not exceed 15 seconds.

- Total-etch technique: Etch dentin and enamel with 37% phosphoric acid gel. Start the etching process with the enamel and then expand this to include the dentin. The contact time of the etching gel on the enamel should be 20 to 60 seconds. The etched enamel should have a chalky white appearance. If this is not the case, etching must be repeated. The contact time on the dentin must not exceed 15 seconds.
- 4. Rinse with water for at least 15 seconds.
- 5. After the rinsing process an even, moist layer must remain on the dentin surface to be treated. If necessary, wet the dentin surface using a lightly moistened micro brush.

Apply LuxaBond Universal as an adhesive to the tooth structure

1. Mix 1 to 2 drops each of bond A and bond B in a trough in the mixing pad in a 1:1 ratio for approx. 5 seconds.
2. Apply a layer of the bond mix onto the preparation using a micro brush and work into the tooth structure for 20 seconds.
3. Remove excess material carefully.
4. Direct a stream of oil-free and water-free air at the material for at least 10 seconds until the surface has an evenly moist and shiny appearance and is firm.
5. If the mixed material in the mixing palette is still liquid, use the liquid material. If the mixed material in the mixing pad has already hardened, remix 1 drop each of bond A and bond B in a trough in the mixing pad in a 1:1 ratio for approx. 5 seconds.
6. Apply another layer of the bond mix onto the preparation using a micro brush and work into the tooth structure for 20 seconds.
7. Remove excess material carefully.
8. Direct a stream of oil-free and water-free air at the material for at least 10 seconds until the surface has an evenly moist and shiny appearance and is firm.
9. Cure the material as follows:
 - For core build-up materials, light-curing filling materials: Light-cure LuxaBond Universal for at least 20 seconds with a suitable light-curing unit.
 - For use when cementing root posts, inlays, onlays, crowns, bridges and veneers: Allow LuxaBond Universal to set. After the working time is completed, the material is sufficiently cured. In addition, the material can be light-cured with a suitable light-curing unit for at least 20 seconds.

Note: LuxaCore Z Dual is recommended for core build-up and for root post cementation. PermaCem Universal is recommended for cementing indirect restorations.

Apply LuxaBond Universal as a primer for indirect restorations made of metals (non-noble metal, precious metal), oxide and silicate ceramics or composites

1. Prepare the clean and dry restoration in accordance with manufacturer's instructions. If not otherwise described, sandblast the surfaces to be bonded with aluminum oxide ($\leq 50 \mu\text{m}$).

Note: Phosphoric acid can hinder an optimal bond with zircon oxide aluminum oxide, and metal restorations.

2. Clean the restoration with alcohol and dry with air that is free of oil and water. Etch the silicate ceramic in advance with 5% hydrofluoric acid in accordance with the manufacturer's instructions and rinse thoroughly with water spray. Dry the restoration with oil-free and water-free air.
3. Mix 1 to 2 drops each of bond A and bond B in a trough in the mixing pad in a 1:1 ratio for approx. 5 seconds.
4. Apply the bond mix onto the restoration using a micro-brush and work in for 20 seconds.
5. Remove excess bond carefully from the restoration.
6. Direct a stream of oil-free and water-free air at the material for at least 10 seconds until the surface has an evenly moist and shiny appearance and is firm.
7. Additional light-curing is optional.

Note: The bonding values can be increased by the optional light-curing.

Timing

Optional: Contact time of the etching gel on enamel	20 – 60 s
Optional: Contact time of the etching gel on dentin	max. 15 s
Working time when keeping protected from light and incl. mixing time depending on the ambient temperature	max. 60 s
Mixing time	5 s
Work-in time as adhesive on tooth: 2 layers	Massage in for 2x 20 s + Blow with compressed air for 10 s
Working time as a primer on restorations: 1 layer	Massage in for 1x 20 s + Blow with compressed air for 10 s
Optional: Light-curing time	min. 20 s

Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.
- The working time of the luting cement or the core build-up material may be reduced if there is contact with LuxaBond Universal.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

SBS bonding values on etched and non-etched tooth substance measured in accordance with ISO29022:2013 are $\geq 7 \text{ MPa}$

Storage/disposal

- Store in a dry place in the refrigerator at 2–8 °C/36–46 °F.
- Close bottles immediately after use.
- Do not use the material past the stipulated shelf life or if the fluid assumes a gelatinous consistency.
- Ensure that the lids of bond bottles are not switched in order to avoid any cross contamination.

Composition

Bond A: HEMA, bis-GMA, MMHE, BPO, additives.

Bond B: Ethanol, water, Na-toluenesulfinate, DABE, CO₂, additives.

Contains: Ethyl 4-(dimethylamino)benzoate

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