Optimized  
Radiopacity

# LuxaPost

Glasfaserverstärkter Composite Wurzelstift · Glass fibre composite post  
Pivot radiculaire en matériau composite renforcé par des fibres de verre  
Perno endocanalare in fibra di vetro · Poste de resina composite con fibra de vidrio  
Pino de raiz em compósito reforçado a fibra de vidro · Glasvezelversterkte composiet-  
stompobouw · Glasfiberforstærket komposit-rodstift · Glasfiberförestärkt komposit-  
rotstift · Kompozytowy sztyft korzeniowy wzmocniony włóknem szklanym  
Усиленный композитом стекловолоконный штифт

## **Inhalt | Contents**

<b>DE</b> – Gebrauchsinformation	4
<b>EN</b> – Instructions for use	13
<b>FR</b> – Mode d'emploi	21
<b>IT</b> – Istruzioni d'uso	30
<b>ES</b> – Instrucciones de empleo	39
<b>PT</b> – Instruções de uso	48
<b>NL</b> – Gebruiksaanwijzing	57
<b>DA</b> – Brugsanvisning	66
<b>SV</b> – Användarinstruktioner	74
<b>PL</b> – Instrukcja użycia	83
<b>RU</b> – Информация по применению	92

## **Produktbeschreibung**

---

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## **Product Description**

---

Fiberglass-reinforced, conical composite post, x-ray opaque and pre-silanized.

## **Indications**

---

- Following endodontic treatment, in order to establish reliable retention of the core build-up where insufficient coronal tooth structure exists.
- Only for use with dual or self-curing luting cements.

## **Contraindications**

---

The material must not be used in conjunction with light-curing luting cements.

## **Basic safety instructions**

---

- ▶ For dental use only!
- ▶ Keep out of reach of children!

## **Interactions**

---

- ▶ Cements containing eugenol can affect the curing of luting cements. Avoid the use of materials containing eugenol.
- ▶ Do not use chlorhexidine preparations when rinsing the root canal.
- ▶ With phosphate and glass ionomer cements exhibiting lower compression strengths, it is possible that extreme loads together with conical posts may result in the post pressing into the cement and potentially causing a springing of the tooth (wedging or shim effect). Avoid the use of phosphate and glass ionomer cements.
- ▶ Moisture or chemical vapors (e.g. from acids) can cause corrosion. Protect the material against moisture and chemical vapors.

## **Instructions for use**

---

- ▶ In order to improve the fracture resistance of core build-up and restoration,

ensure that a minimum 2 mm of tooth structure remains above the gum after preparation, in order that this can be covered by the crown itself.

- ▶ In order to ensure the enhanced transfer of axial forces, ensure that the minimum wall thickness of the dentin is 1 mm.
- ▶ When preparing the root canal always use a LuxaPost drill of the appropriate size. This guarantees an adequate cement gap.
- ▶ Use LuxaPost drills no more than 10 times.
- ▶ Sterilize LuxaPost drills after every use.
- ▶ Light-curing units should have an output of 450 mW and should be checked regularly. The light intensity should be a minimum of 400 mW/cm<sup>2</sup>. Place the light as close as possible to the material.
- ▶ Observe the respective manufacturer's instructions when using any other materials.

In order to generate an homogenous monoblock, for the cementing and core build-up

DMG recommends using the dual-curing LuxaCore Z, or LuxaCore in combination with the dual-curing adhesive system LuxaBond-Total Etch.

DMG recommends the use of a rubber dam.

### **Prior to every use**

---

- ▶ LuxaPost drills are not sterile upon delivery. Disinfect LuxaPost drills in a disinfecting bath. Use a disinfection bath with a screen insert. Only use special combination disinfectants and cleaning solutions, which provide simultaneous protection against corrosion.
- ▶ LuxaPost drills must not be left in the disinfecting solution for extended periods (e.g. overnight).
- ▶ Do not touch contaminated LuxaPost drills with bare hands. Remove LuxaPost drills from the disinfecting solution using tweezers and dry them with highly absorbent paper (paper towels).

- ▶ The root posts are not sterile upon delivery. Prior to insertion, perform cold sterilization of the root posts using glutaraldehyde or phenol. Subsequently rinse and dry the root posts thoroughly.

## Recommended use

---

1. Using periapical x-ray images and a LuxaPost drill of the appropriate size, determine the diameter of the post to be used.

**Please note:** The diameter of the post should be at least as large as the root canal to be filled.

2. Using the LuxaPost drill, determine the depth of the preparation.

**Please note:** 4mm of the root filling should remain at the apical end after preparation and the root post length should correspond to at least the length of the clinical crown.

3. Slide a depth marker onto the LuxaPost drill in accordance with the desired prep-

aration depth, in order to avoid excessive preparation.

**Please note:** In order to simplify specification of the post length, it is possible to use the measuring card enclosed (scale of 1:1 for analogue x-ray images).

**Caution!** Incorrect use of the LuxaPost drill can lead to risks to the patient, the user and to wearing or breakage of the instruments.

- ▶ Slide the LuxaPost drill fully into the drill chuck.
- ▶ Ensure sufficient water spray cooling (> 50 ml/min).
- ▶ Observe the speeds and contact pressures.
- ▶ Avoid tilting, bending or levering.

4. Remove gutta-percha and prepare the root canal. In order to do so, use a LuxaPost drill with an appropriate size for the root post (see color coding).

**Please note:** Prior to preparation of the root canal, gutta-percha can be removed to the

desired depth using a Gates Glidden drill, a Peeso Reamer and/or a heated instrument.

5. Check the seating of the depth marker through comparison with the x-ray image.
6. Rinse the root canal and dry using an air jet and paper tips.
7. Tentatively insert the root post into the root canal.

**Please note:** In order to ensure sufficient retention, the post should be coronally roughly the length of the crown and should reach roughly twice the coronal portion into the tooth root.

**Caution!** When shortening the post, sanding dust can get into the eyes and respiratory tract.

- ▶ Wear mouth protection and goggles.
- ▶ Do not inhale the sanded dust and use an extractor.

8. Shorten the post to the desired length extraorally, using a water-cooled, fast-action diamond grinder or a diamond cutting wheel. The post must be shortened from the occlusal end.
9. Clean the post using alcohol and dry with oil-free air.

**Please note:** The LuxaPost is pre-silanized. Additional silanization is unnecessary. If additional silanization is however required, it is possible to carry out the silanization process using DMG Silane or another suitable silane. Observe the manufacturer's instructions for the material used.

10. If the bonding agent used requires an additional etching step then etch inside the canal in accordance with the etching gel manufacturer's instructions.
11. Apply bonding agent in accordance with the manufacturer's instructions. If an adhesive system with a light-curing component is used, remove the excess from the canal prior to the application of light,

e.g. using a paper tip. Otherwise the lumen may be blocked. If a core build-up is planned using a composite material, it is also necessary to bond the coronal portion of the preparation.

12. Coat the walls of the root canal with a resin-based attaching cement, e.g. using a lentulo spiral.

Please note: LuxaCore Z-Dual and LuxaCore-Dual can be applied to the root canal using special Endo tips. If desired, it is additionally possible to coat the post with material.

13. Carefully slide the post into the root canal. Allow the excess cement to flow out and remove.

Please note: When using LuxaCore Z-Dual or LuxaCore-Dual for the attachment and core build-up, it is not necessary to remove the excess.

14. Start the bonding of a dual-hardening material by light-curing the upper end of

the canal. Observe the manufacturer's instructions for the cement used.

15. In order to carry out a subsequent composite core build-up, bond the post and the remaining dentin.

Please note: In order to create the shape it is possible to use a matrix band or post form. The use of LuxaCore Z-Dual or LuxaCore-Dual results in an homogenous monoblock.

**Caution!** When processing the core build-up, sanding dust can get into the eyes and respiratory tract.

- ▶ Wear mouth protection and goggles.
- ▶ Do not inhale the sanded dust and use an extractor.

16. After curing is complete, prepare the core build-up for receiving the permanent restoration using a diamond or hard metal grinder.

## Removing the root post

If renewed treatment of the tooth is required, it is possible to remove the set post as follows:

**Caution!** Incorrect use of the LuxaPost drill can lead to risks to the patient, the user and to wearing or breakage of the instruments.

- ▶ Slide the LuxaPost drill fully into the drill chuck.
  - ▶ Ensure sufficient water spray cooling (> 50 ml/min).
  - ▶ Observe the speeds and contact pressures.
  - ▶ Avoid tilting, bending or levering.
1. Drill out the inserted post using a new LuxaPost drill of the appropriate size. When removing root posts with larger diameters, initially start with the smallest LuxaPost drill and move up through the larger LuxaPost drills.
  2. Perform follow-up X-rays as needed.

## Disinfection

---

- Strictly observe the user information, reaction time and suitability of disinfectants and cleaning substances for certain types of instruments, as detailed in the respective manufacturer's instructions. DMG refers in this regard to the list of disinfection processes tested per the Guidelines for the testing of chemical disinfectants and found to be effective by the German association for hygiene and microbiology (DGHM), and to the list published by the Robert Koch Institute (RKI).
- Information regarding suitable and validated disinfectants and disinfection processes (including all parameters) can also be found on the internet at [www.rki.de](http://www.rki.de) and [www.dghm.de](http://www.dghm.de).
- One of the recommended validated processes for instrument disinfection is the use of lysetol (4% dilution, reaction time 1 h) from Schülke & Mayr GmbH.

- The plastic packaging can be disinfected using a standard surface disinfectant.

## Ultrasonic disinfection

---

The use of ultrasound reduces the disinfection time to 15 minutes. The LuxaPost drill bits must not come into contact with hard objects during this process. Cleaning in a thermal disinfectant is not possible.

## Checking/secondary cleaning

---

- After disinfection, check LuxaPost drills for contamination. If necessary, carry out secondary cleaning using an ultrasonic cleaner. If necessary, carry out secondary cleaning manually using a soft brush.
- Check the functionality. LuxaPost drills that have oxidized, become blunt or exhibit irregular running must be separated and no longer used.

## Sterilization

---

Autoclave (if possible use pre and post-vacuum types) at 121 °C or 134 °C. Avoid steri-

lization temperatures above 180 °C. Exceeding this can lead to the drill losing hardness.

## Storage and shelf life

---

- ▶ Store in a dust-free and dry place at 2 to 25 °C.
- ▶ Protect against chemical vapors.
- ▶ Do not use after the expiration date.
- ▶ Observe the specifications of the hygiene plan and the respective manufacturer's instructions for devices and materials used.

## Ingredients

---

Glass fibers in a Bis-GMA based resin matrix (LuxaPost).

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DMG Chemisch-Pharmazeutische Fabrik GmbH