LuxaPost

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English

Product description

Glass fiber-reinforced, conical composite root post, radiopaque and pre-silanized.

Intended purpose

Retention of core build-ups in endodontically treated teeth

Limitations of use

- Do not use phosphate or glass ionomer cements for luting the post, because these cements have lower compressive strengths.
- Only for use with dual or self-curing luting cements.
- Do not use with light-curing luting cements.

Indication

Insufficient coronal tooth structure for the reconstruction of the tooth crown

Contraindication

- Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

People to be treated in the course of a dental procedure.

Intended user

Dentist

Notes for use

- In order to improve the fracture resistance of the core build-up and restoration, ensure that a minimum of 2 mm of tooth structure remains above the gum after preparation, so that it can be encased by the definitive crown.
- For enhanced transfer of axial forces, ensure that the minimum wall thickness of the dentin is 1 mm.
- When preparing the root canal, always use a LuxaPost drill of the corresponding size. This ensures that a space of adequate size is created for the cement.
- Use LuxaPost drills no more than 10 times.
- Sterilize LuxaPost drills after every use.
- Light-curing units should emit at 450 nm and should be checked regularly. The light intensity should be a minimum of 600 mW/cm². Place the light guide as close as possible to the material.
- Follow the respective manufacturer's instructions when using any additional materials.
- In order to create a homogeneous system, DMG recommends using the dual-curing LuxaCore Z in combination with the dual-curing adhesive system LuxaBond Total Etch or LuxaBond Universal for post cementation and core build-up.
- DMG recommends the use of a rubber dam.

Before first use

LuxaPost drills

LuxaPost drills are supplied non-sterile. Disinfect LuxaPost drill in a tray bath, then autoclave. Use a disinfection bath with a screen insert. Only use special combined cleaning and disinfection solutions that provide simultaneous protection against corrosion. LuxaPost drills must not be left in the disinfectant solution for prolonged periods (e.g. overnight).

Use tweezers to remove disinfected LuxaPost drills from the disinfection solution and dry them with highly absorbent paper (paper towels). Do not touch LuxaPost drills with bare hands.

LuxaPost root posts

Root posts are supplied non-sterile. Before using the root post, disinfect it with alcohol.

Recommended use

1. Use periapical X-ray images and a LuxaPost drill of the appropriate size to determine the diameter of the root post to be used.

Note: The diameter of the root post should be at least as large as the root canal to be filled.

2. Using LuxaPost drill, determine the depth of the preparation.

Note: After preparation, 4 mm of the root canal filling must be retained at the apical end and the post length must be at least equal to the length of the clinical crown.

3. Slide a depth marker onto the LuxaPost drill to the required depth to avoid over-preparation.

Note: In order to simplify the determination of the root post length, it is possible to use the enclosed measuring card (scale of 1:1 for analog X-ray images).

Caution! Improper use of the LuxaPost drill can lead to risks to the patient and the user, and to excessive wear or breakage of the instrument.

- ▶ Insert LuxaPost drill fully into the handpiece chuck.
- ▶ Ensure adequate water spray cooling (> 50 ml/min).
- ▶ Observe the limits for rotational speed and contact pressure. Avoid tilting, bending or levering.

4. Remove gutta-percha and prepare the root canal. In order to do so, use a LuxaPost drill of the appropriate size for the root post (see color coding).

Note: Prior to preparation of the root canal, the gutta-percha can be removed from the desired depth using a Gates Glidden drill, a Peeso reamer, and/or a heated instrument.

5. Check the seating of the depth marker by comparing it with the X-ray image.
6. Rinse the root canal and dry using an air jet and paper points.
7. Insert the root post into the root canal for verification.

Note: To ensure sufficient retention, the coronal part of the root post should be roughly the length of the crown and the root part of the post should be roughly twice the length of the coronal portion.

Caution! When shortening the root post, grinding dust can enter the eyes and respiratory tract.

- ▶ Wear a dust mask and protective goggles.
- ▶ Use suction for grinding dust and do not inhale it.

8. Shorten the root post to the desired length extra-orally using a water-cooled, fast-action diamond grinder or a diamond cutting wheel. The root post must be shortened from the occlusal end.

Note: Do not use scissors or pliers to shorten the root post.

9. Clean the root post using alcohol and dry with oil-free air.

Note: LuxaPost is pre-silanized. Additional silanization is unnecessary. If additional silanization is desired, it is possible to carry out the silanization process using a suitable silane. Observe the manufacturer's instructions for the material used.

10. If the bonding agent to be used requires an additional etching step, then etch inside the canal in accordance with the etching gel manufacturer's instructions.

11. Apply the bonding agent in accordance with the manufacturer's instructions. If an adhesive system with a light-curing component is used, remove the excess from the canal prior to the exposure to light, e.g. using a paper point. Otherwise the lumen may be blocked. If a composite core build-up is planned, the coronal portion of the preparation must also be bonded.

12. Coat the walls of the root canal with a resin-based luting cement.

Note: LuxaCore Z Dual can be applied to the root canal using special Endo Tips. If desired, it is additionally possible to coat the root post with material.

13. Glass fiber posts must be fully embedded in cement or composite. Exposed fibers absorb water, which can alter the mechanical properties of root posts.

14. Carefully insert the root post into the root canal. Allow excess cement to flow out and then remove the excess.

Note: When using LuxaCore Z Dual for the luting and core build-up, it is not necessary to remove the excess.

15. Start the bonding of a dual-curing material by light-curing the upper end of the canal. Observe the manufacturer's instructions for the cement used.

16. For a subsequent composite core build-up, bond the root post and the remaining dentin.

Note: In order to shape the core, it is possible to use a matrix band or core form. The use of LuxaCore Z Dual results in a homogeneous system.

Caution! When finishing the core build-up, grinding dust can enter the eyes and respiratory tract. Wear a dust mask and protective goggles. Use suction for grinding dust and do not inhale it.

17. After curing is complete, prepare the core build-up in preparation for the permanent restoration using a diamond or carbide bur.

Removing the root post

If retreatment of the tooth is required, it is possible to remove the cemented root post as follows:

Caution! Improper use of the LuxaPost drill can lead to risks to the patient and the user, and to excessive wear or breakage of the instruments.

- ▶ Insert LuxaPost drill fully into the handpiece chuck.
- ▶ Ensure adequate water spray cooling (> 50 ml/min).
- ▶ Observe the limits for rotational speed and contact pressure.
- ▶ Avoid tilting, bending or levering.

1. Drill out the inserted root post using a new LuxaPost drill of the corresponding size. When removing root posts with larger diameters, start with the smallest LuxaPost drill and continue with the larger LuxaPost drills.
2. If necessary, make a control radiograph.

Disinfecting LuxaPost drill

- Strictly follow the user information, exposure time and suitability of disinfectants and cleaning agents for certain types of instruments, as specified in the respective manufacturer's instructions. In this regard, DMG refers to the list of disinfection processes tested and found to be effective according to the guidelines for the testing of chemical disinfectant solutions by the German Association for Hygiene and Microbiology (DGHM), and to the list published by the Robert Koch Institute (RKI).
- Information regarding suitable and validated disinfectant solutions and disinfection processes (including all parameters) can also be found on the Internet at www.rki.de and www.dghm.de.
- The plastic packaging can be disinfected using a standard surface disinfectant.

Disinfecting LuxaPost drill in ultrasound

Using ultrasound reduces the disinfection time to 15 minutes. The cutting edges of the LuxaPost drills must not come into contact with any hard objects during this process. Cleaning in a thermal disinfectant is not possible.

Control/subsequent cleaning

- After disinfection, check the LuxaPost drills for contamination. If necessary, carry out a subsequent cleaning using an ultrasonic cleaner. If necessary, manually clean the drill using a soft brush.
- Check the functionality. LuxaPost drills that have become oxidized, blunt or exhibit irregular rotation must be discarded and no longer used.

Sterilizing the LuxaPost drill

Autoclave (ideally use pre- and post-evacuating types) at 121 °C (249 °F) or 134 °C (273 °F). Exceeding the temperature limit can cause the drill to lose hardness.

Recommended speed

Approx. 1000 rpm

Interactions

- Cements containing eugenol can affect the curing of luting cements.
- Avoid the use of materials containing eugenol.
- Do not use chlorhexidine preparations when rinsing the root canal.
- Moisture or chemical vapors (e.g. from acids) can cause corrosion.
- Protect the material against moisture and chemical vapors.

Residual risks/side effects

- We place great importance on the quality of our products. However, please check the product for damage and visible contamination before use.
- In the case of extreme loads, when phosphate or glass ionomer cements are used with tapered root posts, the root post can be pressed into the cement and the tooth can fracture (wedging effect). Using these cement classes in combination with LuxaPost should fundamentally be avoided. Glass fiber root posts should be luted exclusively with composite cement.
- A low residual risk of unknown side effects cannot be ruled out.
- Hypersensitivity and/or allergic reactions cannot be ruled out in individuals who are correspondingly predisposed.
- Improper use of the product may lead to excessive loss of tooth substance as well as to injuries to the skin or mucous membranes, eye irritation, obstructed or irritated airways, choking, and infections.

Warnings/precautions

- For dental use only.
- Avoid contact with the skin! In the event of accidental contact with the skin, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer (productsafety@mtcompanies.com) and to the responsible registration authorities.

Product performance characteristics

Flexural strength: ≥ 1200 MPa, flexural modulus: ≤ 50 GPa

Storage/disposal

- Store at 2 – 25 °C (36 – 77 °F) in a dry place, free from dust.
- Protect from chemical vapors.
- Do not use after the expiration date.
- Observe the requirements of the hygiene plan and the respective manufacturer's instructions for devices and materials used.
- Disposal must comply with national regulations.

Composition

Posts: Ceramic fibers, glass fillers, Bis-GMA, UDMA, HDDMA, additives.

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