LuxaPrint Ortho

DMG Digital Enterprises SE
Elbgaustraße 248 · 22547 Hamburg · Germany · www.dmg-dental.com
092864/2023-12

Gebrauchsinformation

Deutsch

Produktbeschreibung

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Instructions for use English

Product description

Light-curing resin based on (meth)acrylate for the production of dental devices. For DLP / SLA printers that work at wavelengths of 385 nm or 405 nm.

Devices with material parameters already defined in the system are recommended. You will find a list of validated devices on our website. Select material variants according to the wavelength used by the printer (see »Packaging«).

Intended purpose

Manufacturing of surgical templates and orthodontic devices with a wearing period of up to 24h by means of 3D printing.

Limitations of use

Only insert LuxaPrint Ortho intra-orally in a fully polymerized state.

Indication

- Surgical template: Tooth loss, navigational implant placement
- Orthodontic appliances: Stabilization of the tooth's position after the completion of active orthodontic treatment

Contraindication

Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended user

Dentist, dental technician

Notes for use

- During handling (e.g. transferring), the material should be exposed to daylight / ambient light only for as short a time as possible, since lengthy exposure to light can lead to unwanted curing of the material.
- Wear nitrile gloves, protective clothing and goggles when working with the material/printed object.

For optimum properties of the end product:

- Ensure that the devices used during the production process are in proper condition and comply with machine and material parameters.
- Ensure clean working procedures. Residues on the devices can cause defects on the printed object.
- Use a cleaning unit with material parameters defined in the system. You will find a list of validated devices on our website.
- Allow objects to cure completely.
- Use a light-curing unit with material parameters defined in the system. You will find a list of validated devices on our website.
- Minimum material thickness for structures: 1.5 mm.
- Maximum curing depth during post-curing: For solid objects and light-curing on both sides, the material can be up to 7 mm thick (at a curing depth of 3.5 mm).
- Thinner layers increase the accuracy of the printed object and increase printing time.
- Position the printed object as flat as possible on the printer's build platform. A flatter positioning increases accuracy and reduces printing time.

Recommended use

Note: Observe the device and material parameters.

Note: Maintain a processing temperature of 23 °C ± 2 °C (73.4 °F ± 3.6 °F).

Note: For optimal process reliability, scan the material's RFID code.

Finishing the printed object

Production steps	
Recommended drip-off time after raising the platform: 10 minutes	
Note: Finishing should take place immediately after the construction process where possible.	
1.	Carefully remove the printed object from the build platform.
2.	Cleaning: Clean the printed object with isopropanol (approx. 99 %) using an appropriate cleaning device and specified program.
2. (Alternatively)	1. Pre-clean the printed object with ethanol (approx. 96%) or isopropyl alcohol (approx. 99%) for max. 3 min in an ultrasonic cleaner (with a brush if necessary). 2. Optional: Clean openings, holes and gap areas with compressed air. Clean the printed object in a separate container with fresh ethanol (approx. 96%) or isopropyl alcohol (approx. 99%) for max. 2 min.
3.	Post-curing: Post-cure the printed object using the appropriate light-curing unit and the specified program.
4.	After removing the support, finish the surface of the printed object using the appropriate instruments and polishing paste.

Preparation prior to use on the patient

Disinfection

According to the manufacturer's specifications, the following disinfectants may be used:

- PrintoSept-ID (on the basis of quart. ammonium salts)
- SprayActiv, alcoholic disinfectant spray (also contains didecyldimethyl-N-chloride)
- Dentavon (solution prepared from a granulate; contains penta-potassium-bis (peroxymonosulfate)-bis(sulfate), anionic surfactants, non-ionic surfactants, soap, phosphonate)

Autoclavability

Before use on the patient, the printed object can be autoclaved once.

Autoclaving parameters for steam sterilization:

Temperature: 134 °C / 273 °F at a pressure of 2 bar and a duration of 5 min.

Interactions

Residues of methacrylate-based materials (uncured material residues) may influence the setting behavior of silicone impression materials.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid skin contact with the liquid material and components prior to post-curing! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

Cured material	
Flexural strength	≥ 70 MPa
Transverse modulus	≥ 1 GPa
Shore-D hardness	≥ 60

Uncured material	
Dynamic viscosity at 50 Pa shear stress	≤ 2.5 Pa s
Light sensitivity to ambient light	≤ 5 min

Storage/disposal

- Store in a dry place at room temperature (15°–25 °C/59°–77 °F) and protected from light!
- Even low exposure to light can trigger polymerization.
- Always keep the container tightly sealed; carefully seal immediately after every use.
- Do not use after the expiration date!
- The cured material can be discarded together with household waste following consultation with the disposal company and in compliance with regulations.
- Incinerate the resin mixture as hazardous waste in accordance with official regulations.

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