Luxatemp Glaze & Bond

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English

Product description

Luxatemp Glaze & Bond is a light-curing, 1-component material based on multi-functional methacrylates for surface sealing and repair of temporary crowns and bridges, as well as for the pre-treatment of indirect plastic restorations.

Intended purpose

- Coating and repair of provisional restorations
- Primer/bonding agent for indirect restorations of PEEK, composite, and PMMA

Limitations of use

Do not use Luxatemp Glaze & Bond as a dentin or enamel adhesive

Indication

Restoration of oral function and esthetics:

- Major tooth defects that cannot be restored by direct restoration, e.g. extensive carious lesions, traumatic injuries or serious dental wear
- Endodontically-treated teeth that are susceptible to fractures
- Tooth loss
- Dental misalignments and deformities, e.g. individual or multiple diastemas and internal discoloration
- It is not possible to predict the functional or esthetic outcomes of direct restorations
- Occlusal modification

Contraindication

Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

People treated in the course of a dental procedure.

Intended user

Dentist, dental assistant staff and dental technician

Notes for use

- Light-curing units should have an output of 450 nm and should be checked regularly. Place the light as close as possible to the material. The light intensity should be at least 600 mW/cm².
- Luxatemp Glaze & Bond binds to teeth and films made of polycarbonate or polyacetate.
- Do not use Luxatemp Glaze & Bond as a dentin or enamel adhesive.
- Avoid contact with materials containing eugenol.
- Close the bottle immediately after use.

Recommended use

Surface coating

1. Apply a thin layer of the material without entrapped air bubbles to the required surfaces.
2. Allow the material to act for 20 s and then cure each unit for 40 s at a light intensity of ≥ 600 mW/cm² or for 20 s at a light intensity of ≥ 1000 mW/cm² with a halogen or LED lamp. Alternatively cure the entire restoration for at least 90 seconds in the light oven.

Note: To ensure that no oxygen inhibition layer is formed, the material should be light-cured within 2 min after being applied.

Adhesion for repair purposes

1. Roughen the surface to be repaired with a cutter.
2. Apply a thin layer of the material without entrapped air bubbles to the required surfaces.
3. Cure for 40 s at a light intensity of ≥ 600 mW/cm² or for 20 s at a light intensity of ≥ 1000 mW/cm² with a halogen or LED lamp.
4. Apply freshly mixed temporary material, e.g. Luxatemp Star, or a light-curing composite, e.g. LuxaFlow Star, to the relevant surfaces. Observe the manufacturer's instructions for use for the material used.

Bonding agent for restorations

1. Prepare the restoration according to manufacturer's specifications.
2. Apply a thin layer of the material without entrapped air bubbles to the required surfaces. Cure for 40 s at a light intensity of ≥ 600 mW/cm² or for 20 s at a light intensity of ≥ 1000 mW/cm² with a halogen or LED lamp.
3. Lute restoration with a suitable adhesive dual- or self-curing composite luting cement according to the manufacturer's instructions.

Interactions

- Leftover methacrylate-based materials may alter the setting behavior of the silicone impression and bite registration materials.
- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- Avoid contact with materials containing eugenol.
- Keep out of the reach of children!
- For dental use only!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.
- Do not inhale vapors.

Product performance characteristics

Generates a protective acrylate surface film

Bonding on composite, PMMA and PEEK surfaces: SBS bonding properties ≥ 10 MPa

Storage/disposal

- Store at 2 - 25 °C (36 - 77 °F) in a dry place.
- If stored in the refrigerator (2-8 °C/36 - 46 °F), prior to applying bring the material to room temperature (15 - 25 °C/59 - 77 °F).
- Do not use after the expiration date.
- Close the bottle immediately after each use and store out of the way of direct sunlight.
- Dispose of in accordance with national regulations.

Composition

MMA, DPPA, vinyl ester resin, BAPO, MMHE, polymeric aminobenzoate derivative, additives.

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