PermaCem

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English

Product description

PermaCem is a universal, adhesive, self-curing compomer luting cement for the cementation of indirect restorations in combination with an enamel-dentine bonding agent such as LuxaBond Total Etch (Total-Etch bonding system), or LuxaBond Universal and a metal or ceramic bonding agent.

Intended purpose

Definitive cementation of crowns, bridges, inlays and onlays made of metal, ceramics or composites.

Limitations of use

Do not use PermaCem with single-bottle dentine bonding agents.

Indication

Cementation of indirect restorations and prostheses for the treatment of large dental defects, tooth loss, malocclusions and malformations of the teeth.

Contraindication

- Do not apply PermaCem directly to the exposed dental pulp.
- Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended user

Dentist

Notes for use

- With the Smartmix version the material is applied directly without a dispensing device.
- With the Automix version the material is applied using a 1:1 dispenser.
- Cementation of combined restorations such as with telescopic units: Using a suitable separating agent, ensure that PermaCem cannot flow between primary and secondary parts.

Using the Automix Cartridge (for image, see »Handling«)

Use of the Automix Cartridge

- Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] all the way.

- Lift up the plastic lever [C], insert the Automix Cartridge and lock it using the plastic lever [C].

Note: Make sure that the notches on the Automix Cartridge line up with those on the Automix Dispenser.

Attaching a mixing tip

- Remove the cartridge cap [D] or used mixing tip after turning it 90° counterclockwise, and discard it.
- Attach a new mixing tip [E].

Note: In order to achieve an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available in refill packs.

Note: Make sure that the notches on the mixing tip line up with those on the Automix Cartridge.

- Twist the mixing tip clockwise by 90° to lock it.

Applying the material

- ▶ When using an Automix Cartridge for the first time, squeeze out a pea-sized amount of material and discard it.
- ▶ Mix the material in the mixing tip using the lever [F] on the Automix Dispenser and apply directly.

Note: After the application, leave the used mixing tip on the Automix Cartridge as a seal!

Removal of the Automix Cartridge

- Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] all the way.
- Lift up the plastic lever [C] and remove the Automix Cartridge.

Using the Smartmix Syringe (for image, see »Handling«)

- Before attaching the mixing tip, remove the cap [A] or used mixing tip after turning it 90° counterclockwise, and discard it.
- Attach a new mixing tip [B].

Note: In order to achieve an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available in refill packs.

Note: Make sure that the notches on the mixing tip line up with those on the Smartmix Syringe [C].

- Twist the mixing tip clockwise by 90° to lock it.
- The material is mixed in the mixing tip when ejected and can then be applied directly.

Note: When using a Smartmix Syringe for the first time, squeeze out a pea-sized amount of material and discard it!

Note: After application, leave the used mixing tip on the Smartmix Syringe as a seal!

Disinfection of the Automix Dispenser

- The Automix Dispenser can be reused after disinfection, although it should be replaced if there are visible signs of wear. Disinfection can be carried out using a commercially available dipping disinfectant. The Automix Dispenser is also autoclavable.
- Dirty or contaminated mixing tips used as seals can be cleaned with disinfectants normally used in the practice.
- DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute).

Recommended use

- Prepare the clean and dry restoration as follows:
 - Etchable silicate ceramics: Etch the restoration with 5 % hydrofluoric acid in accordance with the instructions from the manufacturer and rinse thoroughly with water spray. Subsequently dry with oil-free and water-free air, and pre-treat the surface using a suitable silanizing agent (e.g. Vitique Silane or LuxaBond Universal, DMG) in accordance with the manufacturer's instructions.
 - Metals and metal alloys: Pretreat in accordance with the manufacturer's instructions. Unless described otherwise, sandblast the inside of the restoration with aluminum oxide, clean with alcohol, dry with oil-free and water-free air and apply a suitable metal bonding agent (e.g. LuxaBond Universal).
 - Oxide ceramics: Pretreat in accordance with the manufacturer's instructions. Ensure that there is sufficient parallel-wall retention. Unless described otherwise, sandblast the inside of the restoration with aluminum oxide, clean with alcohol, dry with oil-free and water-free air and apply a suitable ceramic bonding agent (e.g. LuxaBond Universal) according to the manufacturer's instructions.
- Remove any residue of the temporary cement from the prepared tooth, clean thoroughly with water spray and dry with oil-free and water-free air. In doing so, ensure that the dentin retains a little residual moisture and does not become too dry.

Note: It is advisable to use a rubber dam.

- Apply a suitable bonding agent to the prepared tooth according to manufacturer's instructions. DMG recommends LuxaBond Total Etch (Total-Etch-System) or LuxaBond Universal.
- Apply PermaCem in the pre-treated restoration or on the prepared tooth.
- Insert the restoration within 60 seconds after mixing has begun and attach to the dry, prepared tooth using slight pressure.
- Instruct the patient to adopt habitual occlusion.
- Remove cement residues within 60 to 90 seconds after inserting the restoration with the aid of a microbrush, brush, foam pellet or scaler. Remove interdental excess carefully using dental floss.
- Allow the material to cure. The restoration can support a full load after 360 seconds.

Timing

Working time (including mixing time)*	≤ 60 s
Maximum intraoral curing time (including mixing time)	360 s
Removal of excess material	60 - 90 s

* The specified times are valid at a room temperature of 23 °C and a normal relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times.

Interactions

- Residues of methacrylate-based materials may alter the setting behavior of silicone impression and bite registration materials.
- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

Complies with ISO 4049:2019, Type 2, Class 1 (including shade stability and radiopacity of ≥ 190 % Al).

Storage/disposal

- Store in a dry place at a room temperature of 2 - 25 °C (36 - 77 °F).
- After application, leave the used mixing tip on the Automix Cartridge / Smartmix Syringe as a seal!
- Do not use after the expiration date.
- Disposal must comply with national regulations.

Composition

Dental glass, bis-GMA, EBPADMA, TEDMA, TMPTMA, MMHE, SiO₂, additives. Inorganic filling materials: approx. 49 vol% (0.02 – 4 µm).

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