# Ecosite One

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## Product description

Ecosite One is a light-curing, radiopaque composite for direct fillings, especially suitable for use in the posterior region.

## Intended purpose

- Fillings in the posterior region (Classes I and II, including individual cusps)
- Restorations of class V
- Core build-up
- Extended fissure sealing
- Restorations in deciduous teeth

## Limitations of use

Application of the material is contraindicated if dry isolation or the recommended application technique is not possible.

## Indication

- Loss of tooth substance due to caries, fractures, dental wear or developmental disorders

## Contraindication

- Do not use directly on the exposed pulp.
- Do not use the material where there are known allergies to any of the components or contact allergies.

## Patient target group

Persons receiving treatment as part of a dental procedure.

## Intended user

Dentist

## Notes for use

- For optimum polymerization, position the light-curing unit as close to the restorative material as possible and observe the recommendations relating to the light-curing time and the light intensity (see table). Observe the manufacturer's instructions for the light-curing unit used.
- Follow the manufacturer's instructions for other materials that are used with Ecosite One.
- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

## Recommended use

1. Clean the tooth. Remove all cleaning residue with water spray.

**Note:** DMG recommends the use of a rubber dam.

2. Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Ensure that the areas close to the pulp are protected.
3. Condition the dental enamel and apply the bonding agent according to the manufacturer's instructions for the material used. DMG advises using the light-curing bonding agent Ecosite Bond.
4. If necessary, apply a base of low-viscosity composite or compomer. We recommend using a matrix.
5. Insert Ecosite One with a maximum layer thickness of 3 mm into the cavity. Apply taller restorations one layer at a time. Light-cure each layer separately.

Light intensity	≥ 1000 mW/cm <sup>2</sup>	≥ 600 mW/cm <sup>2</sup>
Exposure time	10 s	20 s
Max. layer thickness	3 mm	3 mm
Wavelength	420 – 480 nm	420 – 480 nm

**Note:** Additional buccal or lingual light-curing is advisable. Remove metal matrices beforehand.

6. Removal of any excess material, contouring, and finishing are done using diamond finishing tools and flexible polishing disks of decreasing grain size. Ecosite One can be given a high shine using polishing paste together with polishing cups or disks.

## Handling of the Safetip Dispenser Pro

1. Pull the quick-lock mechanism all the way back.
2. Insert the Safetip into the Safetip Dispenser Pro, and then guide the quick-lock mechanism back into place until the Safetip is clearly engaged.
3. Turn the Safetip dispensing nozzle in the desired direction, remove the cap from the Safetip and dispense the material by pressing evenly. Do not exert too much force when pressing.
4. After dispensing, pull the quick-lock mechanism back again and remove the Safetip.
5. The Safetip Dispenser Pro should be cleaned and disinfected after use.

## Disinfection of the Safetip Dispenser Pro

The Safetip Dispenser Pro is reusable after disinfection, although it should be replaced if visible signs of wear appear.

Disinfection can be carried out using a commercially available wiping or immersion disinfectant. DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute).

## Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid these conditions.
- The use of oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

## Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

## Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Reuse: Standard surface disinfection as required. No specific maintenance required.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

## Product performance characteristics

Ecosite One complies with ISO 4049:2019 (Type 1, Class 2, Group 1); polymerization depth: ≥ 3 mm by using a light-curing time of 20 seconds with a light intensity of ≥ 600 mW/cm<sup>2</sup> or by using a light-curing time of 10 seconds with a light intensity of ≥ 1000 mW/cm<sup>2</sup>

## Storage/disposal

- Store in a dry place at 2 to 25°C (36 to 77°F)!
- Dispose of in accordance with national regulations.

## Composition

Dental glass, EBPADMA, SiO<sub>2</sub>, neopentyl glycol propoxylate diacrylate, Bis-GMA, UDMA, additives. Inorganic filling materials: approx. 65 vol%, 0.02 – 0.7 µm.

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