# LuxaPrint Ortho Comfort

DMG Digital Enterprises SE  
 Elbgaustraße 248 · 22547 Hamburg · Germany · www.dmg-dental.com  
 092878/2024-04

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### Produktbeschreibung

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IBOMA, alip. Urethanacrylat, alip. Urethanmethacrylat, EBPADMA, TPO, Additive

## Instructions for use

English

### Product description

Light-curing resin based on (meth)acrylate for the production of dental devices. For DLP/SLA printers that work at a wavelength of 385 nm.

Devices with material parameters already defined in their systems are recommended. You will find a list of validated devices on our website. Select material variants according to the wavelength used by the printer (see »Packaging«).

### Intended purpose

Manufacturing of occlusal splints using 3D printing.

### Limitations of use

Only insert LuxaPrint Ortho Comfort intraorally in a fully polymerized state.

### Indication

Bruxism and craniomandibular dysfunction (CMD)

### Contraindication

Do not use the material where there are known allergies to any of the components or contact allergies.

### Patient target group

People to be treated in the course of a dental procedure.

### Intended user

Dentist, dental technician

### Notes for use

- During handling (e.g. transferring), the material should be exposed to daylight/ambient light only for as short a time as possible, since lengthy exposure to light can lead to unwanted curing of the material.
- Wear nitrile gloves, protective clothing and goggles when working with the material/printed object.

### For optimum properties of the end product:

- Ensure that the devices used during the production process are in proper condition and comply with machine and material parameters.
- Ensure clean working procedures. Residues on the devices can cause defects on the printed object.
- Use a cleaning unit with material parameters defined in the system. You will find a list of validated devices on our website.
- Allow objects to cure completely.
- Use a light-curing unit with material parameters defined in the system. You will find a list of validated devices on our website.
- Maximum curing depth during post-curing: For solid objects and light-curing on both sides, the material can be up to 7 mm thick (at a curing depth of 3.5 mm).
- Thinner layers increase the accuracy of the printed object and increase printing time.
- Position the printed object as flat as possible on the printer's build platform. A flatter positioning increases accuracy and reduces printing time.

### Recommended use

**Note:** Observe the machine and material parameters.

**Note:** Maintain a processing temperature of 23 °C ± 2 °C (73.4 °F ± 3.6 °F).

**Note:** Scan the RFID code of the material for a higher level of process security.

### Finishing the printed object

Production steps	
Recommended drip-off time after raising the platform: 10 min	
<b>Note:</b> Finishing should take place immediately after the construction process where possible.	
1.	Carefully remove the printed object from the build platform.
2.	Cleaning: Clean the printed object with isopropanol (approx. 99 %) using an appropriate cleaning device and specified program.
2. (Alternatively)	1. Pre-clean the printed object with isopropyl alcohol (approx. 99 %) for max. 3 min in an ultrasonic cleaner (with a brush if necessary). 2. Optional: Clean openings, boreholes and gap areas using pressurized air. Clean printed object in a separate container using clean isopropanol (approx. 99 %) for max. 2 min.
3.	Post-curing: Post-cure the printed object using the appropriate light-curing unit and the specified program.
4.	After removing the support, finish the surface of the printed object using the appropriate instruments and polishing paste.

### Interactions

Residues of methacrylate-based materials (uncured material residues) may influence the setting behavior of silicone impression materials.

### Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

### Warnings/precautions

- Keep out of the reach of children!
- Avoid skin contact with the liquid material and components prior to post-curing! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Avoid breathing in the fumes.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.
- Do not use heat-based disinfection or sterilization methods. Otherwise, the workpiece may be deformed.

### Product performance characteristics

Cured material	
Shore-D hardness	≤ 72

Uncured material	
Dynamic viscosity at 50 Pa shear stress	≤ 3.0 Pa s
Light sensitivity to ambient light	≤ 5 min

### Storage/disposal

- Store in a dry place at room temperature (15 – 25 °C/59 – 77 °F) and protected from light!
- Even low exposure to light can trigger polymerization.
- Always keep the container tightly sealed; carefully seal immediately after every use.
- Do not use after the expiration date!
- The cured material can be discarded together with household waste following consultation with the disposal company and in compliance with regulations.
- Incinerate the resin mixture as hazardous waste in accordance with official regulations.

### Composition

IBOMA, aliphatic urethane acrylate, aliphatic urethane methacrylate, EBPADMA, TPO, additives

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