Flairesse Bleaching Gel CP 5%

Flairesse Bleaching Gel CP 10%

DMG Chemisch-Pharmazeutische Fabrik GmbH
Elbgaustraße 248 · 22547 Hamburg · Germany · www.dmg-dental.com
092791/2025-02

Gebrauchsinformation

Deutsch

Produktbeschreibung

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English

Product description

Flairesse Bleaching Gel CP is a gel containing 5 % or 10 % carbamide peroxide (equivalent to 1.8 % or 3.6 % hydrogen peroxide) for the gentle bleaching of teeth. After the dentist has given instructions, the product can be used independently by the patient, under the dentist's supervision.

Intended purpose

- Whitening of intrinsic and extrinsic discolorations due to various causes.
- Treatment of developmental enamel defects in minors.

Limitations of use

- Whitening treatment is not effective against discolorations caused by metals (e.g. amalgam or silver points).
- The material is not suitable for heavy smokers.
- The material is not suitable for patients with a high alcohol intake.

Indication

- Intrinsic and extrinsic discolorations.
- Developmental enamel defects.

Contraindication

- Do not use the material where there are known allergies to any of the components or contact allergies.
- Do not use the material with patients who are pregnant or breastfeeding.
- Do not use the material in cases of periodontal disease.
- Do not use the material in cases of exposed dentin, defective filling margins or untreated caries.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended user

Initial treatment is carried out by the dentist at the dental practice. The patient is then able to use the bleaching product independently under the dentist's supervision.

Notes for use

- Demineralized areas (white spots) may become more visible at the start of treatment. Any difference in shade compared to the healthy tooth enamel will even out during the course of treatment.
- Treatment for developmental enamel defects may require a longer treatment period.

Using the syringe

- The gel is supplied in a 5 ml syringe.
- The syringe must be resealed with its tip cap after each application.
- Only use in combination with a bleaching tray.

Recommended use

Pretreatment

1. Thorough examination of the patient's mouth to identify any defective (provisional) restorations, untreated caries or exposed dentin. These areas must be excluded from the bleaching treatment or treated accordingly beforehand, as the patient may otherwise experience hypersensitivity.
2. Professional tooth cleaning to remove extrinsic stains and calculus prior to treatment.

Note: If irritation of the gingiva occurs, treatment must be delayed until the irritation subsides to avoid further sensitization.

3. Assessment of the type and extent of the discoloration involved and whether or not tooth bleaching is indicated.
4. Assessment of the initial tooth shade using a suitable shade guide (e.g. Vita classical®). Photo documentation of the initial tooth shade can help patients to visualize their treatment success.

Note: Patients should be informed that the tooth shade of existing restorations will not change and that a shade difference to the surrounding enamel may occur.

5. If the patient is unable to follow the treatment instructions competently, the bleaching treatment should not be recommended.

Bleaching treatment

1. A bleaching tray is made for this purpose using suitable materials and in accordance with the relevant manufacturer's instructions.
2. Open the syringe by unscrewing the cap. Keep the cap in a safe place.

Note: The syringe must be resealed with its cap after use.

3. Dispense a pinhead-sized amount into each reservoir of the tray for the areas to be treated.

Note: Check whether the amount is correct by looking at how the material is distributed in the reservoir after the tray has been inserted. It may be necessary to adjust the amount for the next application.

4. Insert the filled bleaching tray. Remove any excess material with a soft cloth or a soft toothbrush.
5. Wear the tray for up to 8 hours, depending on the dentist's specific recommendation. Repeat the treatment daily until the desired tooth shade is achieved. Avoid eating and drinking during treatment.

Note: It is possible to pause the treatment for several days without jeopardizing subsequent treatment.

Note: It is important to avoid staining foods and drinks (e.g. wine, coffee, cola) as well as acidic foods (e.g. citrus fruits) and drinks (e.g. lemonade) throughout the entire treatment cycle.

6. Remove the tray and rinse the mouth several times with water. Do not swallow the water! Any gel remaining on the teeth can be carefully removed with a soft toothbrush.
7. Clean the tray with cold water and a soft toothbrush. Store the tray in the tray box provided.

Interactions

Carbamide peroxide can have a negative effect on the radical polymerization of resins. Therefore, adhesive dental materials must not be used for at least two weeks after the bleaching treatment has ended.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Children and adolescents should only use the product under adult supervision.
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! Risk of loss of sight! In case of accidental contact with eyes, immediately and carefully rinse with water for several minutes. If contact lenses are worn, remove these if possible. Continue rinsing.
- If exposure has occurred or is suspected, seek medical advice/attention.
- Contains alkali.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.
- Patient should be told to avoid swallowing large amounts of gel.

- If there is risk of swallowing or inhaling the bleaching tray contents, treatment should not take place.
- Patient should be advised to make an appointment if any symptoms occur.

Product performance characteristics

Complies with ISO 28339:2021.

Storage/disposal

Store in a dry place at 2 to 25 °C (36 to 77 °F).

Composition

Flaïresse Bleaching Gel CP 5 %: Water, glycerin, sodium hydroxide, carbamide peroxide, polyacrylic acid (cross-linked), xylitol, potassium nitrate, additives.

Flaïresse Bleaching Gel CP 10 %: Water, glycerin, sodium hydroxide, carbamide peroxide, polyacrylic acid (cross-linked), xylitol, potassium nitrate, additives.

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