LuxaPrint Ortho Flex NX

DMG Digital Enterprises SE

Elbgaustraße 248 · 22547 Hamburg · Germany · www.dmg-dental.com
093196/2025-06

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Produktbeschreibung

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Aliphatisches Urethanacrylat, IBOMA, TPO-L, Additive.

Instructions for useEnglish

Product description

Light-curing resin based on (meth)acrylate for the production of dental devices. For DLP / SLA printers that work at wavelengths of 385 nm and 405 nm. Devices with material parameters already defined in their systems are recommended. A list of validated devices can be found on our website (Digital Workflow—Resins—LuxaPrint: Validated Systems).

Intended purpose

Manufacturing of bleaching trays and transfer trays (indirect bracket bonding technique) using 3D printing.

Limitations of use

Only insert LuxaPrint Ortho Flex NX intraorally in a fully polymerized state.

Indication

- Bleaching trays: Tooth discoloration
- Transfer splints: Indirect bracket bonding in orthodontic treatment.

Contraindication

Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

People to be treated in the course of a dental procedure.

Intended user

Dentist, dental technician

Notes for use

- During handling (e.g. transferring), the material should be exposed to daylight / ambient light only for as short a time as possible, since lengthy exposure to light can lead to unwanted curing of the material.
- Wear nitrile gloves, protective clothing and goggles when working with the material/printed object.

For optimum properties of the end product:

- Ensure that the devices used during the production process are in proper condition and comply with machine and material parameters.
- Ensure clean working procedures. Residues on the devices can cause defects on the printed object.
- Use a cleaning unit with material parameters defined in the system. A list of validated devices can be found on our website (Digital Workflow—Resins—LuxaPrint: Validated Systems).
- Allow objects to cure completely.
- Use a light-curing unit with material parameters defined in the system. A list of validated devices can be found on our website (Digital Workflow—Resins—LuxaPrint: Validated Systems).
- Thinner layers increase the accuracy of the printed object and increase printing time.
- Position the printed object as flat as possible on the printer’s building plate. A flatter positioning increases accuracy and reduces printing time.
- Maximum curing depth during post-curing: For solid objects and light-curing on both sides, the material can be up to 7 mm thick (at a curing depth of 3.5 mm).

Recommended use

- Note: Observe the device and material parameters.
- Note: Maintain a processing temperature of 23 °C ± 2 °C (73.4 °F ± 3.6 °F).
- Note: Scan the RFID code of the material for a higher level of process security.

Finishing the printed object

Working steps	
Recommended drip-off time after raising the platform: 10 min	
Note: Finishing should take place immediately after the construction process where possible.	
1.	Carefully remove the printed object from the build platform.
2.	Cleaning: Clean the printed object with isopropanol (approx. 99 %) using an appropriate cleaning device and specified program.
2. (Alternatively)	1. Pre-clean the printed object with ethanol (approx. 96%) or isopropyl alcohol (approx. 99%) for max. 3 min in an ultrasonic cleaner (with a brush if necessary). 2. Optional: Clean openings, holes and gap areas with compressed air. 3. Clean the printed object in a separate container with fresh ethanol (approx. 96%) or isopropyl alcohol (approx. 99%) for max. 2 min.
3.	Post-curing: Post-cure the printed object using the appropriate light-curing unit and the specified program.
4.	After removing the support, finish the surface.

Interactions

Residues of methacrylate-based materials (uncured material residues) may influence the setting behavior of silicone impression and bite registration materials.

Residual risks/side effects

- We place great importance on the quality of our products. However, please check the product for damage and visible contamination before use.
- There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- Avoid skin contact with the liquid material and components prior to post-curing! In the event of accidental contact with the skin, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Avoid breathing in the fumes.
- Follow the manufacturer’s instructions for other products that are used with the material/ materials.
- Observe the safety data sheet.
- Do not use the material where there are known allergies to any of the components or contact allergies.
- Serious incidents involving this product must be reported to the manufacturer (productsafety@mtcompanies.com) and to the responsible registration authorities.
- Do not mix products with different lot numbers.

Product performance characteristics

Cured material	
Shore-D hardness	≤ 60

Uncured material	
Dynamic viscosity at 50 Pa shear stress	≤ 4.0 Pa s
Light sensitivity to ambient light (500 lux)	≤ 5 min

Storage/disposal

- Store in a dry place at room temperature (15 – 25 °C/59 – 77 °F) and protected from light!
- Even low exposure to light can trigger polymerization.

- Always keep the container tightly sealed; carefully seal immediately after every use.
- Do not use after the expiration date!
- The cured material can be discarded together with household waste following consultation with the disposal company and in compliance with regulations.
- Incinerate the resin mixture as hazardous waste in accordance with official regulations.

Composition

Aliphatic urethane acrylate, IBOMA, TPO-L, additives.

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