Constic

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English

Product description

Constic is a self-etching, self-adhesive, radiopaque, flowable composite which is cured with light. The composite is immediately ready for use because the preparatory steps of etching, priming and bonding are contained in Constic.

Intended purpose

- Small restorations of class I and III
- Restorations of class V
- Base material for class I and II
- Fissure sealing
- Extended fissure sealing
- Repairs of composite restorations
- Modifications to temporaries and long-term temporaries
- Blocking out and filling of undercuts in cavities
- Fillings for small occlusal deciduous tooth cavities

Limitations of use

Application of the material is contraindicated if dry isolation or the recommended application technique is not possible. Do not use Constic on etched dentin.

Indication

- Loss of tooth substance due to caries, trauma or development disorders
- Use in molars and premolars with pronounced deep fissures and increased caries risk: Fissure sealing

Contraindication

- Do not apply Constic to the exposed pulp.
- Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended user

Dentist

Notes for use

- Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 600 mW/cm². Observe the manufacturer's instructions for the light-curing unit.
- Position the light-curing unit as close to the material as possible.
- Do not remove the oxygen inhibited layer because this is required for the bond with the next layer.
- When using Constic in conjunction with other materials, follow the applicable manufacturer's instructions.
- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.
- Disinfection: The brush holder is reusable after disinfection, although it should be replaced if visible signs of wear, e.g. scratches, appear. Disinfection can be carried out using a commercially available wiping or dipping disinfectant. DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute).

Recommended use

Small restorations of class I, III and V and base material for class I and II

1. Determine tooth color and select the appropriate shade.
2. Prepare the cavity in accordance with the general rules for adhesive technique. Bevel the enamel edges if necessary. Protect pulp in deep cavities with appropriate preparations.

Note: DMG recommends the use of a rubber dam.

3. Clean tooth, remove all cleaning residue with water spray.

Note: Cleaning with alcohol may cause overdrying of the dentin.

4. Dry the tooth using water-free and oil-free air in order to avoid overdrying the dentin. A moist layer must remain on the surface of the tooth.
5. Apply the material onto the cavity surface with the aid of the Luer-Lock Tips by pressing the syringe and massage a thin layer (= 0.5 mm) into the entire surface of the cavity wall for 25 seconds using the brush.
6. Remove any excess, if necessary.
7. Light-cure for 20 seconds.
8. Layer the material as follows:
 - Small restorations of class I, III and V: Build up layers with a maximum thickness of 2 mm with the material and light-cure each layer for 20 seconds at a light intensity of at least 1000 mW/cm² or 40 seconds at a light intensity of at least 600 mW/cm².
 - Base material for class I and II: Use the material to build up 1 layer with a maximum thickness of 1 mm and light-cure for 20 seconds.

Note: It is not necessary to massage the individual layers more than once.

9. Remove excess, contour the filling surface and finish using suitable instruments. Check the occlusion and correct if necessary.
10. Use a polishing paste to polish the filling surface and neighboring surfaces. Use a polishing cup or polishing disc for polishing work.

Fissure sealing and extended fissure sealing

1. For extended fissure sealing, widen particularly narrow and vial-shaped fissures.
2. Clean enamel with fluoride-free paste, remove all cleaning residue with water spray and dry the tooth with water-free and oil-free air. In doing so, leave the enamel surface moist.
3. Determine tooth color and select the appropriate shade.

Note: DMG recommends enamel etching in addition to cleaning. To do this, etch the enamel for 15 – 30 seconds with 37% phosphoric acid gel, rinse for 15 seconds with water and dry with oil-free and water-free air. In doing so, leave the enamel surface moist.

4. Press the syringe to apply the material onto the surface of the fissure or cavity wall with the aid of the Luer-Lock Tips and use the brush to massage a thin layer (= 0.5 mm) into the entire surface for 25 seconds.
5. Remove any excess, if necessary.
6. Light-cure for 20 seconds.
7. If necessary, use the material to build up layers with a maximum thickness of 2 mm and light-cure each layer for 20 seconds at a light intensity of at least 1000 mW/cm² or 40 seconds at a light intensity of at least 600 mW/cm².

Note: It is not necessary to massage the individual layers more than once.

8. Check the occlusion and correct if necessary.

Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid materials containing eugenol, moisture and air containing oil!
- The use of mouthwashes, plaque indicators and chlorhexidine may lead to discoloration.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

The product meets the requirements for mechanical properties according to DIN EN ISO 4049:2019, Type 1, Class 2, Group 1.

Polymerization depth ≥ 2 mm with a 20 second light-curing time when the light intensity of the lamp is ≥ 1000 mW/cm²

Polymerization depth ≥ 2 mm with a 40 second light-curing time when the light intensity of the lamp is ≥ 600 mW/cm²

Storage/disposal

- Store in a dry place at 2 to 25°C (36 to 77°F)!
- Disposal must comply with national regulations.

Composition

Dental glass, bis-GMA, EBPADMA, alip. EDA, UDMA, MDP, HEMA, TEDMA, additives. Inorganic filling materials: approx. 43 vol%, 0.02 – 2.3 μ m.

 Contains: 2-ethylhexyl 4-(dimethylamino)benzoate

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