# DMG Etching Gel

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English

### Product description

DMG Etching Gel is a universal etching gel (37 % phosphoric acid gel) with optimum flowability and thixotropy. DMG Etching Gel delivers reliable conditioning results and, due to its blue coloring, contrasts with the tooth structure. The Jumbo Syringe (50 ml) allows an empty 2 ml Luer-Lock Syringe to be refilled.

### Intended purpose

Etching of enamel and dentin

### Limitations of use

DMG Etching Gel is not suitable for etching ceramics.

### Indication

Loss of tooth substance due to caries, fractures, tooth abrasion or developmental disorders.

## Contraindication

Do not use the material where there are known allergies to any of the components or contact allergies.

## Patient target group

Persons receiving treatment as part of a dental procedure.

## Intended user

Dentist

## Notes for use

- Before continuing the dental procedure – e.g. by application of a bonding agent – the etched surface should not be touched and should not be contaminated, for example, with saliva. Should contamination occur after drying, the surface should be re-etched for approximately 10 seconds.
- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

## Recommended use

1. Clean the tooth prior to preparation. Remove all cleaning residue with water spray. Prepare the cavity in accordance with the general rules for adhesive technique. Bevel the enamel edges. It is advisable to use a rubber dam. In deep cavities in areas close to pulp provide protection in the form of a suitable liner, e.g. by using a small quantity of calcium hydroxide and a thin layer of glass ionomer cement.
2. Apply DMG Etching Gel to the enamel edges of the cavity and allow it to act for 20 – 60 seconds. When using the Total-Etch technique start the etching process with the enamel and then expand this to the dentin, so that the application time on the enamel is between 20 and 60 seconds, whilst the application time on the dentin does not exceed 15 seconds. The etched enamel should have a chalky white appearance. If this is not the case, etching must be repeated. Please observe the instructions for use applicable to the materials used with the DMG Etching Gel.
3. Rinse with water for a minimum of 15 seconds. Then dry using oil-free and water-free air.

## Filling the 2 ml Luer-Lock Syringe

► **Caution!** Risk of chemical burns. Wear gloves, goggles and a face mask.

1. Remove the sealing caps from the 50 ml Jumbo Syringe and the 2 ml Luer-Lock Syringe.
2. Join the two syringes by turning them in a clockwise direction. Ensure that the two syringes are firmly connected.

**Note:** Do not overtighten the syringes.

► **Caution!** Risk of chemical burns. Do not exceed the plunger stop of the 2 ml Luer-Lock Syringe. Ensure that the syringes are not pointed at people during filling. Do not dispense DMG Etching Gel from the 2 ml Luer-Lock Syringe back into the 50 ml Jumbo Syringe.

3. Slowly and smoothly retract the plunger of the 2 ml Luer-Lock Syringe to fill the 2 ml Luer-Lock Syringe. At the same time, carefully press the plunger of the 50 ml Jumbo Syringe with the other hand.

**Note:** Do not overfill the 2 ml Luer-Lock Syringe.

4. Turn syringes counterclockwise to separate them.
5. Place the sealing cap on the 50 ml Jumbo Syringe.
6. Place the Luer-Lock Tip on the 2 ml Luer-Lock Syringe.
7. Apply the DMG Etching Gel as described in the section »Recommended use«.

## Interactions

DMG Etching Gel should not be applied on acid-soluble linings.

## Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out

## Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

## Product performance characteristics

Etches dentin and enamel optimally because of its acid content and thixotropic behavior.

## Storage / disposal

- Store in a dry place at room temperature 15-25 °C (59-77 °F).
- Do not use after the expiration date.
- Disposal must comply with national regulations.

## Composition

Water, phosphoric acid, SiO<sub>2</sub>, additives.

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