# Ionosit Baseliner

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English

### Product description

Ionosit Baseliner is a light-curing, radiopaque, compomer base liner. The controlled expansion of Ionosit Baseliner counteracts the shrinkage of a composite filling.

### Intended purpose

- Base liner for composite restorations
- Blocking of undercuts

### Limitations of use

Application of the material is contraindicated if dry isolation or the recommended application technique is not possible.

### Indication

Loss of tooth substance due to caries, trauma, tooth attrition/abrasion or development disorders

### Contraindication

- Do not use directly on the exposed pulp.
- Do not use the material where there are known allergies to any of the components or contact allergies.

### Patient target group

Persons receiving treatment as part of a dental procedure.

## Intended user

---

Dentist

## Notes for use

---

- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.
- Light-curing units should emit at 450 nm and should be checked regularly.
- The light intensity should be at least 600 mW/cm<sup>2</sup>. Place the light-curing unit as close as possible to the material.
- Ionosit Baseliner adheres to metal instruments. Therefore, clean the instruments with a paper towel before the material sets.
- Ionosit Baseliner is light-sensitive. Protect the material from premature activation by ambient light and close the syringe immediately after use.
- Do not remove the oxygen inhibition layer because it is required for bonding to the next layer.
- The use of a matrix is recommended. Follow the manufacturers' instructions for other materials used.

## Recommended use

---

1. Before preparation, clean the tooth and remove all cleaning residues with water spray. It is recommended to use a rubber dam. Avoid contamination with saliva during application.
2. Prepare the cavity in accordance with the general rules for an adhesive technique. Protect pulp in deep cavities with suitable materials.
3. Pre-treat the cavity according to the preferred technique (e.g. Total Etch). Observe the manufacturer's instructions for the material used.
4. Apply a bonding agent (e.g. Ecosite Bond) in accordance with the manufacturer's instructions.
5. Ionosit Baseliner can be used both as a conventional base liner (coating the floor of the cavity) and as a coating for the entire cavity (stress-breaking function). Apply Ionosit Baseliner with a brush or ball-tipped instrument in a thickness of about 1 mm.
6. Cure 20 seconds using at least 1000 mW/cm<sup>2</sup> or 40 seconds using at least 600 mW/cm<sup>2</sup> with a suitable light unit.
7. Place the composite filling material according to the manufacturer's instructions. DMG recommends the use of Ecosite Elements.

## Interactions

---

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid materials containing eugenol, moisture and air containing oil!
- The use of cationic oral rinses, plaque indicators and chlorhexidine may lead to discoloration. Leftover methacrylate-based materials may influence the setting behavior of silicone materials.

## Residual risks/side effects

---

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

## Warnings/precautions

---

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

## Product performance characteristics

---

Compressive strength  $\geq$  100 MPa

## Storage/disposal

---

- Store in a dry place at room temperature 15-25 °C (59-77 °F).
- Do not use after the expiration date.
- Disposal must comply with national regulations.

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