Luxatemp Solar

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English

Product description

Luxatemp Solar is a light-curing composite with initial chemical curing for the fabrication of provisional crowns, bridges, inlays, onlays, partial crowns and veneers. The 2-component material is based on multi-functional methacrylates and it is automatically mixed.

Intended purpose

Fabrication of temporaries and long-term temporaries: crowns, bridges, inlays, onlays, partial crowns and veneers

Limitations of use

Do not use the material for permanent restorations

Indication

Restoration of oral function and esthetics:

- Major tooth defects that cannot be restored by direct restoration, e.g. extensive carious lesions, traumatic injuries or extensive dental wear
- Endodontically treated teeth that are susceptible to fractures
- Tooth loss
- Dental misalignments and deformities, e.g. individual or multiple diastemas and internal discoloration
- Cases for which it is not possible to predict the functional or esthetic outcomes of direct restorations
- Occlusal modification

Contraindication

Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended user

Dentist, dental assistant staff and dental technicians

Notes for use

- The material is applied using a 10:1 dispenser.
- When using an Automix Cartridge for the first time, squeeze out a pea-sized amount of material and discard it.

Recommended use

1. Before preparing the crown or bridge or prior to extraction of a tooth, take an impression using alginate, silicone or a thermoplastic impression material.

Alternatively, the temporary can be made using a laboratory-fabricated vacuum-formed stent or, in the case of single crowns, with a strip crown. Adhere precisely to the relevant manufacturer's instructions!

Note: In a silicone impression, undercuts should be compensated for and, if necessary, relief channels should be cut. To create a connection with the abutment teeth in the case of a gap in the molar region, cut a groove into the casting.

1. Dry the prepared teeth and compensate for any undercuts present in the preparation. Then lightly lubricate the cores and surrounding tissue plus any composite core build-ups using petroleum jelly or a similar separating medium.
2. Insert the Luxatemp Solar cartridge into the Automix Dispenser (see "Using the Automix Cartridge").

Note: To avoid bubbles, it is important that the end of the mixing tip always remains immersed in the material during application.

3. Dispense Luxatemp Solar into the impression or vacuum-formed stent. Then dispense the material onto the occlusal surfaces of the impression or vacuum-formed stent and fill up to the gingival areas, only overbuilding slightly.
4. Place the impression or vacuum-formed stent filled with Luxatemp Solar on the prepared teeth 1:10 min at the latest after mixing has begun and hold in place with moderate pressure.

Note: The setting process must be monitored intraorally (e.g. using a probe), as the temperature in the mouth has a critical effect on curing, and the provisional can only be removed during its elastic phase.

5. The provisional must be removed from the mouth during the elastic phase, i.e. during the period from 2 to 4 min after mixing has begun.

Note: It is essential that the provisional be light-cured for it to achieve its final hardness.

6. Light-curing: Expose each unit of the provisional to a halogen lamp for 20 seconds each or cure for at least 1 min in a polymerization unit.
7. Before finishing, remove the oxidation prevention layer from the surface of the provisional (e.g. with alcohol).
8. Prepare the provisional with cross-cut burs, flexible disks etc.
9. Check occlusion and adjust if necessary using suitable instruments.
10. If necessary, coat the surface of the provisional with a light-curing high-gloss glaze (e.g. Luxatemp Glaze & Bond) or manually polish it with a rubber polisher or small goat's hair brush.
11. Clean the core, dry it and thoroughly remove the separating medium.
12. Attach the provisional.

Note: Any temporary cement can be used for fixing the provisional. DMG recommends the use of automatically mixing luting cement. Cements that contain eugenol may prevent the polymerization of resin-based luting cements during the definitive luting

13. Check the occlusion and adjust accordingly.

Customizing and repairing a provisional

Note: When customizing and repairing a freshly-made provisional, it does not need to be roughened and a bonding agent is not required (continue with step 2).

1. Repairing a provisional that is being worn: Roughen the provisional slightly using a grinder or a sandblaster at the appropriate points. Afterwards apply a bonding agent (e.g. Luxatemp Glaze & Bond) in accordance with the manufacturer's instructions.
2. Apply Luxatemp Solar to all surfaces.
3. Press the parts together.
4. Light-cure the relevant part of the provisional for 20 seconds with a curing light and then process as normal.

Note: Alternatively, the provisional can also be customized or repaired using the light-curing, easy-flowing composite LuxaFlow Star or Luxatemp Star. Follow the corresponding instructions for use!

Timing

| | |
|-----------------------|---------------|
| Place into the mouth | 0:00–1:10 min |
| Remove from the mouth | 2:00–4:00 min |

Note: The indicated times apply at a room temperature of 23 °C/73 °F and a normal relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times.

| Light-curing | |
|---------------------|---------------------|
| Hand lamp | 20 seconds per unit |
| Polymerization unit | ≥ 1:00 min |

Using the Automix Cartridge (for image, see "Handling")

Use of the Automix Cartridge

1. Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] all the way.
2. Lift up the plastic lever [C], insert the Automix Cartridge and lock it using the plastic lever [C].

Note: Make sure that the notches on the Automix Cartridge line up with those on the Automix Dispenser.

Attaching a mixing tip

1. Remove the cartridge cap [D] or used mixing tip after turning it 90 ° counterclockwise, and discard it.
2. Attach a new mixing tip [E].

Note: Make sure that the notches on the mixing tip line up with those on the Automix Cartridge.

3. Twist the mixing tip clockwise by 90 ° to lock it.

Applying the material

- ▶ When using an Automix Cartridge for the first time, squeeze out a pea-sized amount of material and discard it.
- ▶ Mix the material in the mixing tip using the lever [F] on the Automix Dispenser and apply directly.

Note: After the application, leave the used mixing tip on the Automix Cartridge as a seal!

Remove the Automix Cartridge

1. Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] all the way.
2. Lift up the plastic lever [C] and remove the Automix Cartridge.

Disinfection of the Automix Dispenser

- The Automix Dispenser can be reused after disinfection, although it should be replaced if there are visible signs of wear.
- Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute). The Automix Dispenser is also autoclavable

Interactions

- Residues of methacrylate-based materials may alter the setting behavior of the silicone impression and bite registration materials.
- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- For dental use only!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Reuse: Standard surface disinfection as required. No specific maintenance required.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.
- When finishing the provisional, sanding dust can get into the eyes and respiratory tract.
- Wear a face mask and protective goggles.
- Protect against grinding dust by proper ventilation – do not inhale the dust.

Product performance characteristics

Compressive strength ≥ 250 MPa, flexural strength ≥ 60 MPa.

Storage/disposal

- Store in a dry place at 15 to 25 °C/59 to 77 °F! Do not store in the refrigerator!
- Protect against exposure to direct sunlight!
- After the application, leave the used mixing tip on the Automix Cartridge as a seal!
- Do not use after the expiration date!
- Disposal must comply with national regulations.

Composition

Dental glass, EBPADMA, aliph. polyester diurethanes, unsaturated polyester resin, SiO₂, additives.

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