# PermaCem Dual

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English

### Product description

PermaCem Dual is a universal, adhesive, dual-curing compomer luting cement for cementing indirect restorations in combination with an enamel-dentine bonding agent such as LuxaBond Total Etch (Total-Etch bonding system), or LuxaBond Universal and a metal or ceramic bonding agent, e.g. LuxaBond Universal.

### Intended purpose

Definitive cementation of crowns, bridges, inlays and onlays made of metal, ceramics or composites.

### Limitations of use

Do not use PermaCem Dual with single-bottle dentine bonding agents.

### Indication

Cementation of indirect restorations and prostheses for the treatment of large dental defects, tooth loss, malocclusions and malformations of the teeth.

### Contraindication

- Do not apply PermaCem Dual directly to the exposed dental pulp.
- Do not use the material where there are known allergies to any of the components or contact allergies.

### Patient target group

Persons receiving treatment as part of a dental procedure.

### Intended user

Dentist

### Notes for use

- With the Smartmix version the material is applied directly without a dispensing device.

- With the Automix version the material is applied using a 1:1 dispenser.
- Cementation of combined restorations such as with telescopic units: Using a suitable separating agent, ensure that PermaCem Dual cannot flow between primary and secondary parts.
- Light-curing units should emit at 450 nm and should be checked regularly. The light intensity should be a minimum of 600 mW/cm<sup>2</sup>. Place the light-curing unit as close as possible to the material. Observe the manufacturer's instructions for the light-curing unit used.

## Using the Automix Cartridge (for image, see »Handling«)

### Use of the Automix Cartridge

1. Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] all the way.
2. Lift up the plastic lever [C], insert the Automix Cartridge and lock it using the plastic lever [C].

**Note:** Make sure that the notches on the Automix Cartridge line up with those on the Automix Dispenser.

### Attaching a mixing tip

1. Remove the cartridge cap [D] or used mixing tip after turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [E].

**Note:** In order to achieve an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available in refill packs.

**Note:** Make sure that the notches on the mixing tip line up with those on the Automix Cartridge.

3. Twist the mixing tip clockwise by 90° to lock it.

### Applying the material

- ▶ When using an Automix Cartridge for the first time, squeeze out a pea-sized amount of material and discard it.
- ▶ Mix the material in the mixing tip using the lever [F] on the Automix Dispenser and apply directly.

**Note:** After the application, leave the used mixing tip on the Automix Cartridge as a seal!

### Removal of the Automix Cartridge

1. Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] all the way.
2. Lift up the plastic lever [C] and remove the Automix Cartridge.

## Using the Smartmix Syringe (for image, see »Handling«)

1. Before attaching the mixing tip, remove the cap [A] or used mixing tip after turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [B].

**Note:** In order to achieve an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available in refill packs.

**Note:** Make sure that the notches on the mixing tip line up with those on the Smartmix Syringe [C].

3. Twist the mixing tip clockwise by 90° to lock it.
4. The material is mixed in the mixing tip when ejected and can then be applied directly.

**Note:** When using a Smartmix Syringe for the first time, squeeze out a pea-sized amount of material and discard it!

**Note:** After application, leave the used mixing tip on the Smartmix Syringe as a seal!

## Disinfection of the Automix Dispenser

- The Automix Dispenser can be reused after disinfection, although it should be replaced if there are visible signs of wear. Disinfection can be carried out using a commercially available dipping disinfectant. The Automix Dispenser is also autoclavable.
- Dirty or contaminated mixing tips used as seals can be cleaned with disinfectants normally used in the practice.
- DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute).

## Recommended use

1. Prepare the clean and dry restoration as follows:
  - Etchable silicate ceramics: Etch the restoration with 5 % hydrofluoric acid in accordance with the instructions from the manufacturer of the etching gel and rinse thoroughly using a water spray. Subsequently dry with air that is free of oil and water, and prepare the surface using a silanizing agent (e.g. Vitique Silane or LuxaBond Universal) in accordance with the manufacturer's instructions.
  - Metals and metal alloys: Pretreat in accordance with the manufacturer's instructions. Unless described otherwise, sandblast the inside of the restoration with aluminum oxide, clean with alcohol, dry with oil-free and water-free air and, if necessary, apply a metal bonding agent according to the manufacturer's instructions.
  - Oxide ceramics: Pretreat in accordance with the manufacturer's instructions. Ensure that there is sufficient parallel-wall retention. Unless described otherwise, sandblast the inside of the restoration with aluminum oxide, clean with alcohol, dry with oil-free and water-free air and apply a suitable ceramic bonding agent (e.g. LuxaBond Universal) according to the manufacturer's instructions.
2. Thoroughly clean the prepared tooth of any remnants of provisional cement with water spray and dry with oil-free and water-free air. In doing so, ensure that the dentin retains a little residual moisture and does not become too dry.

**Note:** It is advisable to use a rubber dam.

3. Apply a suitable bonding agent to the prepared tooth according to manufacturer's instructions. DMG recommends LuxaBond Total Etch (Total-Etch-System) or LuxaBond Universal.
4. Apply PermaCem Dual in the pretreated restoration or on the prepared tooth.
5. Insert the restoration within 60 seconds after mixing has begun and attach to the dry, prepared tooth using slight pressure.
6. Instruct the patient to adopt habitual occlusion.
7. Remove cement residues within 60 to 90 seconds after inserting the restoration with the aid of a microbrush, brush, foam pellet or scaler. Remove interdental excess carefully using dental floss.
8. Cure the material with a suitable light-curing device (e.g. LED or halogen lamp) for at least 20 seconds. The restoration is then initially fixed. The restoration can support a full load after 360 seconds.

## Timing

Working time (including mixing time)*	≤ 60 s
Maximum intraoral curing (including mixing time)	360 s
Light exposure time for curing (optional) min. 20 s	20 seconds with a light intensity of ≥ 1000 W/cm <sup>2</sup> or 40 seconds with a light intensity of ≥ 600 mW/cm <sup>2</sup>
Removal of excess material (without light-curing)	60 - 90 s

\* The specified times are valid at a room temperature of 23 °C and a normal relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times.

## Interactions

- Residues of methacrylate-based materials may alter the setting behavior of silicone impression and bite registration materials.
- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.

## Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

## Warnings/precautions

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

## Product performance characteristics

Complies with ISO 4049:2019, Type 2, Class 3 (including shade stability and radiopacity of ≥ 190 % Al).

## Storage/disposal

- Store in a dry place at room temperature (2 - 25 °C/ 36 - 77 °F)!

- After application, leave the used mixing tip on the Automix Cartridge / Smartmix Syringe as a seal
- Do not use after the expiration date.

## Composition

Dental glass, bis-GMA, EBPADMA, TEDMA, DMAEMA, TMPTMA, MMHE, SiO<sub>2</sub>, additives. Inorganic filling materials: approx. 50 vol% (0.02 – 4 µm).

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