PermaCem Dual

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English

Product description

PermaCem Dual is a universal, adhesive, dual-curing compomer luting cement for cementing indirect restorations in combination with an enamel-dentin bonding agent such as LuxaBond Total Etch (total-etch system), or a metal or ceramic bonding agent, e.g. LuxaBond Universal.

Intended purpose

Definitive cementation of crowns, bridges, inlays and onlays made of metal, ceramics or composites.

Limitations of use

Do not use PermaCem Dual with single-bottle dentin bonding agents.

Indication

Cementation of indirect restorations and prostheses for the treatment of large dental defects, tooth loss, malocclusions and malformations of the teeth.

Contraindication

- Do not apply PermaCem Dual directly to the exposed dental pulp.
- Do not use the material if there are known allergies to any of the components or contact allergies.

Patient target group

People who are receiving dental treatment.

Intended user

Dentist

Notes for use

- If the material is applied directly in the patient's mouth using a single-use applicator, the single-use applicator must be used only on this one patient for hygienic reasons.
- With the Smartmix version, the material is applied directly without a dispensing unit.
- With the Automix version, the material is applied using a 1:1 dispenser.
- Cementation of combined restorations such as with telescopic restorations: Using a suitable separating agent, ensure that PermaCem Dual cannot flow between primary and secondary parts.
- Light-curing units should emit at 450 nm and should be checked regularly. The light intensity should be a minimum of 600 mW/cm². Place the light-curing unit as close as possible to the material. Observe the manufacturer's instructions for the light-curing unit used.
- Dirty or contaminated mixing tips used as caps can be cleaned with standard disinfectant solutions.

Handling the Automix Cartridge (Fig. see "Handling")

Inserting the Automix Cartridge

1. Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] all the way.
2. Lift up the plastic latch [C], insert the Automix Cartridge and lock it using the plastic latch [C].

Note: Make sure that the notches on the Automix Cartridge line up with those on the Automix Dispenser.

Attaching a mixing tip

1. Remove the cartridge cap [D] or used mixing tip after turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [E].

Note: To achieve an optimum mixing result, use the mixing tips available from DMG. All mixing tips are also available in refill packs.

Note: Make sure that the notches on the mixing tip line up with those on the Automix Cartridge.

3. Turn the mixing tip 90° clockwise to lock it.

Applying the material

- ▶ When using an Automix Cartridge for the first time, extrude a pea-sized amount of material and discard it.
- ▶ Mix the material in the mixing tip using the lever [F] on the Automix Dispenser and apply directly.

Note: After use, clean the used mixing tip with a standard disinfectant solution and leave it on the Automix Cartridge to act as a cap!

Removing the Automix Cartridge

1. Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] all the way.
2. Lift up the plastic latch [C] and remove the Automix Cartridge.

Handling the Smartmix Syringe (Fig., see "Handling")

1. Before attaching the mixing tip, remove the cap [A] or used mixing tip after turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [B].

Note: To achieve an optimum mixing result, use the mixing tips available from DMG. All mixing tips are also available in refill packs.

Note: Make sure that the notches on the mixing tip line up with those on the Smartmix Syringe [C].

3. Turn the mixing tip 90° clockwise to lock it.
4. The material is mixed in the mixing tip when ejected and can then be applied directly.

Note: When using a Smartmix Syringe for the first time, apply a pea-sized amount of material and discard it!

Note: After use, clean the used mixing tip with a standard disinfectant solution and leave it on the Smartmix Syringe to act as a cap!

Disinfecting the Automix Dispenser

- The Automix Dispenser can be reused after disinfection, although it should be replaced if there are visible signs of wear. Disinfection can be carried out using a commercially available immersion disinfectant.
- DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute).

Recommended use

1. Prepare the clean and dry restoration as follows:
 - Etchable silicate ceramics: Etch the restoration with 5 % hydrofluoric acid in accordance with the instructions from the manufacturer of the etching gel and rinse thoroughly using a water spray. Then dry with air that is free of oil and water, and prepare the surface using a silanizing agent (e.g. Vitique Silane or LuxaBond Universal) in accordance with the manufacturer's instructions.
 - Metals and metal alloys: Pretreat in accordance with the manufacturer's instructions. Unless described otherwise, sandblast the inside of the restoration with aluminum oxide, clean with alcohol, dry with oil-free and water-free air and, if necessary, apply a metal bonding agent according to the manufacturer's instructions.
 - Oxide ceramics: Pretreat in accordance with the manufacturer's instructions. Ensure that there is sufficient parallel-wall retention. Unless described otherwise, sandblast the inside of the restoration with aluminum oxide, clean with alcohol, dry with oil-free and water-free air and apply a suitable ceramic bonding agent (e.g. LuxaBond Universal) according to the manufacturer's instructions.
2. Thoroughly clean the prepared tooth of any remnants of provisional cement with water spray and dry with oil-free and water-free air. In doing so, ensure that the dentin retains a little residual moisture and does not become too dry.

Note: It is advisable to use a rubber dam.

3. Apply a suitable bonding agent to the prepared tooth according to manufacturer's instructions. DMG recommends LuxaBond Total Etch (total-etch system) or LuxaBond Universal.
4. Apply PermaCem Dual in the pretreated restoration or on the prepared tooth.
5. Insert the restoration within 60 s after mixing has begun and attach to the prepared tooth using slight pressure.
6. Instruct the patient to adopt their habitual occlusion.
7. Remove any excess cement within 60 to 90 s after inserting the restoration with the aid of a microbrush, brush, foam pellet or scaler. Remove interdental excess carefully using dental floss.
8. Cure the material with a suitable light-curing unit (e.g. LED or halogen lamp) for at least 20 s. The restoration is then initially fixed. The restoration can support a full load after 360 s.

Timing

Working time including mixing time*	≤ 60 s
Maximum intraoral curing (including mixing time)	≥ 360 s
Exposure time for curing (optional) min. 20 s	20 s with a light intensity of ≥ 1000 mW/cm ² or 40 s with a light intensity of ≥ 600 mW/cm ²
Removal of excess material (without light-curing)	60 – 90 s

* The specified times are valid at a room temperature of 23 °C (73 °F) and a normal relative humidity of 50 %. Higher temperatures shorten these times; lower temperatures lengthen these times.

Interactions

- Residues of methacrylate-based materials may alter the setting behavior of the silicone impression and bite registration materials.
- Materials containing eugenol, moisture, and oily air may inhibit polymerization at the contact site.

Residual risks/side effects

- Please check the product for damage and visible contamination before use.
- Hypersensitivity and/or allergic reactions cannot be ruled out in individuals who are predisposed.
- Incorrect use of the product may lead to excessive loss of dental hard tissues, as well as respiratory tract irritation, irritation and/or injury to the skin or mucous membranes, accidental ingestion and infections.
- A low residual risk of unknown side effects cannot be ruled out.

Warnings/precautions

- Observe the safety data sheet.
- For dental use only!
- Avoid contact with the skin! In the event of accidental contact with the skin, immediately wash the affected area thoroughly with soap and water.

- Avoid contact with the eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products used with the material/materials.
- The manufacturer (productsafety@mtcompanies.com) and the responsible registration authorities must be notified of any serious incidents involving this product.

Product performance characteristics

Complies with ISO 4049:2019, Type 2, Class 3 (including shade stability and radiopacity of $\geq 190\%$ Al).

Storage/disposal

- Disposal must comply with national regulations.
- Store at 2 – 25 °C (36 – 77 °F) in a dry place!
- After use, clean the used mixing tip with a standard disinfectant solution and leave it on the Automix Cartridge/Smartmix Syringe to act as a cap!
- Do not use after the expiration date!

Composition

Dental glass, bis-GMA, EBPADMA, TEDMA, DMAEMA, TMPTMA, MMHE, SiO₂, additives
Inorganic fillers: approx. 50 vol% (0.02 – 4 µm).

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