Silagum Mono

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Instructions for use English

Product description

Monophase impression material based on addition-curing silicones. The material is available in the Automix Cartridge.

Intended purpose

- Impressions for the fabrication of indirect restorations, such as crowns, bridges, inlays and onlays
- Functional impressions

Limitations of use

There are no known limitations of use.

Indication

All dental situations requiring final impressions, e.g. tooth loss, caries and trauma.

Contraindication

Do not use the material if there are known allergies to any of the components or contact allergies.

Patient target group

People who are receiving dental treatment.

Intended user

Dentist, dental assistant staff

Notes for use

- The recommended times can be found in the section headed “Timing.”
- If the material is applied directly in the patient’s mouth using a single-use applicator, the single-use applicator must be used only on this one patient for hygienic reasons.
- To achieve an optimum mixing result, use the mixing tips available from DMG. All mixing tips are available as refill packs.
- The Intraoral Tips are to be used for direct, intrasulcular application around the preparation.
- In order to improve the adhesion of the impression material to metal or plastic trays, a tray adhesive for A-silicones should be applied (e.g. DMG Tray Adhesive).
- Silagum Mono can be used in all common non-perforated impression trays.

Handling the Automix Cartridge (Fig. see “Handling”)

Inserting the Automix Cartridge

1. Push up the lever [A] at the back of the dispensing unit and pull back the slide [B] completely.
2. Lift up the plastic latch [C], insert the Automix Cartridge and lock it using the plastic latch [C].

Note: Make sure that the notches on the Automix Cartridge line up with those on the dispensing unit.

Attaching a mixing tip

1. Remove the cartridge cap [D] or used mixing tip after turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [E].

Note: Make sure that the notches on the mixing tip line up with those on the Automix Cartridge.

3. Turn the mixing tip 90° clockwise to lock it.
4. Attach an Intraoral Tip [F], if needed.

Applying the material

- When using an Automix Cartridge for the first time, squeeze out a pea-sized amount of material and discard it.
- Mix the material in the mixing tip using the lever [G] on the dispensing unit and apply directly.

Note: After use, clean the used mixing tip with a standard disinfectant solution and leave it on the Automix Cartridge to act as a cap!

Removing the Automix Cartridge

1. Push up the lever [A] at the back of the dispensing unit and pull back the slide [B] completely.
2. Lift up the plastic latch [C] and remove the cartridge.

Disinfecting the Automix Dispenser

- The Automix Dispenser can be reused after disinfection, although it should be replaced if there are visible signs of wear.
- Disinfection can be carried out using a commercially available immersion disinfectant. DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute).

Recommended use

1. Fill the prepared impression tray with Silagum Mono. When doing so, hold the impression tray at a slight angle beneath the mixing tip.

Note: To avoid air bubbles, the end of the mixing tip should remain immersed in the material throughout application.

2. Fill a small amount of material into a commercially available elastomer syringe or use the 50 ml cartridge directly and inject around the preparation using an Intraoral Tip.

Note: Filling of the impression tray and injection should be concluded at the same time, as the higher temperature in the mouth causes the injected material to set more quickly than the material in the impression tray.

3. Position the filled impression tray slowly and gently in the mouth within the working time and take the impression of the preparation.
4. Remove the impression from the mouth once the setting time has ended (see "Removal").

Removal

- In the area where the teeth have been prepared, remove the impression parallel to the longitudinal axis of the teeth.
- After removal, proceed as described in the section "Cleaning and disinfection of the final impression."

Cleaning and disinfection of the final impression

- Disinfect the final impression using immersion solutions only.
- DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute) or the EPA (Environmental Protection Agency). Observe the manufacturer's instructions!
- Ensure that the disinfectant solution used is indicated for use with addition-curing impression silicones.
- Leave the final impression to dry and store until casting (see "Storage of the impression").
- Further information is also available in the DGZMK's statement on the disinfection of dental impression materials: www.dgzmk.de

Casting and galvanization

- After removing the impression from the mouth, wait at least 30 min before casting. DMG recommends making the cast from hard modeling plaster (ISO 6873:2013; Type 3, e.g. Moldano, Bayer) and super-hard plaster for molds (ISO 6873:2013; Type 4, e.g. Fujirock II, GC).
- The impression can be copper or silver-plated.

Storing the impression

- Store the impression at 15 – 25 °C (59 – 77 °F) in a dry place!
- Protect from exposure to direct sunlight!

Timing

Silagum Mono	
Mixing time [s]	
Working time including mixing time* [s]	≤ 90
Setting time in the mouth [s]	≥ 150

* The specified times are valid at a room temperature of 23 °C (73 °F) and a normal relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times. The working time describes the period from when mixing begins until the latest point at which the tray must be positioned in the patient's mouth.

Interactions

- Rinse off any residual retraction solution thoroughly with water, otherwise this may lengthen the setting time in some instances.
- Never allow the impression to come into contact with organic solvents or liquids containing solvents as this may result in a swelling of the material and may therefore lead to inaccuracies in the cast.
- Do not allow the impression material to come into contact with latex gloves. In the event of prior contact with latex gloves, wash hands carefully. DMG recommends the use of PVC or PE-based gloves.
- Do not mix Silagum Mono with condensation curing silicone or polyether impression materials.
- Residues of methacrylate-based materials may alter the setting behavior of the silicone impression materials. If a restoration was prepared beforehand, carefully clean the preparation and the surrounding tissue with alcohol-soaked gauze and water spray before taking the impression.

Residual risks/side effects

- Please check the product for damage and visible contamination before use.
- Hypersensitivity and/or allergic reactions cannot be ruled out in individuals who are predisposed.
- A low residual risk of unknown side effects cannot be ruled out.

Warnings/precautions

- For dental use only!
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products used with the material/materials.
- Observe the safety data sheet.
- The manufacturer (productsafety@mtcompanies.com) and the responsible registration authorities must be notified of any serious incidents involving this product.

Product performance characteristics

Meets the requirements of ISO 4823:2025, Type 2: medium-bodied consistency

Storage/disposal

- Store at 15 – 25 °C (59 – 77 °F) in a dry place.
- Protect from exposure to direct sunlight!
- After use, clean the used mixing tip with a standard disinfectant solution and leave it on the Automix Cartridge to act as a cap!
- Do not use after the expiration date!
- Disposal must comply with national regulations.

Composition

SiO₂, vinyl silicones/siloxanes, methyl silicones/siloxanes, additives, pigments, catalysts, hydrogen silicone

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