VITIQUE

Silane

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English

Indication

Silanization of siliceous ceramics for improved adhesion of resin based materials, Improved bond between the luting cement and the glass fiber enforced composite post, LuxaPost

Recommended application

1. Prepare ceramic restoration with hydrofluoric acid. Carefully follow the directions for using hydrofluoric acid. Clean root post with solvents such as acetone or ethanol.
2. Mix 1 to 2 drops each of Vitique Silane Adhesive and Vitique Silane Activator in a ratio of 1:1 for approx. 15 seconds. Use silane mixture immediately.
3. Apply a thin layer of the silane mixture with a brush onto the prepared ceramic surface and let set for approx. 10 seconds.
4. Blow off the solvent carefully with a gentle stream of air.
5. Now the tooth can be cemented.

Please note

Do not apply the material to exposed dentine or enamel.

Avoid contact with skin and mucous membranes. In case of contact, wash immediately with water and soap. In case of contact with eyes, rinse immediately with a lot of water and contact a physician.

So far no side effects are known. In individual cases hypersensitivity to components of the material cannot be precluded.

Hold the bottle vertically to draw out material.

Composition

Silane; activator: distilled water, ethanol, phosphoric acid and catalyst

Storage

Store at room temperature not exceeding 25°C/77°F. Recap bottles after use. Do not use after expiration date.

Available product sizes

1 x 3 ml/2.7 g bottle of Vitique Silane Adhesive	REF 213118
1 x 3 ml/2.7 g bottle of Vitique Silane Activator	
1 brush holder + 50 brushes	

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

Keep away from children! For dental use only!

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