# DeltaFil Conditioner

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## Instructions for use English

### Product description

DeltaFil Conditioner is an aqueous polyacrylic acid solution for treating the smear layer after preparation and for conditioning the cavity, so as to optimize the bond between DeltaFil and the tooth structure.

### Intended purpose

- Smear layer treatment process
- Cavity conditioning

### Limitations of use

Use only in combination with DeltaFil.

### Indication

- Loss of dental enamel due to caries, trauma, tooth wear or developmental disorders
- Secondary prevention of caries

### Contraindication

- Do not use the material where there are known allergies to any of the components or contact allergies.
- Do not use directly on the exposed pulp.

### Patient target group

People to be treated in the course of a dental procedure.

### Intended user

Dentist

### Notes for use

- Observe the manufacturer's instructions for all materials that are used with the product.
- Close the bottle carefully after use and store it in a place protected from light.

- In deep cavities in areas close to pulp, provide protection in the form of a suitable liner, e.g. by using a small quantity of a strongly adhesive calcium hydroxide preparation.

## Recommended use

1. Pour DeltaFil Conditioner into a clean mixing well or onto a clean mixing pad.
2. Use a cotton pellet or a mini applicator sponge to absorb the DeltaFil Conditioner, apply it to the prepared area and allow to take effect for 10 s.
3. Then rinse thoroughly with water. Remove any excess water (e.g. with a cotton pellet or by blow-drying it lightly with a gentle flow of air).
4. The prepared surface should still look moist for the next steps. Do not dry the tooth structure too much!

## Interactions

DeltaFil Conditioner should not be applied to acid-soluble liners.

## Residual risks/side effects

- We place great importance on the quality of our products. However, please check the product for damage and visible contamination before use.
- A low residual risk of unknown side effects cannot be ruled out.
- Hypersensitivity and/or allergic reactions cannot be ruled out in individuals who are correspondingly predisposed.
- If DeltaFil Conditioner comes into contact with the mucous membranes, it may cause skin irritation. Rinse the affected skin area thoroughly with water.
- Improper use of the product can lead to excessive loss of tooth substance, eye irritation, respiratory tract irritation, skin or mucous membrane irritation and/or damage, choking and infection.

## Warnings/precautions

- For dental use only!
- Make sure the product is not swallowed! If it is accidentally swallowed, rinse the mouth thoroughly with water to eliminate any remaining material from the oral cavity. Do not induce vomiting!
- Avoid contact with the skin! In the event of accidental contact with the skin, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer (productsafety@mtcompanies.com) and to the responsible registration authorities.

## Product performance characteristics

Improves the bond between the tooth substance and DeltaFil.

## Storage/disposal

- Store at 15 – 25 °C (59 – 77 °F) in a dry and dark place!
- Do not use after the expiration date!
- Disposal must comply with national regulations.

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