
DeltaFil

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Instructions for use — English

Product description

DeltaFil is a classic, radiopaque glass ionomer cement in capsule form for restorative treatment. DeltaFil bonds chemically to enamel and dentin and releases fluoride ions. It is available in the shades White and Dark as well as A1, A2, A3, and A4 of the Vita tooth shade guide. The extrudable volume per capsule is at least 0.13 ml.

Intended purpose

- Class I and II permanent restorations: areas not subject to occlusal load; also for sandwich technique
- Class I, II*, and V permanent restorations; also for sandwich technique
- Class I, II, III, and V temporary restorations
- Extended fissure sealing
- Base/liner
- Class I, II, III, and V deciduous teeth restorations.

* Also for areas subject to occlusal load, provided the isthmus measures less than half of the intercuspal area

Limitations of use

Application is contraindicated if a relative dry isolation or the recommended application technique is not possible.

Indication

- Loss of dental hard tissue due to caries, fractures, tooth abrasion or developmental disorders
- Secondary prevention of caries

Contraindication

- Do not use directly on the exposed pulp.
- Do not use the material if there are known allergies to any of the components or contact allergies.

Patient target group

People who are receiving dental treatment.

Intended user

Dentist

Notes for use

- Rinse instruments in water immediately after use.
- Observe the manufacturer's instructions for all materials that are used with the product.
- In deep cavities in areas close to pulp, provide protection in the form of a suitable liner, e.g., apply a small quantity of a strongly adhesive calcium hydroxide preparation.
- Keep the times between activation, mixing, and application of the capsule contents short to avoid premature curing of the material (see "Timing").
- While applying, keep the tip of the capsule immersed in the material to avoid the formation of air bubbles.
- Do not use DeltaFil without DeltaFil Conditioner.
- For optimum results, use compatible dispensers as specified in the enclosed compatibility table (Short guide).
- Once the DeltaFil restoration has been placed and finished, the use of DeltaFil Coat is optional.

Handling the DeltaFil Capsule (Fig., see "Handling")

Activating and mixing the DeltaFil Capsule

1. Shake the powder in the capsule before inserting it into a suitable dispenser (as per the compatibility table in the enclosed Short guide). Do not activate by tapping on a table or any other surface! (Fig. 1)
2. Insert the capsule by fitting the groove into the dispenser holder (Fig. 2) and activate it by firmly (or quickly) squeezing the dispenser once. This should cause the retaining ring to release (Fig. 3).
3. Remove the capsule (Fig. 4), insert it into the mixer fork of a capsule mixer (e.g., CapMix, Silamat S6, or RotoMix), and mix at 4000 to 4500 rpm for 10 s (CapMix, Silamat S6) or for 8 s (RotoMix) (Fig. 5).

4. Reinsert the capsule into the dispenser and remove the cap (Fig. 6). The material is dispensed by pressing repeatedly and can be applied directly. Apply within 10 s (Fig. 7).
5. After pressing the ejection button, the capsule can be removed from the dispenser (Fig. 8).

Recommended use

1. Prepare the tooth as usual. When shaping the cavity it should be ensured that the material layer is not thinner than 0.5 mm at any point.
2. Then rinse the cavity thoroughly with water, remove excess water (e.g., with a cotton pellet or blow dry lightly with a gentle air flow). The prepared surface should still be slightly damp. The dentin must not become too dry.
3. Protect areas near the pulp.
4. Using DeltaFil Conditioner (observe the relevant instructions for use!): apply DeltaFil Conditioner with a cotton pellet or mini applicator sponge to the prepared area and leave in place for 10 s. Then rinse thoroughly with water. Remove excess water as described above. The surface should still look moist for the next steps.

Note: DeltaFil Conditioner is a polyacrylic acid solution in water, intended for treating the smear layer and conditioning the enamel and dentin to optimize the bond. Apply as described above. Carefully seal the bottle of the DeltaFil Conditioner after use.

5. If necessary, place a matrix.
6. Apply the mixed material directly from the capsule into the cavity, compress it, adapt it to the edges of the cavity, and model the surface. DeltaFil can be molded for up to 120 s from the start of mixing.
7. Immediately after modeling, care should be taken to prevent moisture contamination and excessive drying. If this cannot be guaranteed, apply DeltaFil Coat immediately and light-cure it. Follow the relevant instructions for use!
8. Then let the material cure for approx. 260 s.
9. The restoration can be finished as usual, e.g., with abrasive tooth grinders and finishing burs, 360 s after mixing has been started. Grinding can be performed with or without water spray.
10. DeltaFil Coat can be optionally applied and cured to protect the restoration. Follow the relevant instructions for use!
11. Disinfect the dispenser as indicated in the manufacturer's instructions for use.

Timing

Mixing time	10 s (CapMix/ Silamat S6) or 8 s (RotoMix)
Apply the capsule content after mixing	≤ 10 s
Working time (23 °C/73 °F) from the start of mixing*	≤ 120 s
Net setting time (37 °C/99 °F)	≥ 260 s
Begin the final polishing from the start of mixing	≥ 360 s

* The specified times are valid at a room temperature of 23 °C (73 °F) and a normal relative humidity of 50%. Higher temperatures shorten these times; lower temperatures lengthen these times.

Interactions

No interactions are currently known.

Residual risks/side effects

- Please check the product for damage and visible contamination before use.
- Hypersensitivity and/or allergic reactions cannot be ruled out in individuals who are predisposed.
- Improper use of the product can lead to excessive loss of dental hard tissue, eye irritation, respiratory tract irritation, skin or mucous membrane irritation and/or damage, choking, and infection.
- A low residual risk of unknown side effects cannot be ruled out.

Warnings/precautions

- For dental use only!
- Single use only! DeltaFil Capsule is intended for single use, and each capsule may only be used on one patient for hygienic reasons.
- Avoid contact with the skin! In the event of accidental contact with the skin, immediately wash the affected area thoroughly with soap and water.
- Wear safety goggles and avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products used with the material/materials.
- Observe the safety data sheet.
- The manufacturer (productsafety@mtcompanies.com) and the responsible registration authorities must be notified of any serious incidents involving this product.

Product performance characteristics

Complies with ISO 9917-1:2025 chemical type C and application type C

Storage/disposal

- Store at 15-25 °C (59-77 °F) in a dry place!
- Do not use after the expiration date!
- Disposal must comply with national regulations.

Composition

Ionomer glass, polyacrylic acid, water, tartaric acid, additives.

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