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Instructions for use — English

Product description

Flairesse Bleaching Gel CP is a gel containing 5 % or 10 % carbamide peroxide (equivalent to 1.8 % or 3.6 % hydrogen peroxide) for the gentle bleaching of teeth. After the dentist has given instructions, the product can be used independently by the patient, under the dentist's supervision.

Intended purpose

- Treatment of developmental enamel opacities in the permanent teeth of patients under the age of 18.
- Removal of medically-induced discoloration.

Limitations of use

- Whitening treatment is not effective against discolorations caused by metals (e.g. amalgam or silver points).
- The material is not suitable for heavy smokers.
- The material is not suitable for patients with a high alcohol intake.

Indication

- Developmental enamel opacities.
- Medically-induced discoloration.

Contraindication

- Do not use the material where there are known allergies to any of the components or contact allergies.
- Do not use the material with patients who are pregnant or breastfeeding.
- Do not use the material in cases of periodontal disease.
- Do not use the material in cases of exposed dentin, defective filling margins or untreated caries.

Patient target group

People to be treated in the course of a dental procedure.

Intended user

Initial treatment is carried out by the dentist at the dental practice. The patient is then able to use the bleaching product independently under the dentist's supervision.

Notes for use

- Demineralized areas (white spots) may become more visible at the start of treatment. Any difference in shade compared to the healthy tooth enamel will even out during the course of treatment.
- In cases of developmental enamel opacities, treatment may take longer.

Using the syringe

- The gel comes in a 5 ml syringe.
- The syringe must be resealed with its tip cap after each application.
- Only use it in combination with a bleaching tray.

Recommended use

Pretreatment

1. Thorough examination of the patient's mouth to identify any defective (provisional) restorations, untreated caries or exposed dentin. These areas must be excluded from the bleaching treatment or treated accordingly beforehand, as the patient may otherwise experience hypersensitivity.
2. Removal of external stains and tartar by having the teeth professionally cleaned.

Note: If irritation of the gingiva occurs, treatment must be delayed until the irritation subsides.

3. Assessment of the type and extent of the discoloration involved and whether or not tooth bleaching is indicated.
4. Assessment of the initial tooth shade using a suitable shade guide (e.g. Vita classical®). Photo documentation of the initial tooth shade can help patients to visualize their treatment success.

Note: Patients should be informed that the tooth shade of existing restorations will not change and that a shade difference to the surrounding enamel may occur.

5. If the patient is unable to follow the treatment instructions competently, the bleaching treatment should not be recommended.

Bleaching treatment

1. To manufacture the bleaching tray, suitable materials must be used in accordance with the manufacturer's instructions.
2. Open the syringe by removing the cap. Keep the cap in a safe place.

Note: The syringe must be resealed with its cap after use.

3. Dispense a pinhead-sized amount into each reservoir of the tray for the areas to be treated.

Note: Check whether the amount is correct by looking at how the material is distributed in the reservoir after the tray has been inserted. It may be necessary to adjust the amount for the next application.

4. Insert the filled bleaching tray. Remove any excess material with a soft cloth or a soft toothbrush.
5. Wear the tray for up to 8 hours, depending on the dentist's specific recommendation. Repeat the treatment daily until the desired tooth shade is achieved. Do not eat or drink while wearing the filled bleaching tray.

Note: If irritation of the gingiva occurs, treatment must be delayed until the irritation subsides.

Note: It is possible to pause the treatment for several days without jeopardizing its subsequent treatment.

Note: It is important to avoid staining foods and drinks (e.g. wine, coffee, cola) as well as acidic foods (e.g. citrus fruits) and drinks (e.g. lemonade) throughout the entire treatment cycle.

6. Remove the tray and rinse the mouth several times with water. Do not swallow the water! Any gel remaining on the teeth can be carefully removed with a soft toothbrush.
7. Clean the tray with cold water and a soft toothbrush. Store the tray in tray box provided.

Interactions

Carbamide peroxide can have a negative effect on the radical polymerization of plastics. The bleaching treatment must therefore not be carried out in the two weeks prior to or immediately after an adhesive restoration.

Residual risks/side effects

- We place great importance on the quality of our products. However, please check the product for damage and visible contamination before use.

- During use, increased tooth sensitivity may occur temporarily. Contact with oral mucosa may cause irritation. These reactions are transient and subside after the treatment is completed.
- Tooth sensitivity may occur during use.
- Hypersensitivity and/or allergic reactions cannot be ruled out in individuals who are correspondingly predisposed.

Warnings/precautions

- Keep out of the reach of children!
- Children and adolescents should only use the product under adult supervision.
- Avoid contact with the skin! In the event of accidental contact with the skin, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with mucous membranes! In the event of accidental contact with the skin, immediately wash the affected area thoroughly with water.
- Avoid contact with eyes! Risk of loss of sight! In case of accidental contact with eyes, immediately and carefully rinse with water for several minutes. If contact lenses are worn, remove these if possible. Continue rinsing.
- If exposure has occurred or is suspected, seek medical advice/attention.
- Contains alkali.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Patient should be told to avoid swallowing large amounts of gel.
- If there is a risk that the bleaching tray could be swallowed or aspirated, avoid the treatment.
- Patient should be advised to make an appointment if any symptoms occur.
- Serious incidents involving this product must be reported to the manufacturer (productsafety@mtcompanies.com) and to the responsible registration authorities.

Product performance characteristics

Complies with ISO 28399:2021.

Storage/disposal

Store in a dry place at 2 to 25 °C (36 to 77 °F).

Composition

Flairesse Bleaching Gel CP 5%: Water, glycerin, sodium hydroxide, carbamide peroxide, polyacrylic acid (cross-linked), xylitol, potassium nitrate, additives.
Flairesse Bleaching Gel CP 10%: Water, glycerin, sodium hydroxide, carbamide peroxide, polyacrylic acid (cross-linked), xylitol, potassium nitrate, additives.

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