



PRODUKTÜBERSICHT

Retraction Paste

Packaging: 25 Compules @ 0,37 g

REF 919381



New!

Retraction Paste

Your uncomplicated impression aide.

DMG
Chemisch-Pharmazeutische Fabrik GmbH
Elbgaustraße 248 22547 Hamburg Germany
Fon: +49. (0) 40.84006-0 Fax: +49. (0) 40.84006-222
info@dmg-dental.com www.dmg-dental.com
www.facebook.com/dmgdental





Your uncomplicated
impression aide.
DMG Retraction Paste

Every accurate restoration relies on a precise impression. DMG Retraction Paste is particularly helpful here:

Temporary sulcus widening has never been this easy or this effective – and it includes an astringent effect on the gingiva.

The newly developed compule made of special resin can be bent just as desired. You can be shaped to your preferred bend and work well everywhere in the patient's mouth.

The homogeneous consistency of DMG Retraction Paste is optimally adjusted. Thanks to its high-level stability, the paste remains in the spot where it was applied and once there, it takes effect optimally.

Indications

Temporary marginal gingival retraction and establishment of a dry sulcus environment

GOOD REASONS:

- ↗ Flexible cannula for simplest application
- ↗ Use of standard composite dispensers
- ↗ Remains secure in the sulcus thanks to optimum consistency
- ↗ Made from a base of natural alumina
- ↗ Good color contrast



Compule tip inserted into the sulcus.



Apply the retraction paste.



Leave the paste to work into the sulcus, to enable the gingival retraction.



The tissue is retracted mechanically.