

# Dental solutions for practices and laboratories



Icon Constic Luxatemp Star Silagum O-Bite LuxaPost **TEcosite T**LuxaCore Z LuxaBond PermaCem Luxatemp MaxProtect Vitique **¬** Ionosit Baseliner **¬** LuxaPrint MiniDam Tlairesse 

**DentaMile** 







# There's no doubt about it: dentistry and dental technology are fascinating and diverse fields to work in – and extremely demanding ones

As a dentist, your commitment goes far beyond the treatment chair. Whether it is as someone for your patients and staff to turn to, as a listening ear or as a manager and business person. In the dental laboratory, it is important to remain faithful to traditional craftsmanship while incorporating the use of cutting-edge technology. Constantly having to tackle new challenges is a fundamental part of being a dentist or a dental technician.

At DMG, we give our all to help lighten your load.

And we have been doing so passionately for over 60 years.

## A smile ahead together

## Welcome to the DMG community

#### Knowledge for all - the DMG training portfolio

At DMG, we support you with a wide range of training courses as well as high-quality products. Whether it is talks, workshops or live treatment - with us, you can learn online, in our training centre or on your own premises.

#### On-site for you

Throughout the year, we have a presence at numerous trade fairs and congresses worldwide. We look forward to your visit and to having the chance to share ideas with you.

#### You can find dates and information here:

#### **DMG Academy**

Visit our digital platform, the DMG Academy, to complete training on various topics at any time. The courses are free of charge and are available to dentists, dental assistants and dental technicians.















#### Interdentale

Interdentale is the event series made especially for young dentists. With innovative lectures and workshops - online and on-site!



my-interdentale.com/en







#### **DMG** newsletter

Register today for special offers, information on events, new products and much more!





#### Social media

Follow us on Facebook, Instagram or LinkedIn under 'dmgdental' to be kept up to date at all times. Our YouTube channel offers a large number of helpful videos.









## Digital workflow

# Contents

# The DMG products at a glance

#### Software DentaMile connect 13 DentaMile Lab5 18 DentaMile Lab5 Pro 18 Hardware DentaMile Desk MC-5 14 Resins DentaMile Wash MC 14 LuxaPrint Ortho Comfort 19 DentaMile Cure MC 15 LuxaPrint Ortho Flex 19 DentaMile Desk RS-1 15 LuxaPrint Ortho Plus 20 DentaMile Desk RS-5 16 LuxaPrint Ortho 20 DentaMile Wash RS 16 LuxaPrint Model 21 DentaMile Cure RS 17 LuxaPrint Tray 21 DentaMile Cure RS-XL 17

## Prevention and early intervention

Bleaching		Prophylaxis		Infiltration	
Flairesse		Flairesse Paste	24	Icon Proximal	26
Bleaching Gel CP 5%	23	Flairesse Gel/Foam	24	Icon Vestibular	26
Bleaching Gel CP 10%	23	Flairesse Varnish	25	DMG MiniDam	27
Bleaching Gel CP 16%	23				

## Direct composite fillings

Composite		Underfilling material		Core build-up material	
Ecosite Elements	30	Ionosit Baseliner	35	LuxaCore Z Dual	38
Ecosite One	31			LuxaCore Z Post System	38
Ecosite Bulk Fill	31	Bonding agent			
Constic	32	LuxaBond Universal	36	Root posts	
EcuSphere Flow	32	LuxaBond Total Etch	36	LuxaPost	39
EcuSphere Carat	33	Ecosite Bond	37		
EcuSphere Shape	33	Contax	37	Accessories	
				Safetip Dispenser	40
Glass ionomer cement				Safetip Dispenser Pro	40
DeltaFil	34			DMG Etching Gel	41
DeltaFil Coat	34				
DeltaFil Conditioner	35				

### Flairesse Bleaching Gel

A game changer for developmental enamel opacities: An important building block, especially in the treatment of young patients. Find out more on page 23



## Impression

Impression material				Retraction medium	
Honigum Pro Light	43	Silagum Light	47	DMG Retraction Paste	50
Honigum Pro Mono	43	Silagum Medium	47		
Honigum Pro Heavy	44	Silagum Mono	48	Bite registration materia	al
Honigum Pro Putty	44	Silagum Putty	48	O-Bite	51
Honigum Light	45	Silagum Putty Soft	49	O-Bite Scan	51
Honigum Mono	45	StatusBlue	49	LuxaBite	52
Honigum Heavy	46			DMG Tray Adhesive	52
Honigum Putty	46	Mixing unit		Silagum Comfort	53
		MixStar eMotion	50		
				Relining material	
				LuxaPick-up	53

## Temporary prosthetics

High-performance temporary solutions for long-term stability		Surface sealer Luxatemp Glaze & Bond	57	Luxatemp
Luxatemp MaxProtect	55	Temporary cementation		MaxProtect
		TempoCem ID	58	The perfect symbiosis of
Temporaries		TempoCem NE	58	durability and aesthetics.
Luxatemp Star	55	TempoCem	59	Find out more on page 55
Luxatemp Fluorescence	56			
Luxatemp Plus	56	Max	rotect	
Temporary repair		TV - WI		

## Permanent prosthetics

57

LuxaFlow Star

#### **Luting cements** PermaCem Universal 62 PermaCem 2.0 62 PermaCem Dual 63 PermaCem 63 Vitique 64 Vitique Silane 65



## Tables and overviews

## When to use which material


Direct restoration: Restorative materials	68	Prosthetic restoration:	
Prevention and infiltration	69	temporary and long-term stable	76
Impression	70	Automix Dispenser	77
Direct restoration: Bonding agent	71	Application Tips	78
Impression with Honigum	72	MixStar Tips	79
Impression with Honigum Pro	73	Symbols, Index	80
Cements: temporary and permanent	74		



## Flexible, fast and well connected

With DMG DentaMile, you decide for yourself how to use the digital workflow – even without your own 3D printer. Manufacture precise bite splints, dental models or bleaching trays with 3D printing? Fast, easy and economically efficient?

With the DentaMile connect software and the validated workflows, it couldn't be simpler. You can design and print the objects yourself or delegate them to a partner with whom you are networked in the cloud or to DMG's Design & Print Service. The DentaMile network opens up new and economically convincing possibilities.





## **3D** from a single source

An integrated approach: the DMG DentaMile portfolio offers you all the components for modern 3D dental printing from a single source. From the cloud-based software DentaMile connect, via high-performance 3D printers and post-processing units, to high-quality materials for a wide variety of applications.

# Scanned, designed, printed, delivered: **DentaMile Design & Print Service**



With our DentaMile Design & Print Service, you are given tailored solutions to meet your needs. You can use our DentaMile Print Service and have us print out your designs in 3D. Or use the Design Service as well, and we will also design the objects for you.

3D the way I want it.



Would you like more information on 3D printing? www.dentamile.com/en

Don't have a printer of your own?
No problem!
Don't have time to design yourse

Don't have time to design yourself?

No problem!



# The cloud software that overcomes technical barriers

DentaMile connect enables easy networking of all parties involved in the production process: people, instruments and technologies. Without spatial or technical barriers. Simple and efficient.

The cloud-based software developed by DMG simplifies 3D printing so that the technology can easily be completely integrated into everyday work.

DentaMile connect also ensures that the workflow in medical device manufacturing is secure and can be validated. Errors are minimised and new applications are made accessible without requiring any additional technical effort.

## 3D the way I want it: the flexible DentaMile usage options

Printing your own 3D models? Networking with a dentist or a partner laboratory? Just scanning? Or designing digital models as well? No matter how, you can participate in the digital workflow – even without your own 3D printer. The DentaMile connect software's flexible usage model offers you an outstanding range of possibilities, with four different usage options. You decide what you want to do and which option is right for your needs.

- Cloud-based: physical distance is not an issue
- »Smart Connectivity« for networking entire systems
- Control of all processing and fabrication procedures with one software: design, slicing, printing and post-processing
- Intuitive, easy to use
- → Flexible: four usage options, depending on your individual needs



# Validated and efficient: DentaMile workflows for bite splints and bleaching trays

Manufacture precision-fit bite splints or bleaching trays with 3D printing? Fast and efficient? Of course! DentaMile connect and the validated DentaMile workflows help you optimise workflows and achieve even more precise results. The holistic workflows take into account all process steps and components involved. Handling is made as user-friendly as possible. The validated processes guarantee the reliable quality of the printed items.





## **DentaMile workflows:**

concrete advantages for dental practices, laboratories and patients



#### High wearing comfort for the patient

Because of the thin material and optimal fit, the bite splint offers a very high level of wearing comfort



#### **Smarter working**

To the printed item in just a few steps, with simple, intuitive handling. Precise design becomes easier than ever.



#### Highly aesthetic

Barely visible, colourless transparency of the bite splints



#### DentaMile connect



The cloud-based software DentaMile connect facilitates the simple networking of all parties involved in the production process, with no spatial or technical barriers

#### **Product benefits**

- Encryption of patient data with the latest cloud technology from Microsoft Azure
- Consistent logging of workflows
- Software takes into account interaction between materials and manufacturing method.
- Users can get started immediately no major training necessary

- Time savings thanks to helpful default settings and automation
- Control of all processing and fabrication procedures with one software: design, slicing, printing and post-processing
- Consistent, MDR-compliant logging of workflow

#### Overview of options with DentaMile connect

Features included	DMC-ready *	Creator	Creator Plus	Designer	Producer
Creating patient data and uploading STL files	<b>✓</b> **	<b>✓</b>	<b>✓</b>	<b>✓</b>	<b>~</b>
Design workflow (including virtual articulators)	×	×	×	✓	✓
Automated 3D printing	×	×	×	×	<b>✓</b>
Sharing and delegating cases	<b>√</b> **	<b>✓</b>	<b>✓</b>	<b>✓</b>	<b>✓</b>
Connections to other DentaMile connect users	unlimited **	1 (Designer / Producer)	unlimited	unlimited	unlimited
Cloud storage	5 GB	5 GB	5 GB	5 GB	5 GB
Free updates	✓	✓	✓	✓	<b>✓</b>
Service and support	✓	✓	<b>✓</b>	✓	<b>✓</b>

<sup>\*</sup> This usage option is free of charge. Only in conjunction with a printer purchase

#### Packaging

DentaMile connect – Creator DentaMile connect – Creator Plus REF 170904 REF 170905 DentaMile connect – Designer DentaMile connect – Producer REF 170906 REF 170907

<sup>\*\*</sup> Data exchange only

#### DentaMile Desk MC-5



#### Compact 3D printer for easy digital workflow entry

#### Technical data

External dimensions:  $300 \times 400 \times 422.6$  mm; Weight: 16 kg; Build area (W/D/H):  $130 \times 73.1 \times 90$  mm; Technology DLP @385 nm; Native pixel:  $68 \mu m$ ; connections: 2 USB ports, 1 LAN port, 1 mains connection outside; Additional information: 2K light engine, 5" colour touch display, WLAN and LAN, High-precision image alignment and orientation correction process

#### **Product benefits**

Benefits of the 3D printer

- Compact size
- Suitable for a wide spectrum of dental 3D applications
- High precision
- Outstanding price-performance ratio

#### Benefits in use

- · User-friendly and simple user guidance
- Cover open on three sides for a larger working space
- Intelligent networking with DentaMile connect
- Compatible with all standard dental CAD programs (STL files)

#### Benefits in safety

- Reproducible results thanks to a software-supported, validated application
- Reliable service structure
- Developed, quality-tested and calibrated in Germany

#### Packaging

1 DentaMile Desk MC-5

REF 170500

#### DentaMile Wash MC



Simple cleaning unit for the systematic and standardised cleaning of printed objects according to manufacturer specifications

#### Technical data

External dimensions (W/D/H):  $320 \times 330 \times 246$  mm; Volume:  $140 \times 140 \times 120$  mm; Power supply unit: 12 V, 2 A external adapter; Additional information: 2.4" TFT touchscreen

#### **Product benefits**

- Designed for numerous wash cycles
- Removable basket
- Magnetic mixer
- Fast and fully automatic cleaning of the 3D-printed objects
- Preset cleaning programmes
- Colour touch display with simple user guidance
- Guaranteed efficient use of the cleaning fluids

#### Packaging

1 DentaMile Wash MC

#### **DentaMile Cure MC**



## Simple post-curing unit for final curing with validated parameters for a reliable end product

#### Technical data

External dimensions (W/D/H): 320 x 330 x 246 mm; Volume (W/D/H): 145 x 145 x 105 mm; Heating: Eco-passive up to 60 °C; Interface: 2.4" TFT touchscreen connection: 2.4 GHz WiFi; further information: Colour touch display with preset light-curing programs; LED technology (365 nm and 405 nm)

#### **Product benefits**

- UV illumination with exclusive omni 360° technology for uniform curing
- Triple-frequency technology with 3 wavelengths (365 nm, 385 nm, 405 nm)
- Suitable for various light-curing synthetic resins
- Fully automatic material curing programs
- Protection from UV radiation due to opaque walls and safety OFF system
- Network connectivity

#### Packaging

1 DentaMile Cure MC

REF 170502

#### DentaMile Desk RS-1



## The elegant and compact solution for dental 3D manufacturing

#### Technical data

External dimensions (W/D/H):  $320 \times 375 \times 300$  mm (with closed lid),  $320 \times 375 \times 391$  mm (in print mode); Build space:  $133 \times 75 \times 100$  mm; Wavelength: 385 nm; Max. workpiece height: 90 mm; resolution Full HD: 69 µm pixel size; Z resolution: 25-150 µm; Connectivity: WLAN, LAN, USB; Weight approx. 9 kg

#### **Product benefits**

- Force Feedback technology: optimises the printing process by actively managing platform speed for reliable prints
- Certified auto-calibration: Ensures consistent light intensity and print quality
- RFID material identification: Monitors the compatibility of the material and its expiry date
- Compact design: For flexible placements

#### Packaging

1 DentaMile Desk RS-1

#### DentaMile Desk RS-5



High-precision, high-speed 3D printer for a wide range of dental 3D printing applications

#### Technical data

External dimensions (W/D/H):  $335 \times 349 \times 541$  mm; Build area (W/D/H):  $130 \times 75 \times 115$  mm; Max. workpiece height: 115 mm; Native pixel +/- 34 µm; Connections: WLAN, TCP/IP, USB; Weight: approx. 33 kg

#### **Product benefits**

- Suitable for a wide spectrum of dental 3D applications
- Force Feedback option
- Large build area for the simultaneous manufacturing of several workpieces
- · Exceptionally high processing speed
- Precision without compromise
- Simple, material-saving handling and user guidance
- Validated process for reproducible printing results
- Safety via material identification (RFID)
- Intelligent networked workflow in combination with DentaMile Wash RS and DentaMile Cure RS
- Compatible with all standard dental CAD programs (STL files)
- Outstanding price-performance ratio
- Reliable service structure

#### Packaging

1 DentaMile Desk RS-5REF 1704501 DentaMile Desk RS-5 with RS-TurboREF 170451

#### DentaMile Wash RS



Cleaning unit for the systematic and standardised cleaning of printed objects according to the manufacturer's instructions as part of the supported process reliability

#### Technical data

External dimensions (W/D/H): 230 x 450 x 270 mm; Volume: 130 x 75 x 60 mm; Processing time between 4 and 8 minutes, depending on material; cleaning material: isopropyl alcohol; others according to material specifications; connections: WLAN/LAN

#### **Product benefits**

- Fully automatic cleaning in 4-8 minutes
- No contact with sticky resins and cleaning fluids
- Validated workflows
- Touch-display control
- Computer-controlled programs ensure the most effective use of cleaning fluids

#### Packaging

1 DentaMile Wash RS cleaning unit REF 170316

#### Accessories

3 DentaMile Wash RS cleaner swap box REF 170345 6 Empty bottles @ 800 ml for LuxaPrint Clean IPA REF 170343

#### **DentaMile Cure RS**



## Light curing unit for final curing as part of a monitored, controlled process

#### Technical data

External dimensions (W/D/H):  $230 \times 380 \times 270$  mm; Volume:  $130 \times 75 \times 60$  mm; Processing time depending on material; Connections: WLAN/LAN, USB

#### **Product benefits**

- Fully automatic material curing programs
- Integrated touchscreen
- Network connectivity
- 360° lighting
- Can be stacked with DentaMile Wash RS

#### Packaging

1 DentaMile Cure RS post-processing unit

REF 170317

#### **DentaMile Cure RS-XL**



## The high-performance post-curing unit for the series production

#### Technical data

External dimensions (W/D/H):  $650 \times 800 \times 1050$  mm; Volume:  $470 \times 375 \times 280$  mm; Processing time: depending on material; Connections: WLAN/LAN, USB; Weight approx. 150 kg

#### **Product benefits**

- Large volume: Enables the simultaneous post-curing of multiple printed components (up to 36 models at the same time)
- 360° curing: Ensures uniform post-curing of all surfaces of the printed object
- Modern vacuum technology: Shortens post-curing times while maintaining high quality
- Gas free: The use of vacuum means that no gas is required during operation, which significantly reduces operating costs
- Fully automatic material curing programs
- Integrated touchscreen
- Network connectivity
- Easy calibration with ACCS sensor to ensure functionality
- 360° UV lighting for uniform curing

#### Packaging

DentaMile Cure RS-XL 230V
 DentaMile Cure RS-XL 120V

REF 170420 REF 171420

#### **DentaMile Lab5**



## The professional printer for high volumes and efficiency

#### Technical data

External dimensions (W/D/H):  $600 \times 570 \times 1660$  mm; building platform (W/D):  $231 \times 130$  mm; Max. work-piece height: 300 mm; Native pixel +/- 30 µm; Connections: WLAN, TCP/IP, USB; Weight|: 110 kg; Resolution: 4K ( $3840 \times 2160$  pxl) with 385 nm LED

#### **Product benefits**

Superior technology

- Professional high-precision 4K light source
- Print chamber can be heated up to 35 °C
- Integrated force feedback system
- LED status bar
- Motorised hood can be opened and closed via foot sensor

#### Efficient application

- Very large build area for the simultaneous production of several workpieces
- Exceptionally high processing speed
- Precision without compromise

#### Maximum safety

- Software-supported, validated application creates reproducible dental device results
- Safety via material identification (RFID)
- Compatible with all standard dental CAD programs (STL files)

#### Packaging

1 DentaMile Lab5 Printer

REF 170352

#### DentaMile Lab5 Pro



The automated professional printer for high quantities and efficiency with an automatic refill and separation function (ASM) of the finished printed items

#### Technical data

External dimensions (W/D/H):  $600 \times 570 \times 1660$  mm; building platform (W/D):  $231 \times 130$  mm; Max. workpiece height: 300 mm or 100 mm (when using ASM); Native pixel +/- 30 µm; Connections: WLAN, TCP/ IP, USB; 120 kg; Resolution: 4K ( $3840 \times 2160$  pxl) with 385 nm LED

#### Product benefits

The DentaMile Lab5 Pro offers superior technology, efficient application and maximum safety (for details, see DentaMile Lab5)

#### Speciality

Automatic refill and disconnect function (ASM)

#### Packaging

1 DentaMile Lab5 Pro Printer with ASM

#### **LuxaPrint Ortho Comfort**



#### Light-curing 3D print resin

Indication: Additive manufacture of occlusal splints. Classified as a Class IIa medical device

Shade: Transparent

#### **Product benefits**

- Extremely resilient, yet flexible
- Very comfortable to wear
- Material can be printed thinly
- Increased material flexibility in the mouth
- Easy to use even at room temperature
- Tasteless and odourless
- Resistant to discolouration



#### Pack

385 nm

1	Bottle @ 5000 g, shade Transparent	REF 170277
1	Bottle @ 2000 g, shade Transparent	REF 170275
1	Bottle @ 1000 g, shade Transparent	REF 170255
1	Bottle @ 500 g, shade Transparent	REF 170253
1	Bottle @ 200 g, shade Transparent	REF 170251

## **LuxaPrint Ortho Flex**



#### Light-curing 3D print resin

Indications: Additive fabrication of flexible trays in the digital workflow. Classified as a Class I medical device

Shade: Transparent

#### **Product benefits**

- Versatile
- Optimal elasticity
- Highly tear-resistant and unbreakable
- Near-natural transparency
- Easy to clean



#### Pack

385 nm

1 Bottle @ 5000 g, shade Transparent	REF 170284
1 Bottle @ 2000 g, shade Transparent	REF 170281
1 Bottle @ 1000 g, shade Transparent	REF 170252
1 Bottle @ 500 g, shade Transparent	REF 170250
1 Bottle @ 200 g, shade Transparent	REF 170273

## LuxaPrint Ortho Plus



#### Light-curing 3D print resin

Indication: manufacture of splints and orthodontic appliances. Classified as a Class IIa medical device

Shade: Transparent

#### **Product benefits**

- Clear transparent
- Exact fit
- Not brittle



#### Pack

385 nm

1 Bottle @ 5000 g, shade Transparent	REF 170216
1 Bottle @ 2000 g, shade Transparent	REF 170215
1 Bottle @ 1000 g, shade Transparent	REF 170211
1 Bottle @ 500 g, shade Transparent	REF 170210
1 Bottle @ 200 g, shade Transparent	REF 170259
.05 nm	

#### 4

105	nm	
1	Bottle @ 1000 g, shade Transparent	REF 170213
1	Bottle @ 500 g, shade Transparent	REF 170212
1	Bottle @ 200 g, shade Transparent	REF 170260

#### LuxaPrint Ortho



#### Light-curing 3D print resin

Indications: Manufacture of individual surgical templates with maximum precision and orthodontic appliances. Classified as a Class I medical device

Shade: Transparent

#### **Product benefits**

- Exact fit of the drilling sleeve
- Drilling template in accurate position
- Clear transparent, for perfect control
- »Unobstructed View« of the surgical field
- Autoclavable, for high hygiene standards



#### Pack

#### 385 nm

1 Bottle @ 5000 g, sha	de Transparent	REF 17 0218
1 Bottle @ 2000 g, sha	de Transparent	REF 170217
1 Bottle @ 1000 g, sha	de Transparent	REF 170207
1 Bottle @ 500 g, shade	e Transparent	REF 170206
1 Bottle @ 200 g, shade	e Transparent	REF 170256
05 nm		
1 Pottle @ 1000 a shor	da Tuananauant	DEE 1 70000

#### 1 Bottle @ 1000 g, shade Transparent REF 170209 1 Bottle @ 500 g, shade Transparent REF 170208 1 Bottle @ 200 g, shade Transparent REF 170257

## **LuxaPrint Model**



#### Light-curing 3D print resin

Indication: Manufacturing of dental models using 3D printing. The printed model is suitable for producing thermoformed splints

Colours: Beige, Grey, Ivory, Transparent

#### **Product benefits**

- Finest surface structure for perfect identification of details
- Dimensional stability
- Precise



#### Pack

 $385\,\text{nm}\,/\,405\,\text{nm}$ 

1 Bottle @ 2000 g, shade Beige	REF 170265
1 Bottle @ 1000 g, shade Beige	REF 170201
1 Bottle @ 500 g, shade Beige	REF 170200
1 Bottle @ 200 g, shade Beige	REF 170264
1 Bottle @ 2000 g, shade Grey	REF 170263
1 Bottle @ 1000 g, shade Grey	REF 170293
1 Bottle @ 500 g, shade Grey	REF 170292
1 Bottle @ 200 g, shade Grey	REF 170294
1 Bottle @ 2000 g, shade ivory	REF 170279
1 Bottle @ 1000 g, shade ivory	REF 170276
1 Bottle @ 500 g, shade ivory	REF 170274
1 Bottle @ 200 g, shade ivory	REF 170272
385 nm / 405 nm	
1 Bottle @ 5000 g, shade Transparent	REF 170271
1 Bottle @ 2000 g, shade Transparent	REF 170270
1 Bottle @ 1000 g, shade Transparent	REF 170269
=-	

1 Bottle @ 500 g, shade Transparent

1 Bottle @ 200 g, shade Transparent

## **LuxaPrint Tray**



#### Light-curing 3D print resin

Indication: manufacture of individual impression trays. Classified as a Class I medical device

Shade: Turquoise

#### **Product benefits**

- Quick to print
- Extremely stable
- Exact fit



#### Pack

385 nm / 405 nm

1	Bottle @ 5000 g, shade Turquoise	REF 170261
1	Bottle @ 2000 g, shade Turquoise	REF 170258
1	Bottle @ 1000 g, shade Turquoise	REF 170205
1	Bottle @ 500 g, shade Turquoise	REF 170204
1	Bottle @ 200 g, shade Turquoise	REF 170266

REF 170268



# Flairesse Bleaching Gel CP 5% and 10%



# Flairesse Bleaching Gel CP 16%



## A game changer for the treatment of developmental enamel opacities

Indications: Treatment of developmental enamel defects in children, whitening of intrinsic and extrinsic discolouration due to various causes

#### **Product benefits**

- At-home bleaching for use with a custom tray
- Approved medical product
- Suitable for young patients
- Supportive treatment of developmental enamel opacities
- Conditioning of the protein-rich lesions for infiltration treatment
- Especially mild formulation with 5 % and 10 % carbamide peroxide (CP)
- Gentle whitening of opacities with yellowish or brownish opacities
- Bonds to the tray effectively and does not run
- · Contains xylitol and potassium nitrate
- 1-8 hours of wearing per day
- Achieves the material properties required by DIN EN ISO 28399 with daily application over the course of 10 days



#### Starter pack

- 1 Syringe @ 5 ml
- 1 Tray box

Flairesse Bleaching Gel CP 5% REF 220550 Flairesse Bleaching Gel CP 10% REF 220551

#### Refill pack

1 Syringe @ 5 ml

Flairesse Bleaching Gel CP 5% REF 220552 Flairesse Bleaching Gel CP 10% REF 220553

#### Cosmetic bleaching gel for gentle whitening

Indications: Gentle whitening of intrinsic and extrinsic discolouration

#### **Product benefits**

- At-home bleaching for use with a custom tray
- Mild formulation with 16 % carbamide peroxide (CP) for gentle whitening
- Bonds to the tray effectively and does not run
- Contains xylitol and potassium nitrate
- 1-8 hours of wearing per day
- Achieves the material properties required by DIN EN ISO 28399 with daily application over the course of 10 days



#### Starter pack

 $1\,$  Syringe @ 5 ml Flairesse Bleaching Gel CP 16%

1 Tray box REF 310000

#### Refill pack

1 Syringe @ 5 ml Flairesse Bleaching Gel CP 16% REF 310001

#### Flairesse Paste



# Tooth cleaning and polishing paste for professional tooth cleaning (PTC) and polishing of tooth surfaces as part of a prophylactic treatment

Indications: removal of superficial soft and hard plaque, removal of extrinsic stains

#### **Product benefits**

- For a perfect cleaning and polishing effect
- Does not splatter
- 3 degrees of abrasion: fine, medium, coarse
- Contains xylitol and fluoride (12,300 ppm)
- Two fresh flavours: Mint and Melon





#### Pack

200 Single-use unit @ 1.8 g

1 Prophy ring

Mint fine

Mint medium Mint coarse	REF 220415 REF 220416
1 Tube @ 75 ml	
Mint fine	REF 220417
Mint medium	REF 220418
Mint coarse	REF 220419
Melon fine	REF 220423
Melon medium	REF 220424
Melon coarse	REF 220425

#### Flairesse Gel/Foam



#### Prophylaxis gel and foam for final treatment as part of professional tooth cleaning (PTC) and for integration into an individual prevention plan

Indications for gel and foam: promotes remineralisation of the enamel after professional tooth cleaning, remineralisation of incipient caries

#### **Product benefits**

- Increases enamel resistance against acid attacks
- Contains xylitol and fluoride (12,300 ppm)
- Does not drip
- Short application time: only 1 minute
- Can be used in a tray
- Three fresh flavours: Mint, Melon, and Strawberry





#### Pack

RFF 220414

1 Bottle @ 480 ml Gel	
Mint	REF 225411
Melon	REF 225412
Strawberry	REF 225413
1 Dispenser @ 125 g Foam	
Mint	REF 220408
Melon	REF 220409
	220.00

## Flairesse Varnish



#### Moisture and saliva-tolerant prophylaxis varnish

Indications: Treating hyper-sensitive teeth, sealing overly sensitive tooth necks and exposed root dentine surfaces, as a protective coating to prevent caries, remineralisation of incipient caries



#### **Product benefits**

- Reliably eliminates hypersensitivity
- Contains xylitol and fluoride (22,600 ppm)
- Invisible on teeth
- Easy to use
- Two fresh flavours: Mint and Melon



#### Pack

35 Single-use units @ 0.4 ml Mint Melon

REF 220404 REF 220405

#### **Icon Proximal**



#### Caries infiltrant for micro-invasive treatment of carious lesions in the proximal region

Indication: micro-invasive treatment of early, non-cavitated caries in the proximal area (maximal radiological lesion progression into the outer third of dentine)



#### **Product benefits**

- Stops caries at an early stage in the hard-to-reach interdental space
- · Microinvasive, gentle treatment
- Preserves healthy tooth structure, tooth is given a longer life expectancy
- · Comfortable for patients: no need for anaesthesia or drilling



- 1 Treatment unit contains:
- 1 Syringe @ 0.45 ml Icon Etch; 1 Syringe q 0.45 ml Icon Prime;
- 1 Syringe @ 0.45 ml Icon Infiltrant; 6 Icon Proximal Tips; 1 Luer-Lock Tip;
- 4 Icon Wedges

#### Starter pack 1 Treatment unit

1 Treatment unit	REF 220523
Pack	
2 Treatment units	REF 220524

#### Economy pack 7 Treatment units REF 220525

#### Accessories

See right









#### Icon Vestibular



#### Infiltrant for aesthetic treatment for developmental enamel opacities

Indications: minimally invasive treatment of non-cavitated carious enamel lesions, aesthetic treatment of developmental enamel opacities such as caries-related white spots, MIH, fluorosis and trauma





#### **Product benefits**

- Minimally invasive treatment
- Gentle, pain-free and substance-preserving
- · For MIH, caries-related white spots, fluorosis and trauma
- Effectiveness clinically proven in over 400 studies
- A holistic treatment plan: optimise infiltration treatment with medical bleaching and transillumination



- 1 Treatment unit contains:
- 1 Syringe @ 0.45 ml Icon Etch; 1 Syringe @ 0.45 ml Icon Prime;
- 1 Syringe @ 0.45 ml lcon Infiltrant; 6 lcon Vestibular Tips;
- 1 Icon Capillary Tip

50 Icon Wedges

#### Starter pack

1	Treatment unit	REF 220520
Pac	k	
2	Treatment units	REF 220521
Eco	nomy pack	
	Treatment units	REF 220522
Refill pack		
3	Syringes @ 0.45 ml lcon Etch	REF 220530
3	Syringes @ 0.45 ml Icon Prime	REF 220531
3	Syringes @ 0.45 ml lcon Infiltrant	REF 220532
20	Icon Vestibular Tips	REF 220533
20	Icon Proximal Tips	REF 220534
20	Icon Capillary Tips	REF 220535

## **DMG MiniDam**



#### Silicone gingiva protection

Indication: silicone protection to facilitate a dry environment in the proximal area between two teeth to protect the gingiva against dental materials



#### **Product benefits**

- Reliable local protection of the proximal area
- Quick and easy
- Comfortable for patients
- Secure hold without clamps
- Latex-free
- Assists with infiltration treatment

Pack 20 Pcs





## **Ecosite Elements**



## Nano-hybrid composite system for the anterior and posterior region

Indications: loss of tooth substance caused by caries, fractures, dental wear or impaired development; diastema, tooth discolouration or slight misalignment of the front teeth; dental or primary and secondary prevention of caries

#### **Product benefits**

- Innovative and modular shade system
- Exceptional handling thanks to NC1 technology
- Optimum consistency for simple and rapid modelling
- Rapid adaptation to the cavity walls thank to the Push-and-Flow effect
- · Rapid high-gloss polishing for gleaming shine
- · New: the light shades A1 and Bleach Light







(Enamel Light, Enamel Medium)
Refill packs @ 2 g (Incisal, A2 Opaque)
Refill pack Ecosite Bulk Fill, shade (Universal)

Syringes Safetips

	ter pack Safetips @ 0.25 g Ecosite Elements Pure shades (16 of each): B1, A1, A2, A3, A3.5	REF 220861
5	Syringes @ 4 g Ecosite Elements Pure shades: B1, A1, A2, A3, A3.5	REF 220860
64	Safetips @ 0.25 g, Ecosite Elements Layer shades (16 of each): Enamel Bleach, Enamel Light, Enamel Medium, Enamel Dark	REF 220879
4	Syringes @ G 4 g Ecosite Elements Layer shades: Enamel Bleach, Enamel Light, Enamel Medium, Enamel Dark	REF 220874
4	Syringes Ecosite 2 g Ecosite Elements Highlight shades: Incisal, A2 Opaque, White, Brown	REF 220889
1 1	nomy pack – Bundle S Starter pack Ecosite Elements Pure Starter pack Ecosite Elements Layer Starter pack Ecosite Elements Layer Syringes Safetips	REF 111260 REF 111261
1 1 3 1	Nomy pack – Bundle XL Starter pack Ecosite Elements Pure Starter pack Ecosite Elements Layer Refill packs of Ecosite Elements Pure (A1, A2, A3) Ecosite Elements Drawer Insert Shade Guide Layer Syringes Safetips	REF 111262 REF 111263
3	nomy pack – Bundle base Packs of Ecosite Elements Pure (A2, A3, A3.5) Refill packs Ecosite Elements Layer	

#### Refill pack

16	Safetips @ 0.25 g	
	Ecosite Elements Pure, shade A1 Ecosite Elements Pure, shade A2	REF 220851 REF 220852
	Ecosite Elements Pure, shade A3	REF 220853
	Ecosite Elements Pure, shade A3.5	REF 220854
	Ecosite Elements Pure, shade B1	REF 220855 REF 220850
	Ecosite Elements Pure, shade BL1 Ecosite Elements Layer, shade Enamel Bleach	REF 220875
	Ecosite Elements Layer, shade Enamel Light	REF 220876
	Ecosite Elements Layer, shade Enamel Medium	REF 220877
	Ecosite Elements Layer, shade Enamel Dark	REF 220878
1	Syringe @ 4 g	
	Ecosite Elements Pure, shade A1	REF 220841
	Ecosite Elements Pure, shade A2	REF 220842
	Ecosite Elements Pure, shade A3	REF 220843
	Ecosite Elements Pure, shade A3.5 Ecosite Elements Pure, shade B1	REF 220844 REF 220845
	Ecosite Elements Pure, shade BL1	REF 220840
	Ecosite Elements Layer, shade Enamel Bleach	REF 220870
	Ecosite Elements Layer, shade Enamel Light	REF 220871
	Ecosite Elements Layer, shade Enamel Medium	REF 220872
	Ecosite Elements Layer, shade Enamel Dark	REF 220873
1	Syringe @ 2 g	
	Ecosite Elements Highlight, shade Incisal	REF 220885
	Ecosite Elements Highlight, shade A2 Opaque	REF 220886
	Ecosite Elements Highlight, shade White Ecosite Elements Highlight, shade Brown	REF 220887 REF 220888
	0 0 7	NLI 220000
	<b>essories</b> Safetip Dispenser Pro	REF 220829
	Ecosite Elements Shade Guide Layer	REF 220829
	Ecosite Elements Shade Guide BL1	REF 220836
1	Ecosite Elements Drawer Insert	REF 220834
20	Luer-Lock Tips Silver	REF 220712

REF 111270

#### **Ecosite One**



#### **Ecosite Bulk Fill**



## Stopperable one-shade nano-hybrid composite for the posterior region

Indications: fillings in the posterior region (classes I and II including replacement of individual cusps), class V fillings (cervical caries, root erosion, conical defects), core build-ups, deciduous teeth fillings, extended crack fillings on molars and premolars



#### **Product benefits**

- Thanks to special one-shade technology, all common tooth colours can be economically and aesthetically supplied with a universal colour
- Reliable and complete curing to 3 mm
- Exceptional handling thanks to DMG's proprietary NC1 technology
- Easy modeling and precision contact point formation
- In most cases no blocker / opaquer necessary
- Outstanding mechanical properties thanks to special filler technology
- Excellent polishability





#### Pack

1 Syringe @ 4 g 16 Safetips @ 0.25 g REF 220900 REF 220901

#### Packable nano-hybrid composite

Indications: fillings in the posterior region (classes I and II including replacement of individual cusps), class V fillings (cervical caries, root erosion, conical defects), core build-ups, deciduous teeth fillings, extended crack fillings on molars and premolars



#### **Product benefits**

- Can be applied in one layer without an additional top layer
- Reliable and complete curing to 5 mm
- Low shrinkage stress values
- Easy modelling
- Precision contact-point formation
- Outstanding mechanical properties thanks to special filler technology
- Excellent polishability





1 Safetip Dispenser

1 Safetip Dispenser Pro

#### Pacl

ac	K	
16	Safetips @ 0.25 g	
	Shade Light	REF 220755
	Shade Universal	REF 220756
	Shade Contrast	REF 220757
1	Syringe @ 4 g	
	Shade Light	REF 220752
	Shade Universal	REF 220753
	Shade Contrast	REF 220754
Accessories		

REF 213122

#### Constic



## **EcuSphere Flow**



## 3-in-1 flow composite: etching, bonding and filling in one single step

Indications: small class I and III restorations, class V restorations, sealing of fissures, extended sealing of fissures, underfilling of class I-II restorations, blocking out and filling of undercuts in cavities, small occlusal cavities in deciduous teeth

#### **Product benefits**

- · Self-etching and self-adhesive
- Faster, easier and more efficient treatment process
- Radiopacity and tooth-like fluorescence
- Prevents postoperative sensitivities
- Minimises potential sources of mistakes
- Contains MDP



#### Pack

- 2 Syringes @ 2 g
- 1 Brush holder
- 20 Clip-on Brushes
- 20
   Luer-Lock Tips Silver

   Shade A1
   REF 220700

   Shade A2
   REF 220701

   Shade A3
   REF 220702

   Shade Opaque White
   REF 220705

#### Accessories

20 Clip-on Brushes

20 Luer-Lock Tips Silver REF 220712

#### Light-curing flowable universal composite

Indications: small restorations (cavity classes I, III, IV, V), cavities in deciduous teeth, minimally invasive fillings, underfillings; extended fissure sealing, blocking out and filling of undercuts, repair of composite restorations

#### **Product benefits**

- Low-viscosity: fluid during and stable after the application
- Direct application
- Excellent polishing properties
- Natural fluorescence

#### Speciality

Can be applied as a single shade and multi-shade; perfectly compatible with all EcuSphere products



#### Pack

2 Syringes @ 1.5 g 10 Luer-Lock Tips

U	Luoi Look ripo	
	Shade A1	REF 220315
	Shade A2	REF 220314
	Shade A3	REF 220313
	Shade A3.5	RFF 220312

#### Accessories

60 Luer-Lock Tips Black REF 211477

## **EcuSphere Carat**



## **EcuSphere Shape**



## Light-curing packable composite for posterior regions

Indications: restorations of all cavity classes, especially classes I and II; manufacture of direct inlays and onlays

#### **Product benefits**

- · High viscosity: good for packing
- Excellent handling
- Does not stick to the instrument
- Good adaptation to cavity walls
- Excellent stability
- Natural fluorescence

#### Speciality

Recommended for posterior applications





Vdb

raci	N.	
20	Safetips @ 0.3 g	
	Shade A1	REF 220323
	Shade A2	REF 220271
	Shade A3	REF 220322
	Shade A3.5	REF 220270
1	Syringe @ 3 g	
	Shade A1	REF 220287
	Shade A2	REF 220267
	Shade A3	REF 220285
	Shade A3.5	REF 220266
Acc	essories	
1	Safetip Dispenser	REF 213122
1	Safetip Dispenser Pro	REF 220829

#### Light-curing universal composite

Indications: restorations of all cavity classes in the anterior and posterior region

#### **Product benefits**

- Medium viscosity: good modelling properties
- Great variety of shades
- Excellent sculpting properties
- Single and multi-shade layer technique
- Natural fluorescence
- Excellent polishing properties





#### Pack

20	Safetips @ 0.3 g Shade A1 Shade A2 Shade A3 Shade A3.5 Shade A4 Shade B1	REF 220327 REF 220309 REF 220307 REF 220298 REF 220346 REF 220328
1	Syringe @ 3 g Shade A1 Shade A2 Shade A3 Shade A3.5 Shade A4 Shade B1	REF 220334 REF 220308 REF 220299 REF 220297 REF 220347 REF 220335
1	<b>essories</b> Safetip Dispenser Safetip Dispenser Pro	REF 213122 REF 220829

#### DeltaFil



## Conventional, self-adhesive glass ionomer cement with an innovative material concept

Indications: permanent restoration: Class I, II\* and V; including for sandwich technique; temporary restoration: Class I, II, III and V; Deciduous tooth: Class I, II, and III; extended fissure sealing; base/liner; core build-up

\* Suitable for areas subject to occlusal load, provided the isthmus measures less than half of the intercuspal area.

#### **Product benefits**

- Increased flexural strength and fracture toughness thanks to innovative and patented micelle technology
- Ideal for difficult isolation, restricted patient cooperation and a difficult-to-access cavity
- · Easy to use, no absolute drying necessary
- Self- adhesive no bonding necessary
- Range of shades (White and Dark) adapted to suit specific patient group





#### Starter pack

1 Bottle @ 15 g powder

 $1\,$  Bottle @ 8 g liquid, mixing palette, measuring spoon

 Shade A1
 REF 520030
 Shade A4
 REF 520034

 Shade A2
 REF 520031
 Shade White
 REF 520045

 Shade A3
 REF 520032
 Shade Dark
 REF 520046

#### Pack

50~ DeltaFil Capsules @ 370~mg , powder, 75~mg liquid (0.07 ml)

 Shade A1
 REF 520000
 Shade A4
 REF 520004

 Shade A2
 REF 520001
 Shade White
 REF 520005

 Shade A3
 REF 520002
 Shade Dark
 REF 520006

#### Refill pack

1 Bottle @ 15 g powder
Shade A1 REF 520035 Shade A4 REF 520039
Shade A2 REF 520036 Shade White REF 520047
Shade A3 REF 520037 Shade Dark REF 520048
1 Bottle @ 8 g liquid REF 520040

#### Accessories

1 Capsule Dispenser REF 520051 1 Bottle @ 8 g DeltaFil Conditioner REF 520020

#### DeltaFil Coat



## Light-curing, TPO- and MMA-free protective coating for glass ionomer cements

Indications: protection of glass ionomer restorations

#### **Product benefits**

- Protects the glass ionomer cement from drying out and from being contaminated with moisture
- · Ensures optimal surface quality and lustre
- · Reduces the need for fine polishing
- Improves the aesthetics of the restoration



#### Pack

1 Bottle @ 5 g DeltaFil Coat

#### Accessories

50 Single-use Brushes White REF 109320

# **DeltaFil Conditioner**

# **Ionosit Baseliner**





# Polyacrylic acid solution for removing the smear layer after preparation and conditioning the cavity

Indications: conditioning of the cavity and treatment of the smear layer in DeltaFil restorations

### **Product benefits**

- Optimises the bond between DeltaFil and the tooth structure
- Ensures high durability of the restoration



### Pack

1 Bottle @ 8 g

REF 520020

# Light-curing underfilling material for improved longevity of restorations

Indications: base/liner material for composite restorations, and for blocking out undercuts

### **Product benefits**

- Gap-free restorations
- Compensation of polymerisation shrinkage
- Continuous fluoride release
- Antibacterial properties (due to zinc content)
- No complicated handmixing
- Protects the pulp
- Seals the dentine tubules
- Prevents postoperative sensitivities
- Radiopaque

### Speciality

Also available as a Saferinge for easier and faster application





### Pack

2 Coforing co. @ 0.22 g	
3 Saferinges @ 0.33 g REF	220091
2 Syringes @ 1.5 g REF	213002

### Economy pack

20	Saferinges @ 0.33 g	REF 209538
6	Syringes @ 1.5 a	REF 213001

### Accessories

60 Luer-Lock Tips Black REF 211477

### LuxaBond Universal



### Dual curing two-bottle universal adhesive system. Can also be used as a primer for indirect restorations

Indications: As an adhesive: for direct restorations using light-curing composites and compomers, light-curing fissure-sealing materials, build-up materials (light, self, dual-curing), cementing indirect restorations and root posts with composite based luting cements (light, self, dual-curing), cementing veneers with adhesive composites using the melt etching technique

As a primer: for indirect metal restorations (precious and non-precious metals), zirconia/aluminium oxide, silicate ceramics and composites

### **Product benefits**

- Truly comprehensive in every respect
- Contains MDP
- Strong adhesion
- · Minimises sources of error
- Study-demonstrated clinical use with superior adhesion

The ideal partner for optimum bonding: PermaCem Universal (page 62), LuxaCore Z Dual (page 38)



### Pack

- 1 Bottle @ 5 ml Bond A
- 1 Bottle @ 4.5 ml Bond B
- 50 Single-use Brushes White
- 25 Endo Brushes
- 1 Mixing palette REF 213379

### Accessories

50 Single-use Brushes White REF 107291 25 Endo Brushes Black REF 212386

# **LuxaBond Total Etch**







# Dual-curing adhesive system designed for all direct and indirect techniques

Indications: dual-curing bonding agent, in combination with the total etch technique, on enamel and dentine for all direct and indirect techniques, such as challenging adhesive cementations of root posts and core build-ups, inlays, onlays, crowns and bridges

### **Product benefits**

- Self-curing: hardens even without light
- No premature polymerisation of bonding »pools«
- Pre-bond for quick curing of the bonding agent in the dentine: superior adhesion values, longer working time with the bonding agent
- Does not contain acetone

### Speciality

Ideal as a partner of LuxaCore Z Dual (page 38), PermaCem (page 63), PermaCem Dual (page 63)



### Starter pack

- 1 Bottle @ 5 ml Pre-bond
- 1 Bottle @ 5 ml Bond A
- 1 Bottle @ 5 ml Bond B
- 1 Syringe @ 2 ml DMG Etching Gel
- 25 Endo Brushes Red
- 25 Endo Brushes Black
- 3 Luer-Lock Tips White
- 1 Mixing palette REF 212384

### Refill pack

1 Bottle @ 5 ml Pre-bond REF 212389

25	5 Endo Brushes Red	REF 212385
25	5 Endo Brushes Black	REF 212386
50	D Luer-Lock Tips White	REF 109942

### **Ecosite Bond**



### Light-curing universal adhesive

Indication: bonding agent for direct restorations of all cavity classes with light-curing composites, compomers and ormocers

### **Product benefits**

- Freedom of choice when it comes to the etching mode: self-etch, total etch, selective etch
- Contains MDP
- Optimum bonding properties in every situation
- Perfect as a partner for Ecosite One, Ecosite Elements, Ecosite Bulk Fill and Vitique



### Pack

1 Bottle @ 5 ml Ecosite Bond

50 Single-use Brushes White

REF 220758

### Accessories

50 Single-use Brushes White

REF 107291

### Contax







### Self-etching and self-conditioning adhesive system

Indication: light-curing adhesive system for light-curing restorative materials to be used on dentine and enamel

### Contax with activator to mix with Contax Bond

Indication: can be used with all self and dual-curing restoratives

### **Product benefits**

- No additional etching
- Priming and bonding in just two steps
- Secure long-term adhesion because of separate priming step
- Reliable depth of demineralisation
- Water-based primer, thus tolerant as regards technique
- No postoperative sensitivities

### Speciality

Contax must always be light-cured, even when used with the activator



### Starter pack

- 1 Bottle @ 5 ml Primer
- 1 Bottle @ 5 ml Bond
- 50 Single-use Brushes White
- 50 Single-use Brushes Red

1 Mixing palette REF 220263

### Starter pack + activator

- 1 Bottle @ 5 ml Primer
- 1 Bottle @ 5 ml Bond
- 1 Bottle Activator (5 ml)
- 50 Single-use Brushes White
- 50 Single-use Brushes Red

1 Mixing palette REF 220282

### Refill pack

 1 Bottle @ 5 ml Contax Primer
 REF 220281

 1 Bottle @ 5 ml Contax Bond
 REF 220284

 1 Bottle @ 5 ml Contax Activator
 REF 220283

### Accessories

50 Single-use Brushes White REF 107291

### LuxaCore Z Dual



# Dual-curing material for root post cementation and core build-up

Indications: Core build-up, post cementation



### **Product benefits**

- Dentine-like trimmability: score-free preparation
- Particularly high compressive strength: high stability under the crown
- Excellent flow properties: optimum adaptation
- · Dual-cure: flexible working
- Immediate high final hardness
- Direct application with Intraoral Tips and Endo Tips: specially designed for root canals
- · Available in three shades





### Pack

2 Smartmix Syringes @ 9 g, 10 Smartmix Tips Endo, 10 10 Smartmix Tips Combi

	Shade Opaque Light Shade A3 Shade Blue	REF 213336 REF 213334 REF 213335
2	Smartmix Syringes @ 9 g, 10 Smartmix Tips Combi Shade Opaque Light Shade A3 Shade Blue	REF 233336 REF 2333336 REF 2333335
	Cartridge @ 95 g, 25 Automix Tips Yellow 1:1 Intraoral Tips Yellow, 10 Endo Tips Yellow Shade A3 Shade Blue	REF 213430 REF 213431
	Cartridge @ 48 g, 30 Automix Tips Yellow 1:1 Intraoral Tips Yellow, 10 Endo Tips Yellow Shade Opaque Light	REF 213332

### Accessories

See right

Shade A3 Shade Blue

# LuxaCore Z Post System



### System pack for post-endodontic treatment

Indications:

LuxaCore Z Dual, see page 38, left column; LuxaPost, see page 39; LuxaBond Total Etch, see page 36

### **Product benefits**

- Three top products for post-endodontic treatments in one pack
- In addition: CD with application video
- · Handy to use, well arranged











### System pack

LuxaPost

- 5 Posts ø 1.25 mm
- 1 Drill ø 1.25 mm

### LuxaBond Total Etch

- 1 Bottle @ 5 ml Pre-bond
- 1 Bottle @ 5 ml Bond A
- 1 Bottle @ 5 ml Bond B
- 25 Endo Brushes Red
- 25 Endo Brushes Black
- 1 Mixing palette

### DMG Etching Gel

- 1 Syringe @ 2 ml DMG Etching Gel
- 3 Luer-Lock Tips

### LuxaCore Z Dual

1 Smartmix Syringe @ 9 g, shade A3

20 Smartmix Tips Combi REF 110862
Accessories

C33011C3	
Automix Dispenser Type 25 1:1	REF 110253
Automix Tips Yellow 1:1	REF 909200
Automix Tips Endo L	REF 110963
Intraoral Tips Yellow	REF 909222
Endo Tips Yellow	REF 212280
Smartmix Tips Combi	REF 212042
Smartmix Tips Endo	REF 212041
Smartmix Tips Endo L	REF 213083
	Automix Dispenser Type 25 1:1 Automix Tips Yellow 1:1 Automix Tips Endo L Intraoral Tips Yellow Endo Tips Yellow Smartmix Tips Combi Smartmix Tips Endo

REF 213330

# LuxaPost



# Glass fibre-reinforced pre-silanised composite root post with dentine-like flexural modulus

Indication: manufacture of a reliable retention of the core build-up in cases of insufficient coronal tooth substance

### **Product benefits**

- Perfect partner for LuxaCore Z Dual: monoblock
- Dentine-like flexural modulus: reduced risk of root fracture
- Radiopaque
- Pre-silanised
- Gentle on the tooth substance thanks to tapered form
- One effective drill head for each size: no pre-drilling required
- Highly aesthetic



### Starter pack

- $5 \;\; \text{Posts} \, \emptyset \, 1.25 \, \text{mm} \, \text{Black}$
- 5 Posts 1.375 mm Purple
- 5 Posts ø1.5 mm Green
- 1 Drill ø 1.25 mm
- 1 Drill ø 1.375 mm
- $1~\text{Drill}\, \text{\emptyset}\, 1.5~\text{mm}$
- 30 Depth marker

1 Measuring card REF 110780

### Refill pack

en	праск	
5	Posts ø 1.25 mm Black	REF 110781
5	Posts 1.375 mm Purple	REF 110782
5	Posts ø1.5 mm Green	REF 110783
5	Posts 1.75 mm Red	REF 110961
1	Drill ø 1.25 mm	REF 110734
1	Drill ø 1.375 mm	REF 110735
1	Drill ø 1.5 mm	REF 110736
1	Drill ø 1.75 mm	REF 110962

# Safetip Dispenser



Dispenser for direct application of all commercially available composite tips

### **Product benefits**

- Safe loading of Safetips
- Direct application
- Hygienic application
- Autoclavable
- Ergonomic handling

### Speciality

The following products can be used with the Safetip Dispenser:

- Ecosite Elements
- Ecosite Bulk Fill
- Ecosite One
- EcuSphere Carat
- EcuSphere Shape
- DMG Retraction Paste



### Pack

1 Safetip Dispenser

REF 213122

# Safetip Dispenser Pro



High-quality dispenser for ergonomic, direct application of all commercially available composite tips

### **Product benefits**

- Easy dispensing thanks to parallel handles
- Precise treatment through straight-line motion
- Ergonomic handling, even for smaller hands
- High-quality materials and processing
- Secure insertion of Safetips via metal locking mechanism
- Hygienic application
- Autoclavable

### Speciality

The following products can be used with the Safetip Dispenser:

- Ecosite Elements
- Ecosite Bulk Fill
- Ecosite One
- EcuSphere Carat
- EcuSphere Shape
- DMG Retraction Paste



### Pac

1 Safetip Dispenser Pro

# **DMG Etching Gel**



### Universal phosphoric acid gel

Indication: etching of enamel and dentine

### **Product benefits**

- 37 % phosphoric acid gel
- Blue for easy identification
- Ideal retention pattern
- Optimal viscosity
- · Perfect dosing and precise application




### Pack

2 Syringes @ 2 ml 6 Luer-Lock Tips White

REF 212873

### Accessories

50 Luer-Lock Tips White



# **Honigum Pro Light**



# Honigum Pro Mono



# High-precision, VPS-based correction impression material

Indications: crown and bridge impressions, inlay and onlay impressions

### **Product benefits**

- Intelligent flow properties and extraordinary stability
- Good hydrophilicity
- Registers even the finest details
- Extremely high toughness
- Easily readable shade contrast



### Pack

2 Cartridges @ 50 ml

10 Automix Tips Yellow 1:1

10	Intraoral Tips Yellow	
	Honigum Pro Light	REF 929831
	Honigum Pro Light Fast	REF 929832

### Economy pack

8 Cartridges @ 50 ml

40 Automix Tips Yellow 1:1

40 Intraoral Tips Yellow
Honigum Pro Light REF 929835
Honigum Pro Light Fast REF 929836

### Accessories

1 Automix Dispenser Type 50 1:1	REF 909507
50 Automix Tips Yellow 1:1	REF 909200
50 Intraoral Tips Yellow	REF 909222

### Monophase VPS-based impression material

Indications: crown and bridge impressions, inlay, onlay and functional impressions, as well as pressure-spot indications and all types of pickup impressions (e.g. for implants)

### **Product benefits**

- One-step precision
- Trouble-free and safe removal
- Highly hydrophilic and superior flow behaviour
- Precise reproduction of the preparation margin
- Very high tear resistance





### Pack

	· =
1	Cartridge @ 380 ml
ι Λ	MiyCtor Tipo Pluo

TO	Mixotal Tipo bluc	
	Honigum Pro Mono	REF 929568
	Honigum Pro Mono Fast	REF 929978

- 2 Cartridges @ 50 ml
- 4 Automix Tips Green 1:1
- 4 Intraoral Tips Yellow

   Honigum Pro Mono
   REF 929830

   Honigum Pro Mono Fast
   REF 929976

### Economy pack

5 Cartridges @ 380 ml

50	MixStar Tips Blue	
	Honigum Pro Mono	REF 929569
	Honigum Pro Mono Fast	REF 929979

REF 929834

- 8 Cartridges @ 50 ml
- 16 Automix Tips Green 1:1
- 16 Intraoral Tips Yellow Honigum Pro Mono

essories	
MixStar eMotion 230 V Int.	REF 900480
MixStar Tips Blue	REF 903334
Automix Dispenser Type 50 1:1	REF 909507
Automix Tips Green 1:1	REF 909490
Intraoral Tips Transparent	REF 909400
Bottle @ 10 ml DMG Tray Adhesive	REF 909394
	MixStar eMotion 230 V Int. MixStar Tips Blue Automix Dispenser Type 50 1:1 Automix Tips Green 1:1 Intraoral Tips Transparent

# **Honigum Pro Heavy**



# **Honigum Pro Putty**



### VPS-based pre-impression material

Indications: crown and bridge impressions, inlay and onlay impressions, all types of pickup impressions (e.g. for implants)

### **Product benefits**

- High final hardness
- · Convenient insertion resistance
- Trouble-free and safe removal
- Consistently homogeneous, bubble-free material quality





Honigum Pro Heavy

# Pack 1 Cartridge @ 380 ml, 10 MixStar Tips Blue

	Honigum Pro Heavy Fast	REF 929638
2	Cartridges @ 50 ml, 4 Automix Tips Green 1:1 Honigum Pro Heavy Honigum Pro Heavy Fast	REF 929766 REF 929838
	nomy pack	
5	Cartridges @ 380 ml, 50 MixStar Tips Blue Honigum Pro Heavy Honigum Pro Heavy Fast	REF 929538 REF 929637
8	Cartridges @ 50 ml, 16 Automix Tips Green 1:1 Honigum Pro Heavy	REF 929837
Acc	essories	
_	MixStar eMotion 230 V Int. MixStar Tips Blue	REF 900480 REF 903334
	Automix Dispenser Type 50 1:1	REF 903534
	Automix Tips Green 1:1	REF 909490
1	Bottle @ 10 ml DMG Tray Adhesive	REF 909394

# Precision, VPS-based impression material with a putty-like consistency

Indications: crown and bridge impressions, inlay and onlay impressions

### **Product benefits**

- Outstanding recovery after deformation
- Excellent dimensional stability
- Hand-mouldable or putty-like consistency from the DMG MixStar cartridge
- Excellent tear strength





50 MixStar Tips Blue

1 Bottle @ 10 ml DMG Tray Adhesive

# Pack 1 Cartridge @ 380 ml

REF 929537

10	MixStar Tips Blue	REF 929374
2	Jars @ 450 ml Honigum Pro Putty Soft Honigum Pro Putty Soft Fast Honigum Pro Putty Rigid Fast	REF 929364 REF 929365 REF 929366
	nomy pack	
	Cartridges @ 380 ml MixStar Tips Blue	REF 929375
50	MixStal Tips Dide	INLI 929313
8	Jars @ 450 ml Honigum Pro Putty Soft Honigum Pro Putty Soft Fast Honigum Pro Putty Rigid Fast	REF 929369 REF 929370 REF 929371
Acc	essories	
1	MixStar eMotion 230 V Int.	REF 900480

REF 903334

# **Honigum Light**





# Honigum Autor State of the sta

**Honigum Mono** 

### VPS-based correction impression material

Indications: crown and bridge impressions, inlay and onlay impressions

### **Product benefits**

- Excellent flow properties
- Extraordinary stability
- Balanced hydrophilicity
- Registers even the finest details



### Pack

- 2 Cartridges @ 50 ml
- 12 Automix Tips Yellow 1:1
- 12 Intraoral Tips
  Honigum Light REF 909831
  Honigum Light REF 909832

### Economy pack

- 8 Cartridges @ 50 ml
- 50 Automix Tips Yellow 1:1
- 50 Intraoral Tips Yellow Honigum Light Honigum Light Fast

### Accessories

000001100		
1	Automix Dispenser Type 50 1:1	REF 909507
50	Automix Tips Yellow 1:1	REF 909200
50	Intraoral Tips Yellow	REF 909222

### Monophase VPS-based impression material

Indications: crown and bridge impressions, inlay, onlay and functional impressions, as well as pressure-spot indications and all types of pickup impressions

### **Product benefits**

- One-step precision
- Precise reproduction of the preparation line
- Finest surface details
- Lower final hardness
- Good hydrophilicity





### Pack

REF 909835

REF 909836

	Cartridge @ 380 ml MixStar Tips White	REF 909568
4	Cartridges @ 50 ml Automix Tips Green 1:1 Intraoral Tips Transparent	REF 909830
5	nomy pack Cartridges @ 380 ml MixStar Tips White	REF 909569
8	Cartridges @ 50 ml	

### Accessories

20 Automix Tips Green 1:120 Intraoral Tips Transparent

	00001100	
1	MixStar eMotion 230 V Int.	REF 900480
50	MixStar Tips White	REF 909567
1	Automix Dispenser Type 50 1:1	REF 909507
50	Automix Tips Green 1:1	REF 909490
50	Intraoral Tips Transparent	REF 909400
1	Bottle @ 10 ml DMG Tray Adhesive	REF 909394

# **Honigum Heavy**

# **Honigum Putty**





### VPS-based pre-impression material

Indications: crown and bridge impressions, inlay and onlay impressions, all types of pickup impressions (e.g. for implants)

### **Product benefits**

- Extremely high hardness
- · Convenient insertion resistance
- Excellent fixation characteristics
- Trouble-free and safe removal
- Consistently homogeneous, bubble-free material quality





### Pack

1 00	N.	
1	Cartridge @ 380 ml	
10	MixStar Tips White Honigum Heavy Honigum Heavy Fast	REF 909537 REF 909638
2	Cartridges @ 50 ml	
4	Automix Tips Green 1:1	
	Honigum Heavy	REF 909766
	Honigum Heavy Fast	REF 909838
Eco	nomy pack	
5	Cartridges @ 380 ml	
50	MixStar Tips White	
	Honigum Heavy	REF 909538

### 16 Automix Tips Green 1:1 Honigum Heavy Honigum Heavy Fast

### Accessories

Honigum Heavy Fast 8 Cartridges @ 50 ml

1	MixStar eMotion 230 V Int.	REF 900480
50	MixStar Tips White	REF 909567
1	Automix Dispenser Type 50 1:1	REF 909507
50	Automix Tips Green 1:1	REF 909490
1	Bottle @ 10 ml DMG Tray Adhesive	REF 909394

### Precision, VPS-based impression material with a putty-like consistency

Indications: crown and bridge impressions, inlay and onlay impressions

### **Product benefits**

- Real putty-like consistency from the DMG MixStar cartridge or mouldable by hand
- Excellent setting as soon as mixing begins
- Excellent tear strength
- Unique Snap-Set technology (long working time and short time in the mouth)
- Outstanding recovery after deformation
- Excellent dimensional stability





### Pack

REF 909637

REF 909837 REF 909839

	11	
	Cartridge @ 380 ml MixStar Tips Blue	REF 919374
	Jars @ 450 ml	NEI OTOOT I
	Honigum Putty Soft Honigum Putty Soft Fast Honigum Putty Rigid Fast	REF 919364 REF 919365 REF 919366
Eco	nomy pack	

### 5 Cartridges @ 380 ml 50 MixStar Tips Blue

8	Jars @ 450 ml	
	Honigum Putty Soft	REF 919369
	Honigum Putty Soft Fast	REF 919370
	Honigum Putty Rigid Fast	REF 919371

REF 919375

1	MixStar eMotion 230 V Int.	REF 900480
50	MixStar Tips Blue	REF 903334
1	Bottle @ 10 ml DMG Tray Adhesive	REF 909394

# Silagum Light



# Tried and tested correction impression material based on VPS

Indications: crown and bridge impressions, inlay and onlay impressions

### **Product benefits**

- Good hydrophilicity
- Flowable viscosity
- Two speeds: normal and fast-curing



### Pack

2 Cartridges @ 50 ml

12 Automix Tips Yellow 1:1

12 Intraoral Tips Yellow
Silagum Light REF 909713
Silagum Light Fast REF 909714

### Economy pack

8 Cartridges @ 50 ml
50 Automix Tips Yellow 1:1
50 Intraoral Tips Yellow
Silagum Light

Silagum LightREF 909715Silagum Light FastREF 909720

### Accessories

1 Automix Dispenser Type 50 1:1	REF 909507
50 Automix Tips Yellow 1:1	REF 909200
50 Intraoral Tips Yellow	REF 909222

# Silagum Medium



# Medium-viscosity impression material based on VPS

Indications: crown and bridge impressions, inlay and onlay impressions, relining and functional impressions, pressure-spot indication

### **Product benefits**

- Good flow properties
- Good stability
- Normal and fast-curing



### Pacl

2 Cartridges @ 50 ml

12 Automix Tips Yellow 1:1

12 Intraoral Tips Yellow
Silagum Medium REF 909716
Silagum Medium Fast REF 909717

1 Automix Dispenser Type 50 1:1	REF 909507
50 Automix Tips Yellow 1:1	REF 909200
50 Intraoral Tips Yellow	REF 909222

# Silagum Mono



### Monophase VPS-based impression material

Indications: crown and bridge impressions, inlay and onlay impressions, functional impressions

### **Product benefits**

- Reproduces the finest surface details thanks to its exceptional flow properties
- · Easy and precise application with Intraoral Tips
- Remains firm in the tray and on the preparation



### Pack

- 2 Cartridges @ 50 ml
- 4 Automix Tips Yellow 1:1
- 4 Intraoral Tips Yellow REF 909718

### Economy pack

- 8 Cartridges @ 50 ml
- 16 Automix Tips Yellow 1:1
- 16 Intraoral Tips Yellow REF 909719

### Accessories

1 Automix Dispenser Type 50 1:1	REF 909507
50 Automix Tips Yellow 1:1	REF 909200
50 Intraoral Tips Yellow	REF 909222

# Silagum Putty



### VPS-based pre-impression material

Indications: crown and bridge impressions, inlay and onlay impressions

### **Product benefits**

- Easy dosage
- Excellent recovery after deformation
- Easy removal from the mouth
- Exceptional kneadability
- Excellent dimensional stability
- Ideally suited for two-step impressions with Silagum Light/Fast



### Pack

1	. Tub	@	262	ml	base	pas	ste
1	Tub	@	262	ml	cataly	/st	paste

Silagum PuttyREF 909018Silagum Putty FastREF 909037Silagum Putty SoftREF 909019

### Economy pack

4 Tubs @ 262 ml base paste

4 Tubs @ 262 ml catalyst paste
Silagum Putty REF 909405
Silagum Putty Fast REF 909038
Silagum Putty Soft REF 909036

8 Tubs @ 262 ml base paste

8 Tubs @ 262 ml catalyst paste

Silagum Putty REF 909420 Silagum Putty Soft REF 909485

### Accessories

1 Bottle @ 10 ml DMG Tray Adhesive REF 909394

# Silagum Putty Soft

# StatusBlue









### VPS-based pre-impression material

Indications: crown and bridge impressions, inlay and onlay impressions

### **Product benefits**

- Exceptional dimensional stability
- Consistently homogeneous, bubble-free material quality
- Easy and clean handling
- Can be dispensed directly into the tray



### Pack

1	Cartridge @ 380 ml	
10	MixStar Tips White	REF 909534
_		

### Economy pack

5 Cartridges @ 380 ml
50 MixStar Tips White REF 909535

### Accessories

.00	00001100	
1	MixStar eMotion 230 V Int.	REF 900480
50	MixStar Tips White	REF 909567
1	Bottle @ 10 ml DMG Tray Adhesive	REF 909394

### VPS-based situation impression material

Indications: situation impressions, impressions for the manufacture of temporaries or vacuum-formed stents and opposing arch impressions

### **Product benefits**

- Very high dimensional stability
- Very long shelf life
- Repourable
- No time-consuming hand mixing
- Excellent recovery after deformation







REF 909791

### Pack

i ao	ik.	
	Cartridge @ 380 ml	DEE 000000
10	MixStar Tips White	REF 909693
	Cartridges @ 50 ml	DEE 000700
	Automix Tips Pink 4:1	REF 909790
Eco	nomy pack	
5	Cartridges @ 380 ml	
50	MixStar Tips White	REF 909694

### Accessories

8 Cartridges @ 50 ml 16 Automix Tips Pink 4:1

-100	40000001100		
1	MixStar eMotion 230 V Int.	REF 900480	
50	MixStar Tips White	REF 909567	
1	Automix Dispenser Type 50 10:1	REF 110411	
25	Automix Tins Pink 4:1	RFF 909797	

# MixStar eMotion





# Retraction Paste Retraction paste District pasts for graph resident DIMG

**DMG Retraction Paste** 

### Mixing unit

Mixing unit for mixing impression materials

### **Product benefits**

- Fully automatic
- 3 selectable speeds
- 2 timers for working and setting times
- Easy handling thanks to start/stop system
- Suitable for foil bags and cartridges



Packaging

1 MixStar eMotion 230 V Int.

REF 900480

### Astringent retraction paste

Indications: temporary gingival margin retraction and provision of a dry sulcus environment

### **Product benefits**

- Compule can be bent as needed for simplest application
- Optimum consistency optimal results
- Good visibility in the sulcus
- Fast, simple and efficient



### Pack

25 Compules @ 0.37 g REF 919381

### Accessories

1 Safetip Dispenser REF 213122 1 Safetip Dispenser Pro REF 220829

### **O-Bite**







# O-Bite & DMG

O-Bite Scan

### VPS-based bite registration material

Indications: bite registration material, material for application on bite forks for face-bows, key material, e.g. for intraoral needle-point registration

### **Product benefits**

- One of the hardest VPS-based bite registration materials
- Superior fracture resistance
- Fresh fruity orange scent
- Minimal bite resistance
- Fast-curing



### Pack

2 Cartridges @ 50 ml 12 Automix Tips Pink 1:1 REF 909765

### Economy pack

8 Cartridges @ 50 ml 50 Automix Tips Pink 1:1

### Accessories

1 Automix Dispenser Type 50 1:1	REF 909507
50 Automix Tips Pink 1:1	REF 909202
50 Conturation Tips	REF 909203

### Scannable, VPS-based bite registration material

Indications: bite registration material, material for application on bite forks for face-bows, key material, e.g. for intraoral needle-point registration

### **Product benefits**

- Scannable
- Extremely high final hardness and superior fracture stability
- Minimal bite resistance
- · Fresh fruity orange scent
- Can be used easily and precisely for a digital workflow without the need to purchase expensive technology
- · Easy handling, no pretreatment with powder or spray



### Pack

REF 909767

2 Cartridges @ 50 ml 12 Automix Tips Pink 1:1 REF 919017

REF 919018

### Economy pack

8 Cartridges @ 50 ml 50 Automix Tips Pink 1:1

1 Automix Dispenser Type 50 1:1	REF 909507
50 Automix Tips Pink 1:1	REF 909202
50 Conturation Tips	REF 909203

# LuxaBite



### Bite registration material based on bis-acrylates

Indications: Fixation of impression posts for implant impressions, bite registrations



### **Product benefits**

- The hardest bite registration material (Barcol hardness 25 after one hour)
- · No compression or flexing when mounting
- Tasteless and odourless
- Excellent stability
- Easy to trim



### Pack

1 Cartridge @ 50 ml 15 Automix Tips Yellow 10:1

REF 110560

### Economy pack

5 Cartridges @ 50 ml

75 Automix Tips Yellow 10:1 REF 110558

### Accessories

1 Automix Dispenser Type 50 10:1 REF 110411 25 Automix Tips Yellow 10:1 REF 110609

# **DMG Tray Adhesive**



### Tray adhesive for all VPS materials

Indication: tray adhesive for bonding VPS materials with metal and plastic trays

### **Product benefits**

- Promotes the adhesion of VPS-based impression materials to the impression tray
- Easy to apply and handle
- · Short setting time after application
- Blue colour visual control



### Pack

1 Bottle @ 10 ml

REF 909394

# Economy pack

10 Bottles @ 10 ml

# Silagum Comfort



### Permanently soft, VPS-based relining material

Indication: Relining for total, partial and implant-supported dentures



### **Product benefits**

- Excellent bond to the denture
- Permanently soft
- Superior tear strength
- Odourless and tasteless
- Ease of handling
- · Highly aesthetic due to shade stability





### Starter pack

- 1 Cartridge @ 25 ml
- 1 Automix Dispenser
- 5 Automix Tips Yellow 1:1
- 1 Bottle @ 5 ml Primer
- 2 Bottles @ 10 ml Varnish
- 25 Brushes; 5 Trimmers; 1 Holder REF 909087

### Refill pack

1 Cartridge @ 25 ml 5 Automix Tips Yellow 1:1	REF 909088
2 Bottles @ 10 ml Varnish	REF 909089
1 Bottle @ 5 ml Primer	REF 909090
1 Cartridge @ 50 ml 12 Automix Tips Pink 1:1	REF 909951

### Economy pack

8 Cartridges @ 50 ml REF 9099	5	j		)
-------------------------------	---	---	--	---

### Accessories

1	Automix Dispenser Type 25 1:1	REF 110253
1	Automix Dispenser Type 50 1:1	REF 909507
50	Automix Tips Yellow 1:1	REF 909200
50	Automix Tips Pink 1:1	REF 909202
1	Silicone Trimming Set	
0.5	Total and a second seco	

25 Trimmer

1 Holder REF 909047

# LuxaPick-up



### Self-curing relining composite

Indications: direct hard relining of dentures, fixation of secondary components and connection elements in implant-borne dentures

### **Product benefits**

- Direct hard relining of dentures
- Fixation of secondary parts and connection elements in implant-borne dentures
- Ultra-biocompatible: no MMA, no peroxide
- Tension-free, reliable fit
- Direct chairside use
- Ideal for repairing prostheses and customising intermediate prostheses too
- Colour: Pink

### Speciality

Luxatemp Glaze & Bond for an optimum adhesive bond and reliable sealing (see p. 57)



### Pack

1	Cartridge @ 76 g	
15	Automix Tips Blue 10:1	REF 214007

45	Automix Tips Blue 10:1	REF 110409
1	Automix Dispenser Type 50 10:1	REF 110411



# **Luxatemp MaxProtect**



# Temporary crown and bridge material with long-term stability on a bis-acrylate base

Indications: temporary crowns and bridges, inlays and onlays, partial crowns, veneers and long-term temporaries

### **Product benefits**

- Ideal for a longer healing phase
- 1 year wearing time (clinically tested)\*
- For extremely stable long-term provisionals
- Enduring shade stability for excellent aesthetics
- Exceptional stability
- Homogeneous and smooth surface qualities
- Maximum protection
- Highly durable
  - \* Source: clinical study (Reg. No. NCT03160950Bi)



### Pack

1 Cartridge @ 76 g (50 ml)

TO	Automix rips blue 10:1	
	Shade A1	REF 130100
	Shade A2	REF 130101
	Shade A3	REF 130102
	Shade A3.5	REF 130103
	Shade B1	REF 130104
	Shade Bleach Light	REF 130105

### Accessories

50	Smartmix Tips Green-Blue	REF 110990
45	Automix Tips Blue 10:1	REF 110409
1	Automix Dispenser Type 50 10:1	REF 110411

# **Luxatemp Star**







# Highly stable bis-acrylic temporary crown and bridge material

Indications: temporary crowns and bridges, inlays and onlays, partial crowns, veneers and long-term temporaries

### **Product benefits**

- For particularly stable and aesthetic temporaries that meet high standards
- Substantially higher fracture resistance and flexural strength
- Improved initial hardness
- · Faster achievement of final hardness
- High stability and durability
- Excellent aesthetics and long-term shade stability
- Great adaptation and fit
- Ideal for producing mock-ups





### Pack

1	Smartmix Syringe @ 15 g	
10	Smartmix Tips Green-Blue	
	Shade A2	REF 110901
	Shade A3	REF 110902
	Shade B1	REF 110904
	Shade Bleach Light	REF 110905
1	Cartridge @ 76 g	
15	Automix Tips Blue 10:1	
	Shade A1	REF 110906
	Shade A2	REF 110907

REF 110908

REF 110909

REF 110910

REF 110911

### Economy pack

Shade A3

Shade B1

Shade A3.5

Shade Bleach Light

5 Cartridges @ 76 g
75 Automix Tips Blue 10:1
Shade A1 REF 11.0954
Shade A2 REF 11.0912
Shade A3 REF 11.0913

### Accessories

See left

# **Luxatemp Fluorescence**



# **Luxatemp Plus**





# Bis-acrylic temporary crown and bridge material with natural fluorescence

Indications: temporary crowns, bridges, inlays, onlays, partial crowns and veneers

### **Product benefits**

- For naturally aesthetic temporaries
- Low polymerisation shrinkage
- High fracture and abrasion resistance
- Excellent shade stability
- Low allergy potential
- Suitable for long-term temporaries





### Pack

Accessories

50 Smartmix Tips Green-Blue

 $1 \ \ \text{Automix Dispenser Type } 50\ 10{:}1$ 

45 Automix Tips Blue 10:1

	IX.	
1	Smartmix Syringe @ 15 g	
10	Smartmix Tips Green-Blue	
	Shade A1	REF 110718
	Shade A2	REF 110719
	Shade A3	REF 110720
	Shade B1	REF 110722
	Shade Bleach Light	REF 110723
1	Cartridge @ 76 g	
15	Automix Tips Blue 10:1	
	Shade A1	REF 110585
	Shade A2	REF 110586
	Shade A3	REF 110710
	Shade A3.5	REF 110587
	Shade B1	REF 110588
	Shade Bleach Light	REF 110589
Eco	nomy pack	
5	Cartridges @ 76 g	
75	Automix Tips Blue 10:1	
	Shade A2	REF 110712
	Shade A3	REF 110713

### Bis-acrylic temporary crown and bridge material

Indications: temporary crowns, bridges, inlays, onlays, partial crowns and veneers

### **Product benefits**

- An established classic for immediate precision-fit temporaries
- · Simple and quick application into the impression
- Successful temporary restorations from the word go
- · Reliably accurate fit
- Optimised abrasion resistance
- · High flexural strength
- · Very good biocompatibility



### Starter pack

1 Cartridge @ 76 g 1 Automix Dispenser Type 50 10:1 15 Automix Tips Blue 10:1 Shade A2

REF110400

### Pack

1 Cartridge @ 76 g
15 Automix Tips Blue 10:1
Shade A1 REF 110401
Shade A2 REF 110402
Shade A3.5 REF 110403
Shade B1 REF 110404
Shade Bleach Light REF 110368

### Economy pack

5 Cartridges @ 76 g 75 Automix Tips Blue 10:1 Shade A1 Shade A2 Shade A3.5

REF 110405 REF 110406 REF 110407

### Accessories

45 Automix Tips Blue 10:1REF 1104091 Automix Dispenser Type 50 10:1REF 110411

REF 110990

REF 110409

### LuxaFlow Star



# Luxatemp Glaze & Bond



### Flowable light-curing composite

Indications: enhancement and repair of bis-acrylate temporaries; small class III, IV and V cavity fillings; minimally invasive fillings; underfillings



### **Product benefits**

- Natural fluorescence
- Excellent polishing properties
- Luer-Lock Tips for precise and direct applications
- Fluid during the application and stable thereafter
- Radiopaque



### Starter pack

4 Syringes @ 1.5 g 20 Luer-Lock Tips Black

Colours: A1, A2, A3.5, B1 REF 214000

### Refill pack

2 Syringes @ 1.5 g

10 Luer-Lock Tips Black

Shade A1	REF 214001
Shade A2	REF 214002
Shade A3	REF 214003
Shade A3.5	REF 214004
Shade B1	REF 214005

### Accessories

60 Luer-Lock Tips Black REF 211477

# Light-curing varnish for repairing and sealing temporary crowns and bridges

Indications: coating and repairing of provisional crowns and bridges, primer for indirect restorations made of PEEK/composite and PMMA

### **Product benefits**

- High gloss without polishing
- Excellent bond with repair materials
- High abrasion resistance
- Low surface roughness



### Pack

1 Bottle @ 5ml

25 Brushes

1 Mixing palette REF 211927

### Accessories

50 Brushes REF110041

# TempoCem ID



# TempoCem NE



### Aesthetic and transparent provisional luting cement

Indications: temporary luting of temporary and semi-permanent crowns, bridges, inlays, onlays, veneers and implant-supported interim prostheses



### **Product benefits**

- Invisible below the restoration
- Identifiable for easy, reliable excess removal in one go
- Optimum transparency for the highest aesthetic results
- Optimum bonding, easy removability
- User-friendly, easy to apply
- Temporary luting, semi-permanent reliability



### Pack

1 Smartmix Syringe @ 5 ml 20 Smartmix Tips Green

REF 213200

### Accessories

50 Smartmix Tips Green REF 110724

### Eugenol-free, zinc oxide-based provisional luting cement

Indications: temporary luting of provisional and permanent crowns, bridges, inlays and onlays

### **Product benefits**

- Universal use
- No influence on resin-based restoration materials
- Direct application into the temporary
- Good adhesion yet easy removal of temporary restoration
- Easy and clean removal of excess cement
- Exceptional fit through thin film thickness





2 Smartmix Syringes @ 5 ml

50 Automix Tips Blue 1:1

1 Automix Dispenser Type 25 1:1

### Pacl

20 Smartmix Tips Short	REF 212104
1 Cartridge @ 25 ml 35 Automix Tips Blue 1:1	REF 110272
Accessories	
50 Smartmix Tips Short	REF 212040
100 Automix Tips Blue 1:1	REF 110267

REF 110266

# **TempoCem**







### Eugenol zinc oxide-based provisional luting cement

Indications: temporary luting of provisional and permanent crowns, bridges, inlays and onlays as well as implant-borne prostheses

### **Product benefits**

- Direct application into the temporary
- Optimal adhesion and removability
- Easy and clean removal of excess





### Pack

	Smartmix Syringes @ 5 ml Smartmix Tips Short	REF 212129
	Cartridge @ 25 ml Automix Tips Blue 1:1	REF 110247
Accessories		
50	Smartmix Tips Short	REF 212040
100	Automix Tips Blue 1:1	REF 110267
50	Automix Tips Blue 1:1	REF 110266
1	Automix Dispenser Type 25 1:1	REF 110253





### PermaCem Universal



# Adhesive, dual-curing composite luting cement

Indications: permanent luting of inlays, onlays, crowns, bridges, veneers, root posts and restorations on implant abutments

### **Product benefits**

- Universal without compromise
- Universal use versatile with silicate and oxide ceramics, as well as metal and composites
- Clinically proven for all common CAD/CAM materials
- First-class aesthetics in three shades with matching
   Try-in Paste
- Strong adhesion across a broad range
- High abrasion resistance for reliably good marginal quality and longevity of the restoration
- Simple, stress-free removal of excess cement in one step

The ideal partner for optimum bonding: LuxaBond Universal (Page 36)



### Pack

1 Smartmix Syringe @ 9 g 10 Smartmix Tips Short Shade A1 Shade A2.5 Shade Transparent

### Accessories

50 Smartmix Tips Short REF 212040

### PermaCem 2.0



### Self-adhesive, dual-curing composite luting cement

Indications: permanent luting of inlays, onlays, crowns, bridges and root posts, permanent luting of crowns and bridges on implant abutments

### **Product benefits**

- An all-in-one solution
- Quick and simple to use without any etching or bonding
- Contains MDP
- A secure bond with the tooth and restauration
- Outstandingly strong adhesive to zirconia, silicate and oxide ceramics
- Minimal water absorption
- Optimum consistency
- Extraordinary flow under pressure
- Straightforward and stress-free excess removal
- Excellent radiopacity



### Pacl

REF 213401

REF 213402

REF 213404

1 Smartmix Syringe @ 9 g
15 Smartmix Tips Short
5 Smartmix Tips Endo M
Shade Transparent
Shade A 2 Universal

Shade Transparent REF 213368
Shade A2 Universal REF 213366
Shade A3 Opaque REF 213367

### Economy pack

3 Smartmix Syringes @ 9 g 45 Smartmix Tips Short 15 Smartmix Tips Endo M Shade Transparent Shade A2 Universal

Shade TransparentREF 213382Shade A2 UniversalREF 213380Shade A3 OpaqueREF 213381

### Accessories

50Smartmix Tips ShortREF 21204020Smartmix Tips Endo MREF 213373

# PermaCem Dual



# **PermaCem**



### Adhesive, dual-curing luting cement

Indications: final cementation of crowns, bridges, inlays and onlays

### **Product benefits**

- Reduced setting time due to light cure
- Universal use
- Also suitable for all-ceramic work
- Insoluble in a wet environment
- · Low film thickness for perfect fit
- Easy removal of excess





### Pack

Smartmix Syringes @ 10 g Smartmix Tips Short	REF 212035
Cartridge @ 52 g Automix Tips Blue 1:1	REF 110524

### Accessories

REF 212040
REF 110267
REF 110266
REF 110253

### Adhesive, self-curing luting cement

Indications: final cementation of crowns, bridges, inlays and onlays

### **Product benefits**

- Self-curing
- Universal use
- Also suitable for all-ceramic work
- Insoluble in a wet environment
- · Low film thickness for perfect fit
- Easy removal of excess





100 Automix Tips Blue 1:1 50 Automix Tips Blue 1:1

1 Automix Dispenser Type 25 1:1

### Pack

	Smartmix Syringes @ 10 g Smartmix Tips Short	REF 212034
	Cartridge @ 52 g Automix Tips Blue 1:1	REF 110293
	essories Smartmix Tips Short	RFF 21 2040
00	Omartinix ripo onort	INLI ZIZO IO

REF 110267

REF 110266

# Vitique



### Adhesive, light-curing luting composite for veneers

Indication: definitive, highly aesthetic cementation of ceramic or composite veneers

### **Product benefits**

- Superior-quality complete system for veneer cementation
- Outstanding shade reproduction and shade selection
- Exact shade match between Try-in paste and cement
- Ultra-easy removal of excess material
- · High shade stability
- Practical user support

### Speciality

Your adhesive partner for Vitique as a supplementary refill pack: Ecosite Bond (Page 37)



### **Cementation System**

- 3 Syringes @ 6 g Esthetic Resin Cement, shades: A2.5, B1, Transparent
- 3 Syringes @ 3.9 g Try-in Paste, shades: A2.5, B1, Transparent
- $1\,$  Syringe @ 2 ml DMG Etching Gel, 1 bottle @ 5 ml TECO
- 1 Bottle @ 3 ml Vitique Silane Adhesive
- 1 Bottle @ 3 ml Vitique Silane Activator

Accessories REF 213080

### Cementation Kit

- $3\;$  Syringes @ 6 g Esthetic Resin Cement, shades: A2.5, B1, Transparent
- 3 Syringes @ 3.9 g Try-in Paste, shades: A2.5, B1, Transparent
- REF 213081 50 Veneer Tips

### Accessory Kit

- $1\,$  Syringe @ 2 ml DMG Etching GeI, 1 bottle @ 5 ml TECO
- 1 Bottle @ 3 ml Vitique Silane Adhesive
- 1 Bottle @ 3 ml Vitique Silane Activator

REF 213082 Accessories

### Refill pack

1	Syringe @ 6 g Esthetic Resin Cement, 10 Veneer Tips Shade A1 Shade A2.5 Shade A4 Shade B1 Shade Bleach Light Shade Transparent Shade White Shade Pink	REF 213038 REF 213039 REF 213040 REF 213041 REF 213042 REF 213043 REF 213044 REF 213045
1	Syringe @ 3.9 g Try-in Paste, 10 Veneer Tips Shade A1 Shade A2.5 Shade A4 Shade B1 Shade Bleach Light Shade Transparent Shade White Shade Pink	REF 213048 REF 213049 REF 213050 REF 213051 REF 213052 REF 213053 REF 213055
Accessories		

100 Veneer Tips REF 213057

# Vitique Silane



# The bonding agent for ceramic and luting composites

Indication: conditioning of siliceous ceramic substrates

### **Product benefits**

- Reliable bond
- 1:1 mixing ratio
- Fresh mix for each application based on the two components



### Pack

- 1 Bottle @ 3 ml Vitique Silane Adhesive
- 1 Bottle @ 3 ml Vitique Silane Activator
- 50 Brushes
- 1 Brush holder





# Direct restoration: restorative

When to use which material

	Constic	Ecosite One	Ecosite Bulk Fill	Ecosite Elements	EcuSphere Shape	EcuSphere Carat	EcuSphere Flow	lonosit Baseliner	LuxaCore Z Dual	DeltaFil
Class I	•	•		•	•	•	•			*
Class II		•		•	•	•				**
Class III	•			•	•	•	•			***
Class IV				•	•	•	•			
Class V	•	•		•	•	•	•			•
Direct inlays/onlays				•		•				
Direct veneers				•						
Deciduous teeth	•	•		•			•			•
Fissure sealing	•	•		•			•			
Liner	•			•			•	•		•
Root caries	•			•			•			
Tunnel preparation	•			•			•			
Minimally invasive	•			•			•			•
Restoration of temporaries	•			•			•			
Core build-up		•	•	•					•	•
Type of material	Composite	Nano hybrid- composite	Nano hybrid- composite	Nano hybrid- composite	Hybrid- composite	Hybrid- composite	Composite	Compomer	Composite	Glass ionomer cement
Number of shades	6	1	3	14	12	4	4	1	3	6
Shade	A1, A2, A3, A3.5, B1, Opaque White		Light, Universal, Contrast	A1, A2, A3, A3.5, B1 Bleach, Lightf, Enamel Bleach, Enamel Light, Enamel Medium, Enamel Dark, A2 Opque, Incisal, White, Brown	A1, A2, A3, A3.5, A4, B1	A1, A2, A3, A3.5	A1, A2, A3, A3.5	B2	Opaque Light, A3, Blue	A1, A2, A3, A4, White, Dark
Polymerisation time	20 s	10 s	20 s	20 s	20 s	20 s	20 s	20 s	20 / 40* s	
Light-curing	•	•	•	•	•	•	•	•	•	
Self-curing									•	•
Radiopaque	•	•	•	•	•	•	•		•	•
Application	2 g Luer-Lock Syringe	A g Syringe 0.25 g Safetip	4 g Syringe 0.25 g Safetip	4 g Syringe 2 g Syringe 0.25 g Safetip	3 g Syringe 0.3 g Safetip	3 g Syringe 0.3 g Safetip	1.5 g Luer-Lock Syringe	0.33 g Saferinge 1.5 g Luer-Lock Syringe	9 g Smartmix Syringe 48 g Automix Cartridge	370 mg, powder/ 75 mg liquid DeltaFil Capsule

<sup>\*</sup> depending on layer thickness

\*\* suitable for areas subject to occlusal load, provided the isthmus measures less than half of the intercuspal area

<sup>\*\*\*</sup> only temporary

# **Prevention and infiltration**

# When to use which material

		Flairesse Paste	Flairesse Gel	Flairesse Foam	Flairesse Varnish	Icon Proximal	lcon Vestibular
Cleaning		•					
Remineralisation			•	•	•		
Fluoridation			•	•	•		
Dentine sensitivity					•		
Caries – intact surface	E1	•	•	•	•	•	•
	E2	0	0	0	•	•	•
	D1					•	•
Aesthetic -	Fluorosis						•
	Cariogenic white spots						•
	ldiopathically induced spots						•
	Trauma-induced spots						•

- Use highly recommended
- O Use recommended under certain conditions

# **Impression**

# When to use which material

	Honigum Pro Light	Honigum Pro Mono	Honigum Pro Heavy	Honigum Pro Putty	Honigum Light	Honigum Mono	Honigum Heavy	Honigum Putty	Silagum Light	Silagum Medium	Silagum Mono	Silagum Putty	StatusBlue
Crown and bridge impression		•	•		•			•	•	•	•	•	
Inlay and onlay impression		•	•		•			•	•	•	•	•	
Implant impression		•	•					•	•	•	•	•	
Functional impression													
Situation impression		•									•		•
Relining impression		•									•		
Impressions for the fabrication of temporaries		•		•		•		•					•
Impressions for the fabrication of vacuum formed templates		•								•	•		•

# Direct restoration: bonding agent

		LuxaBond Total Etch	LuxaBond Universal	Ecosite Bond	Contax	Contax + Activator
Tochnique	Total-etch	•	•	•		
Technique -	Self-etch		•	•	•	•
	Light-curing restorative materials	•	•	•	•	
For adhesive cementation of	Self and dual-curing restorative materials	•	•			•
	Light curing			•	•	•
Curing mode of the adhesive	Dual curing	•	•			
[	Ecosite One					
	Ecosite Elements	•				
	Ecosite Bulk Fill	•	•		•	
-	EcuSphere	•	•		•	
	Ionosit Baseliner	•	•	•	•	
Usage in combination with DMG products	Vitique	•	•	•	•	
	PermaCem	•	•			•
	PermaCem Dual	•	•			•
	PermaCem Universal	•	•			
	LuxaCore Z Dual	•	•			•
	Core build-up	•	•	•		•
	Direct composite fillings	•	•		•	•
Direct restoration	Luting of veneers		•			
	Root canal	•	•			
	Luting of indirect restorations	•	•		•	•

# Impression with Honigum

	Honigum Light		Honigum Mono		Honigum Heavy				Honigum Putty			
		Fast				Fast		Fast		Soft	Soft Fast	Rigid Fast
ml	50	50	380	50	380	380	50	50	380	450×2	450×2	450×2
Mixing time [min]										00:30	00:30	00:30
Working time [min]	02:15	01:30	02:00	01:40	02:15	01:15	01:45	01:15	01:45	02:45	01:45	01:30
Setting time in the mouth [min]	03:30	02:00	03:15	03:00	03:15	02:00	03:15	02:00	03:15	04:15	02:00	01:45
Dimensional change after 24 hrs [%]	≤ 0.5	≤ 0.5	≤ 0.4	≤ 0.4	≤ 0.4	≤ 0.4	≤ 0.4	≤ 0.4	≤ 0.5	≤ 0.5	≤ 0.5	≤ 0.5
Compression set [%]	≤ 0.4	≤ 0.4	≤ 0.9	≤1.0	≤ 0.6	≤ 0.6	≤ 0.6	≤ 0.6	≤ 0.7	≤ 0.5	≤ 0.5	≤ 0.5
Shore-A-hardness 24 h	≈ 54	≈ 56	≈ 52	≈ 51	≈76	≈76	≈72	≈72	≈73	≈ 68	≈ 68	≈ 75
Strain in compression [%]	≈ 4.2	≈ 4.2	2.9-5.5	2.8-5.5	≈ 1.4	≈ 1.4	≈ 1.4	≈ 1.4	≈ 2.1	≈ 2.2	≈ 2.2	≈ 1.7
Shade	•	•	•	•	•	•	•	•	•	•	•	•

# Impression with Honigum Pro

						Honiaum Pro Liaht			Honigum Pro Mono		Honigum Pro Heavy		Honigum Pro Putty	,
	Vi	scosity				Low viscos	sity	Mediur	n viscosity	/ High	n viscosity	/ Н	igh visco	sity
Si	etting be	haviour		Normal-c		•	1		•		•		•	
				Fast-c	uring	•	1		•					
				Monop	hase				•					
	Tec	hnique		Double-m	nixing	•	)				•		•	
				Corre	ection	•			•		•			
:	Syringe n	naterial				•	)							
Inse	ertion resi	istance		Me	Low	•	1		•		•		•	
Technical data	Honigun	n Pro												
	Light		Mono				Heavy				Putty			
													Soft	Rigid
		Fast		Fast		Fast		Fast		Fast		Soft	Fast	Fast
ml	50	50	380	380	50	50	380	380	50	50	380	450×2	450×2	450×2
Mixing time [min]	0: 25	00.00	04.04	00.05	0:0-							00:30	00:30	00:30
Intraoral working time [min]		≤00:30	≤ 01:00	≤00:30		≤00:30		. 04.00	.00.45	. 04.00	- 01-45	- 00 15	. 01. 15	. 04:00
Working time [min]	≤ 02:45		≤ 02:45	≤ 01:30		5 ≤ 01:30 > 2:00		≤ 01:30	≤ 02:45	≤ 01:30	≤ 01:45	≤ 03:15	≤ 01:45	≤ 01:30
Setting time in the mouth [min]  Dimensional change after 24 hrs [%]	≥ 3:15 ≤ 0.35	≥ 2:00 ≤ 0.35	≥ 3:15 ≤ 0.4	≥ 2:00 ≤ 0.4	≥ 3:15 ≤ 0.4	≥ 2:00 ≤ 0.4	≥ 3:15 ≤ 0.4	≥ 2:00 ≤ 0.4	≥ 3:15 ≤ 0.4	≥ 2:00 ≤ 0.4	≥ 4:15 ≤ 0.5	≥ 2:45 ≤ 0.5	≥ 2:00 ≤ 0.5	≥ 1:45 ≤ 0.5
Compression set [%]	≤ 0.35	≤ 0.35	≤ 0.4 ≤ 0.4	≤ 0.4 ≤ 0.4	≤ 0.4 ≤ 0.4	≤ 0.4	≤ 0.4	≤ 0.4 ≤ 0.5	≤ 0.4	≤ 0.4 ≤ 0.6	≤ 0.5	≤ 0.5	≤ 0.5	≤ 0.5
Shore-A-hardness 24 h	≤ 0.30 ≈ 55	≤ 0.30 ≈ 55	≤ 0.4 ≈ 61	≤ 0.4 ≈ 61	≤ 0.4 ≈ 61	≤ 0.5 ≈ 61	≤ 0.5 ≈ 76	≤ 0.5 ≈ 74	≤ 0.5 ≈ 74	≤ 0.0 ≈ 74	≥ 0.7 ≈ 73	≤ 0.5 ≈ 68	≤ 0.5 ≈ 68	≥ 0.5 ≈ 75
Strain in compression [%]	≈ 4.0	≈ 4.0	≈ 2.8	≈ 2.8	≈ 2.8	≈ 2.8	≈ 1.8	≈ 1.8	≈ 2.0	≈ 2.0	≈ 2.1	≈ 2.2	≈ 2.2	≈ 1.7
Shade	•	•	•	•	•		•	•	•	•	•	•	•	•

# Cements: temporary and permanent

		PermaCem 2.0	PermaCem Universal	PermaCem	PermaCem Dual	Vitique	LuxaCore Z Dual	TempoCem NE	ТетроСет	TempoCem ID
Permanent		•	•	•	•	•	•			
Long-term temporary										•
Temporary								•	•	•
	Crowns/bridges	•	•	•				•	•	•
	Inlays/onlays	•	•	•	•			•	•	•
Indications	Veneers		•			•				•
	Root posts	•	•				•			
	Implant abutments	•	•					•		•
	Adhesive		•	•	•	•	•			
Luting	Self-adhesive	•								
	Retentive							•	•	•
	Composite	•	•			•	•			•
Type of material	Compomer			•	•					
	Zinc oxide							•	•	
	Light-curing	•	•		•	•	•			•
Curing mode -	Self-curing	•	•	•	•		•	•	•	•
Number of shades		3	5	1	1	8	3	1	1	1
SIIdues		A2 Universal, A3, Opaque, Transparent	A1, A2.5, Transparent			A1, A2.5, A4, B1, Bleach Light, Transparent, White, Pink	Opaque Light, A3, Blue			

# Cements: temporary and permanent

	PermaCem 2.0	PermaCem Universal	PermaCem	PermaCem Dual	Vitique	TempoCemNE	ТетроСет	TempoCem ID
Inlays	•	0	0	0		•	•	•
Onlays	•	0	0	0		0	0	•
Crowns/bridges	•	0	0	0		•	•	•
Root posts	•							
On implant abutments	•						0	•
Inlays	0	0	0	0		•	•	•
Onlays	0	0	0	0		0	0	•
Crowns/bridges	•	0	0	0		•	•	•
Veneers		0	0	0	•			•
Crowns/bridges	•	0	0	0		•	•	•
Root posts	•							
On implant abutments	•						0	•
Inlavs	0	0	0	0		•		•
						0	0	
Crowns	•					•	•	•
Root posts	•							
Veneers		0	0	0	•			•
	Onlays Crowns/bridges Root posts On implant abutments  Inlays Onlays Crowns/bridges Veneers  Crowns/bridges Root posts On implant abutments  Inlays Onlays Crowns Root posts	Inlays Onlays Onlays Crowns/bridges Root posts On implant abutments Inlays Onlays Crowns/bridges Veneers  Crowns/bridges Root posts On implant abutments Inlays Onlays Crowns Onlays	Inlays Onlays Onlays Onlays Onlays On implant abutments On implant abutments Onlays On	Inlays Onlays Onlays Onlays Onlays On implant abutments Inlays Onlays On	Inlays	Inlays	Inlays Onlays Onlays Onlays On implant abutments  Inlays On implant abutments  On implant abutments  Inlays Onlays Onlays Onlays Onlays Onlays On implant abutments  Inlays Onlays Onlay	Inlays

Use highly recommended

O Use recommended under certain conditions


# Prosthetic restoration: temporary and long-term stable

	Luxatemp MaxProtect	Luxatemp Star	Luxatemp Fluorescence	Luxatemp Plus	LuxaFlow Star	Luxatemp Glaze & Bond
Temporary crowns and bridges	•	•	•	•		
Temporary inlays and onlays	•	•	•	•		
Temporary partial crowns	•	•	•	•		
Temporary veneers	•	•	•	•		
Long-term temporaries	•	•	•	•		
For long-term stable crowns and bridges	•					
Complex long-term temporary restorations	•					
Enhancement and repair of bis-acrylate temporaries					•	
Small fillings (cavity class III, IV, V)					•	
Minimally invasive fillings					•	
Underfillings					•	
Surface sealer						•
Repair of temporary crowns and bridges					•	•
Type of material	Bis-acrylate	Bis-acrylate	Bis-acrylate	Bis-acrylate	Composite	1-component material, based on multifunctional methacrylates

## **Automix Dispenser**

## **Automix Tips**

#### Dispenser for 25 ml cartridges, mixing ratio 1:1



#### Dispenser for 50 ml cartridges, mixing ratio 1:1



#### Dispenser for 50 ml cartridges, mixing ratio 10:1 / 4:1

Type 50 10:1 REF 110411	Suitable for Luxatemp Plus Luxatemp MaxProtect Luxatemp Star Luxatemp Fluorescence LuxaPick-up	Automix Tips Blue 10:1 45 × REF 110409
	Suitable for LuxaBite	Automix Tips Yellow 10:1 25 × REF 110609
	Suitable for StatusBlue	Automix Tips Pink 4:1 25 × REF 909797

\*Not suitable for Silagum Comfort

## **Application Tips**

## For syringes and Smartmix

#### Smartmix Tips Green-Blue

REF 110990 (50 pcs)



REF 212040 (50 pcs)



Suitable for Luxatemp Star Luxatemp Fluorescence



**Smartmix Tips Endo** 

REF 212041 (50 pcs each)

Suitable for
TempoCem NE
TempoCem
PermaCem Universal
PermaCem 2.0
PermaCem Dual
PermaCem

#### Smartmix Tips Combi

REF 212042 (50 pcs each)





Suitable for LuxaCore Z Dual

#### **Smartmix Tips Green**

REF 110724 (50 pcs)



Suitable for TempoCem ID

#### SmartmixTips Endo L

REF 213083 (20 pcs)



Smartmix Tips Endo M

REF 213373 (20 pcs)

Suitable for PermaCem 2.0



Suitable for LuxaCore Z Dual

#### Veneer Tips

REF 213057 (100 pcs)



REF 109942 (50 pcs)



Suitable for Vitique



Suitable for DMG Etching Gel

#### Luer-Lock Tips Black

REF 211477 (60 pcs)



REF 220712 (20 pcs)



Suitable for LuxaFlow Star Ionosit Baseliner EcuSphere Flow



Suitable for Constic Ecosite Elements Highlight

# MixStar Tips

## **Automix Tips**

#### MixStar eMotion



Suitable for Honigum Pro Mono Honigum Pro Heavy Honigum Pro Heavy Fast Honigum Pro Putty Honigum Putty

MixStar Tips Blue 50 x REF 903334

Suitable for Honigum Mono Honigum Heavy Honigum Heavy Fast Silagum Putty Soft StatusBlue





## DMG symbols for the type of packaging



## Index

Flairesse Bleaching Gel CP 5% 23	LuxaPrint Ortho Plus 20
Flairesse Bleaching Gel CP 10% 23	LuxaPrint Tray 21
Flairesse Bleaching Gel CP 16% 23	Luxatemp Fluorescence 56
Flairesse Gel/Foam 24	Luxatemp Glaze & Bond 57
Flairesse Paste 24	Luxatemp MaxProtect 55
Flairesse Varnish 25	Luxatemp Plus 56
Honigum Heavy 46	Luxatemp Star 55
Honigum Light 45	MixStar eMotion 50
Honigum Mono 45	MixStar Tips 79
Honigum Pro Heavy 44	O-Bite 51
Honigum Pro Light 43	O-Bite Scan 51
Honigum Pro Mono 43	PermaCem 63
Honigum Pro Putty 44	PermaCem 2.0 62
Honigum Putty 46	PermaCem Dual 63
Icon Proximal 26	PermaCem Universal 62
Icon Vestibular 26	Safetip Dispenser 40
Ionosit Baseliner 35	Safetip Dispenser Pro 40
LuxaBite 52	Silagum Comfort 53
LuxaBond Total Etch 36	Silagum Light 47
LuxaBond Universal 36	Silagum Medium 47
LuxaCore Z Dual 38	Silagum Mono 48
LuxaCore Z Post System 38	Silagum Putty 48
LuxaFlow Star 57	Silagum Putty Soft 49
LuxaPick-up 53	StatusBlue 49
LuxaPost 39	TempoCem 59
LuxaPrint Model 21	TempoCem ID 58
LuxaPrint Ortho 20	TempoCem NE 58
LuxaPrint Ortho Comfort 19	Vitique 64
LuxaPrint Ortho Flex 19	Vitique Silane 65
	Flairesse Bleaching Gel CP 10% 23 Flairesse Bleaching Gel CP 16% 23 Flairesse Gel/Foam 24 Flairesse Paste 24 Flairesse Varnish 25 Honigum Heavy 46 Honigum Light 45 Honigum Pro Heavy 44 Honigum Pro Light 43 Honigum Pro Mono 43 Honigum Pro Putty 44 Honigum Pro Putty 46 Icon Proximal 26 Icon Vestibular 26 Ionosit Baseliner 35 LuxaBite 52 LuxaBond Total Etch 36 LuxaCore Z Dual 38 LuxaCore Z Post System 38 LuxaFlow Star 57 LuxaPrint Model 21 LuxaPrint Ortho 20 LuxaPrint Ortho Comfort 19





#### Chemisch-Pharmazeutische Fabrik GmbH

Elbgaustraße 248 22547 Hamburg Germany
Fon: +49. (0) 40. 84 006-0 Fax: +49. (0) 40. 84 006-222
info@dmg-dental.com www.dmg-dental.com
www.facebook.com/dmgdental

