



Looks good – because
you can't see it.

TempoCem ID
Temporary cements



Visible to no-one. Apart from you.

TempoCem ID ensures outstanding results thanks to its visual adaptability. The highly aesthetic temporary cement is invisible when it needs to be, yet is also detectable for dentists when required.

Invisible below the restoration:

Thanks to its optimum transparency, the colour-neutral cement does not affect the shade of the provisional restoration. This makes TempoCem ID the perfect choice for both posterior and anterior restorations.



Benefits:

- Excellent aesthetics thanks to optimum transparency
- Invisible below the restoration
- Yet easy to detect for excess removal
- Optimal bonding properties, easy to remove
- Eugenol-free
- Extremely easy to use
- Outstanding flow and mixing properties
- Perfect for long-lasting provisional restorations



Easy to detect for easy excess removal:

Unlike other aesthetic temporary cements, excess TempoCem ID is easy to detect and can often be removed in one piece. The eugenol-free TempoCem ID has excellent flow and mixing properties. The material is perfect for long-lasting provisional restorations.

See for yourself how TempoCem ID makes aesthetic dreams come true.

Indications:

Temporary and long-term provisional luting of provisional crowns, bridges, inlays, onlays and veneers as well as implant-supported provisional restorations.

PRODUCT OVERVIEW

TempoCem ID

Packaging:

1 Smartmix Syringe @ 5 ml
20 Smartmix Tips

REF 213200



More about
TempoCem ID





DMG

Chemisch-Pharmazeutische Fabrik GmbH

Elbgaustraße 248 22547 Hamburg Germany

Fon: +49, (0) 40. 84 006-0 Fax: +49, (0) 40. 84 006-222

info@dmg-dental.com www.dmg-dental.com

www.facebook.com/dmgdental

