



Ecosite One

**Simply matches.  
Always.**

 **DMG**

# Ecosite One

## Posterior restorations made simple

The innovative Ecosite One layering composite sets new standards in the field of posterior restorations. The Push-and-Flow effect and the One-Shade concept produce outstanding aesthetics and impress all those who opt for universal composite layering technology and want to save valuable time in the process.

### An easy process for dentists, with great results for patients:

- Push-and-Flow effect for excellent handling
- One colour suits all, thanks to the One-Shade concept
- Valuable time savings thanks to 3 mm curing depth and 10 second rapid polymerisation (depending on light intensity\*)



More information on Ecosite One:  
[qr.dmg.dental/ecositeone](http://qr.dmg.dental/ecositeone)

### Ecosite. One collection for everything:

Ecosite One, Ecosite Elements, Ecosite Bulk Fill and Ecosite Bond