



DMG LuxaPrint: Validierte Systeme für unsere hochwertigen 3D-Druck-Resins

Sie haben ein 3D-Drucksystem von DMG oder eines anderen Herstellers und möchten wissen, ob unsere Resins dafür validiert sind? Werfen Sie einfach einen Blick auf die untenstehende Tabelle.

	Drucker Modell	Nachbelichtung	Material								
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRY	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN
3Demax	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Otoflash	✓	✓	✓	✓	✓	✓	✓			✓
3Delite	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Otoflash	✓	✓	✓	✓	✓	✓	✓			✓
DMG	3Decure, RS, Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Desk RS-5	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DentaMile Lab5/Lab5 Pro	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Otoflash	✓	✓	✓	✓	✓	✓	✓			✓
DentaMile Desk RS-1	3Decure, RS, Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Desk MC-5	DentaMile Cure MC	✓	✓	✓	✓	✓	✓	✓	✓	✓
	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rapid Shape	RS D20+	CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓			✓

	Drucker Modell	Nachbelichtung	Material								
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRY	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN
Rapid Shape	RS D10+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
Rapid Shape	PRO 20	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rapid Shape	RS D30+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
Rapid Shape	RS D50+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓		✓
		RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rapid Shape	RS one	CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓		✓
		P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
Straumann	P20+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
		P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
Straumann	P10+	DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓		✓		✓
		P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		PROCURE, DentaMile Cure RS, CURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
Straumann	PRO20	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
		P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
Burms	P50+	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
		P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓		✓
Asiga	Max 2	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
Asiga	Max X	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
Asiga	Ultra	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
Asiga	Pro 4k	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
Asiga	Pro 4k XL	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ackuretta	SOL 405	Curie, Curie Plus					✓	✓		✓	✓
Organical	3D Print X1 NS	Wash&Cure					✓	✓		✓	✓
Shining	AccuFab L4D	Fab Cure 2	✓				✓	✓	✓	✓	✓
Huvitz	Lilivis Print	Lilivis Cure						✓			✓
Rayshape	Edge E2	ShapeCure D									✓
UNIZ	UBEE	UCure									✓