



DMG LuxaPrint: Validierte Systeme für unsere hochwertigen 3D-Druck-Resins

Sie haben ein 3D-Drucksystem von DMG oder eines anderen Herstellers und möchten wissen, ob unsere Resins dafür validiert sind? Werfen Sie einfach einen Blick auf die untenstehende Tabelle.

	Drucker Modell	Nachbelichtung	Material								
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRV	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN
DMG	3Demax	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	3Delite	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	DentaMile Lab5/ Lab5 Pro	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓			✓
DentaMile Desk MC-5	DentaMile Cure MC	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Rapid Shape	RS D20+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	RS D10+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	RS D30+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
RS D50+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Otoflash	✓	✓	✓	✓	✓	✓			✓	
P20+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Otoflash	✓	✓	✓	✓		✓			✓	
Straumann	P10+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	P50+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓			✓



Drucker Modell		Nachbelichtung	Material								
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRY	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN
Asiga	MAX UV	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
	PRO 4K	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ackuretta	SOL 405	Curie					✓	✓		✓	✓
Organical	3D Print X1 NS	Wash&Cure									✓
Shining	AccuFab L4D	Fab Cure 2	✓				✓	✓	✓	✓	
	AccuFab D1s	Fab Cure 2	✓	✓			✓	✓	✓	✓	
Huvis	Lilivis Print	Lilivis Cure						✓			