



## DMG LuxaPrint: validated systems for our high-quality 3D printing resins

Are you using a DMG 3D printing system or one from a different manufacturer and would like to know whether our resins have been validated for it? Have a quick look at the table below.

	Printer model	Post-curing	Material								
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRV	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN
DMG	3Demax	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	3Delite	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	DentaMile Desk RS-5	3Decure, RS, Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
DentaMile Lab5/ Lab5 Pro	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Otoflash	✓	✓	✓	✓	✓	✓			✓	
DentaMile Desk RS-1	3Decure, RS, Cure, P Cure	✓	✓	✓	✓	✓	✓	✓		✓	
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓		✓	
DentaMile Desk MC-5	DentaMile Cure MC	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Rapid Shape	RS D20+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓



	Printer model	Post-curing	Material								
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRY	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN
Rapid Shape	RS D10+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	PRO 20	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
	RS D30+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
RS D50+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Otoflash	✓	✓	✓	✓	✓	✓			✓	
RS one	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓		✓	
	CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓		✓	
Straumann	P20+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	P10+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	PRO20	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		PROCURE, DentaMile Cure RS, CURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
	P50+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
Otoflash		✓	✓	✓	✓	✓	✓			✓	
Burms	MiiCraft Alpha	MiiCure	✓	✓	✓	✓	✓	✓	✓	✓	
Asiga	Max 2	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	
	Max X	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	
	Ultra	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	
	Pro 4k	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	
	Pro 4k XL	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	
Ackuretta	SOL 405	Curie, Curie Plus					✓	✓		✓	
Organical	3D Print X1 NS	Wash&Cure					✓	✓		✓	
Shining	AccuFab L4D	Fab Cure 2	✓				✓	✓	✓	✓	
	AccuFab D1s	Fab Cure 2	✓	✓			✓	✓	✓	✓	
Huvitz	Lilivis Print	Lilivis Cure					✓			✓	
Rayshape	Edge E2	ShapeCure D								✓	
UNIZ	UBEE	UCure								✓	