



## DMG LuxaPrint: validated systems for our high-quality 3D printing resins

Are you using a DMG 3D printing system or one from a different manufacturer and would like to know whether our resins have been validated for it? Have a quick look at the table below.

	Printer model	Post-curing	Material								
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRV	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN
DMG	3Demax	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	3Delite	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	DentaMile Lab5/ Lab5 Pro	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓			✓
DentaMile Desk MC-5	DentaMile Cure MC	✓	✓	✓	✓	✓	✓	✓	✓	✓	
RapidShape	RS D20+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	RS D10+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	RS D30+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
RS D50+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Otoflash	✓	✓	✓	✓	✓	✓	✓		✓	
Straumann	P20+	P cure 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	P10+	P cure 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	P50+	P cure 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓		✓
Asiga	MAX UV	Otoflash	✓	✓		✓	✓	✓	✓	✓	✓
	PRO 4K	Otoflash	✓	✓		✓	✓	✓	✓	✓	✓
Ackuretta	SOL 405	Curie					✓	✓		✓	✓
Organical	3D Print X1 NS	Wash&Cure									✓
Shining	AccuFab L4D	Shining Cure						✓	✓	✓	
	AccuFab D1s	Shining Cure						✓	✓		