



DMG LuxaPrint: 经过验证的系统,适用于我们优质的 3D 打印树脂

您是否正在使用 DMG 3D 打印系统或其他制造商的系统,并想知道我们的树脂是否已针对该系统完成验证?请查看下表。

打印机型号	后固化	材料									
		LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRY	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN	
DMG	3Demax	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	3Delite	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	DentaMile Lab5/ Lab5 Pro	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓			✓
DentaMile Desk MC-5	DentaMile Cure MC	✓	✓	✓	✓	✓	✓	✓	✓	✓	
RapidShape	RS D20+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	RS D10+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	RS D30+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
RS D50+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Otoflash	✓	✓	✓	✓	✓	✓			✓	
Straumann	P20+	P cure 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	P10+	P cure 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	P50+	P cure 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓			✓
Asiga	MAX UV	Otoflash	✓	✓		✓	✓	✓	✓	✓	✓
	PRO 4K	Otoflash	✓	✓		✓	✓	✓	✓	✓	✓
Ackuretta	SOL 405	Curie				✓	✓			✓	✓
Organical	3D Print X1 NS	Wash&Cure									✓
Shining	AccuFab L4D	Shining Cure						✓	✓	✓	
	AccuFab D1s	Shining Cure						✓	✓		