



DMG LuxaPrint: 经过验证的系统,适用于我们优质的 3D 打印树脂

您是否正在使用 DMG 3D 打印系统或其他制造商的系统,并想知道我们的树脂是否已针对该系统完成验证?请查看下表。

	打印机型号	后固化	材料								
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRY	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN
3Demax	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Otoflash	✓	✓	✓	✓	✓	✓	✓			✓
3Delite	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Otoflash	✓	✓	✓	✓	✓	✓	✓			✓
DMG	3Decure, RS, Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Desk RS-5	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DentaMile Lab5/Lab5 Pro	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DentaMile Desk RS-1	Otoflash	✓	✓	✓	✓	✓	✓	✓			✓
	3Decure, RS, Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓		✓
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓		✓
	DentaMile Desk MC-5	DentaMile Cure MC	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rapid Shape	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Otoflash	✓	✓	✓	✓	✓	✓	✓			✓

打印机型号		后固化	材料								
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRY	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN
Rapid Shape	RS D10+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
PRO 20	PRO 20	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
RS D30+	RS D30+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
RS D50+	RS D50+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓			✓
		RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
RS one	RS one	CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
Straumann	P20+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
		P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
Straumann	P10+	DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓		✓		✓
		P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		PROCURE, DentaMile Cure RS, CURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
Straumann	PRO20	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
		P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
Burms	P50+	Otoflash	✓	✓	✓	✓	✓	✓			✓
		MiiCraft Alpha	MiiCure	✓	✓	✓	✓	✓	✓	✓	✓
		Max 2	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓
		Max X	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓
Asiga	Ultra	Ultra	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓
		Pro 4k	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓
		Pro 4k XL	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓
		Ackuretta	SOL 405	Curie, Curie Plus				✓	✓		✓
Organical	3D Print X1 NS	Wash&Cure					✓	✓		✓	✓
		Shining	AccuFab L4D	Fab Cure 2	✓			✓	✓	✓	✓
Huvitz	AccuFab D1s	Fab Cure 2	✓	✓			✓	✓	✓	✓	✓
		Lilivis Print	Lilivis Cure					✓			✓
Rayshape	Edge E2	ShapeCure D									✓
		UNIZ	UBEE	UCure							✓