



## DMG LuxaPrint: 经过验证的系统, 适用于我们优质的 3D 打印树脂

您是否正在使用 DMG 3D 打印系统或其他制造商的系统, 并想知道我们的树脂是否已针对该系统完成验证? 请查看下表。

| 打印机型号               | 后固化                         | 材料                       |                      |                         |                      |                     |                     |                     |                     |                    |   |
|---------------------|-----------------------------|--------------------------|----------------------|-------------------------|----------------------|---------------------|---------------------|---------------------|---------------------|--------------------|---|
|                     |                             | LuxaPrint Ortho          | LuxaPrint Ortho Plus | LuxaPrint Ortho Comfort | LuxaPrint Ortho Flex | LuxaPrint Model BGE | LuxaPrint Model GRY | LuxaPrint Model IVR | LuxaPrint Model TRA | LuxaPrint Tray GRN |   |
| DMG                 | 3Demax                      | 3Decure, RS Cure, P Cure | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  | ✓ |
|                     |                             | DentaMile Cure RS-XL     | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  | ✓ |
|                     |                             | Otoflash                 | ✓                    | ✓                       | ✓                    | ✓                   |                     | ✓                   |                     |                    | ✓ |
|                     | 3Delite                     | 3Decure, RS Cure, P Cure | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  | ✓ |
|                     |                             | DentaMile Cure RS-XL     | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  | ✓ |
|                     |                             | Otoflash                 | ✓                    | ✓                       | ✓                    | ✓                   |                     | ✓                   |                     |                    | ✓ |
|                     | DentaMile Lab5/<br>Lab5 Pro | 3Decure, RS Cure, P Cure | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  | ✓ |
|                     |                             | DentaMile Cure RS-XL     | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  | ✓ |
|                     |                             | Otoflash                 | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   |                     |                    | ✓ |
| DentaMile Desk MC-5 | DentaMile Cure MC           | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
| RS D20+             | RS Cure, 3Decure, P Cure    | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | DentaMile Cure RS-XL        | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | Otoflash                    | ✓                        | ✓                    | ✓                       | ✓                    |                     | ✓                   |                     |                     | ✓                  |   |
| RS D10+             | RS Cure, 3Decure, P Cure    | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | DentaMile Cure RS-XL        | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | Otoflash                    | ✓                        | ✓                    | ✓                       | ✓                    |                     | ✓                   |                     |                     | ✓                  |   |
| Rapid Shape         | RS Cure, 3Decure, P Cure    | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | DentaMile Cure RS-XL        | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | Otoflash                    | ✓                        | ✓                    | ✓                       | ✓                    |                     | ✓                   |                     |                     | ✓                  |   |
| RS D30+             | RS Cure, 3Decure, P Cure    | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | DentaMile Cure RS-XL        | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | Otoflash                    | ✓                        | ✓                    | ✓                       | ✓                    |                     | ✓                   |                     |                     | ✓                  |   |
| RS D50+             | RS Cure, 3Decure, P Cure    | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | DentaMile Cure RS-XL        | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | Otoflash                    | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   |                     |                     | ✓                  |   |
| P20+                | P Cure, 3Decure, RS Cure    | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | DentaMile Cure RS-XL        | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | Otoflash                    | ✓                        | ✓                    | ✓                       | ✓                    |                     | ✓                   |                     |                     | ✓                  |   |
| Straumann           | P Cure, 3Decure, RS Cure    | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | DentaMile Cure RS-XL        | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | Otoflash                    | ✓                        | ✓                    | ✓                       | ✓                    |                     | ✓                   |                     |                     | ✓                  |   |
| P50+                | P Cure, 3Decure, RS Cure    | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | DentaMile Cure RS-XL        | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |   |
|                     | Otoflash                    | ✓                        | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   |                     |                     | ✓                  |   |

| 打印机型号     |                | 后固化          | 材料              |                      |                         |                      |                     |                     |                     |                     |                    |
|-----------|----------------|--------------|-----------------|----------------------|-------------------------|----------------------|---------------------|---------------------|---------------------|---------------------|--------------------|
|           |                |              | LuxaPrint Ortho | LuxaPrint Ortho Plus | LuxaPrint Ortho Comfort | LuxaPrint Ortho Flex | LuxaPrint Model BGE | LuxaPrint Model GRY | LuxaPrint Model IVR | LuxaPrint Model TRA | LuxaPrint Tray GRN |
| Asiga     | MAX UV         | Otoflash     | ✓               | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |
|           | PRO 4K         | Otoflash     | ✓               | ✓                    | ✓                       | ✓                    | ✓                   | ✓                   | ✓                   | ✓                   | ✓                  |
| Ackuretta | SOL 405        | Curie        |                 |                      |                         |                      | ✓                   | ✓                   |                     | ✓                   | ✓                  |
| Organical | 3D Print X1 NS | Wash&Cure    |                 |                      |                         |                      |                     |                     |                     |                     | ✓                  |
| Shining   | AccuFab L4D    | Fab Cure 2   | ✓               |                      |                         |                      | ✓                   | ✓                   | ✓                   | ✓                   |                    |
|           | AccuFab D1s    | Fab Cure 2   | ✓               | ✓                    |                         |                      | ✓                   | ✓                   | ✓                   | ✓                   |                    |
| Huvitz    | Lilivis Print  | Lilivis Cure |                 |                      |                         |                      |                     | ✓                   |                     |                     |                    |