



## DMG LuxaPrint: sistemas validados para nossas resinas de alta qualidade para impressão 3D

Você está usando um sistema de impressão 3D da DMG ou de outro fabricante e gostaria de saber se nossas resinas foram validadas para ele? Dê uma olhada rápida na tabela abaixo.

	Modelo da impressora	Pós-polimerização	Material								
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRY	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN
DMG	3Demax	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	3Delite	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	DentaMile Desk RS-5	3Decure, RS, Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
DentaMile Lab5/ Lab5 Pro	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Otoflash	✓	✓	✓	✓	✓	✓			✓	
DentaMile Desk RS-1	3Decure, RS, Cure, P Cure	✓	✓	✓	✓	✓	✓	✓		✓	
	DentaMile Cure RS, CURE, PROCURE	✓	✓	✓	✓	✓	✓	✓		✓	
DentaMile Desk MC-5	DentaMile Cure MC	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Rapid Shape	RS D20+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓



	Modelo da impressora	Pós-polimerização	Material								
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRY	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN
Rapid Shape	RS D10+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
	PRO 20	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
	RS D30+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓		✓			✓
RS D50+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	RS cure XL, DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Otoflash	✓	✓	✓	✓	✓	✓			✓	
RS one	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓		✓	
	CURE, DentaMile Cure RS, PROCURE	✓	✓	✓	✓	✓	✓	✓		✓	
P20+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Otoflash	✓	✓	✓	✓		✓			✓	
P10+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Otoflash	✓	✓	✓	✓		✓			✓	
Straumann	PRO20	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
		PROCURE, DentaMile Cure RS, CURE	✓	✓	✓	✓	✓	✓	✓	✓	✓
		DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
P50+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	PROCURE, DentaMile Cure RS, CURE	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	DentaMile Cure RS-XL, RS cure XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Otoflash	✓	✓	✓	✓	✓	✓			✓	
Burms	MiiCraft Alpha	MiiCure	✓	✓	✓	✓	✓	✓	✓	✓	
Asiga	Max 2	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	
	Max X	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	
	Ultra	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	
	Pro 4k	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	
	Pro 4k XL	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	
Ackuretta	SOL 405	Curie, Curie Plus					✓	✓	✓	✓	
Organical	3D Print X1 NS	Wash&Cure					✓	✓	✓	✓	
Shining	AccuFab L4D	Fab Cure 2	✓				✓	✓	✓	✓	
	AccuFab D1s	Fab Cure 2	✓	✓			✓	✓	✓	✓	
Huvitz	Lilivis Print	Lilivis Cure					✓			✓	
Rayshape	Edge E2	ShapeCure D								✓	
UNIZ	UBEE	UCure								✓	



# LuxaPrint NX



## DMG LuxaPrint NX: sistemas validados para nossas resinas de alta qualidade para impressão 3D

Você está usando um sistema de impressão 3D da DMG ou de outro fabricante e gostaria de saber se nossas resinas foram validadas para ele? Dê uma olhada rápida na tabela abaixo.

Modelo da impressora	Pós-polimerização	Material			
		LuxaPrint Ortho NX	LuxaPrint Ortho Plus NX	LuxaPrint Ortho Comfort NX	LuxaPrint Ortho Flex NX
3Demax	3Decure, RS Cure, P Cure			✓	✓
	DentaMile Cure RS, CURE, PROCURE			✓	✓
	DentaMile Cure RS-XL, RS cure XL			✓	✓
	Otoflash				✓
3Delite	3Decure, RS Cure, P Cure			✓	✓
	DentaMile Cure RS, CURE, PROCURE			✓	✓
	DentaMile Cure RS-XL, RS cure XL			✓	✓
	Otoflash				✓
DMG DentaMile Desk RS-5	3Decure, RS, Cure, P Cure			✓	✓
	DentaMile Cure RS, CURE, PROCURE			✓	✓
	DentaMile Cure RS-XL, RS cure XL			✓	✓
	Otoflash				✓
DMG DentaMile Lab5/ Lab5 Pro	3Decure, RS Cure, P Cure			✓	✓
	DentaMile Cure RS, CURE, PROCURE			✓	✓
	DentaMile Cure RS-XL, RS cure XL			✓	✓
	Otoflash				✓
DMG DentaMile Desk RS-1	3Decure, RS, Cure, P Cure			✓	
	DentaMile Cure RS, CURE, PROCURE			✓	
DMG DentaMile Desk MC-5	DentaMile Cure MC	✓	✓	✓	✓
Rapid Shape RS D20+	RS Cure, 3Decure, P Cure			✓	✓
	CURE, DentaMile Cure RS, PROCURE			✓	✓
	RS cure XL, DentaMile Cure RS-XL			✓	✓
	Otoflash				✓



Modelo da impressora		Pós-polimerização	Material			
			LuxaPrint Ortho NX	LuxaPrint Ortho Plus NX	LuxaPrint Ortho Comfort NX	LuxaPrint Ortho Flex NX
Rapid Shape	RS D10+	RS Cure, 3Decure, P Cure			✓	✓
		CURE, DentaMile Cure RS, PROCURE			✓	✓
		RS cure XL, DentaMile Cure RS-XL			✓	✓
		Otoflash				✓
	PRO 20	RS Cure, 3Decure, P Cure			✓	✓
		CURE, DentaMile Cure RS, PROCURE			✓	✓
		RS cure XL, DentaMile Cure RS-XL			✓	✓
		Otoflash				✓
	RS D30+	RS Cure, 3Decure, P Cure			✓	✓
		CURE, DentaMile Cure RS, PROCURE			✓	✓
		RS cure XL, DentaMile Cure RS-XL			✓	✓
		Otoflash				✓
	RS D50+	RS Cure, 3Decure, P Cure			✓	✓
		CURE, DentaMile Cure RS, PROCURE			✓	✓
		RS cure XL, DentaMile Cure RS-XL			✓	✓
		Otoflash				✓
RS one	RS Cure, 3Decure, P Cure			✓		
	CURE, DentaMile Cure RS, PROCURE			✓		
Straumann	P20+	P Cure, 3Decure, RS Cure			✓	✓
		DentaMile Cure RS-XL			✓	✓
		Otoflash				✓
	P10+	P Cure, 3Decure, RS Cure			✓	✓
		DentaMile Cure RS-XL, RS cure XL			✓	✓
		Otoflash				✓
	PRO20	P Cure, 3Decure, RS Cure			✓	✓
		PROCURE, DentaMile Cure RS, CURE			✓	✓
		DentaMile Cure RS-XL, RS cure XL			✓	✓
		Otoflash				✓
	P50+	P Cure, 3Decure, RS Cure			✓	✓
		PROCURE, DentaMile Cure RS, CURE			✓	✓
DentaMile Cure RS-XL, RS cure XL				✓	✓	
Otoflash					✓	
Burms	MiiCraft Alpha	MiiCure	✓	✓	✓	✓
Asiga	Max 2	Otoflash			✓	✓
	Max X	Otoflash			✓	✓
	Ultra	Otoflash			✓	✓
	Pro 4k	Otoflash			✓	✓
	Pro 4k XL	Otoflash			✓	✓