



DMG LuxaPrint: sistemas validados para nossas resinas de alta qualidade para impressão 3D

Você está usando um sistema de impressão 3D da DMG ou de outro fabricante e gostaria de saber se nossas resinas foram validadas para ele? Dê uma olhada rápida na tabela abaixo.

	Modelo da impressora	Pós-polimerização	Material									
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRY	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN	
DMG	3Demax	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		Otoflash	✓	✓	✓	✓		✓			✓	
	3Delite	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		Otoflash	✓	✓	✓	✓		✓			✓	
	DentaMile Lab5/Lab5 Pro	3Decure, RS Cure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		Otoflash	✓	✓	✓	✓	✓	✓			✓	
	DentaMile Desk MC-5	DentaMile Cure MC	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Rapid Shape	RS D20+	RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓
			DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓
Otoflash			✓	✓	✓	✓		✓			✓	
RS D10+		RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		Otoflash	✓	✓	✓	✓		✓			✓	
RS D30+		RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		Otoflash	✓	✓	✓	✓		✓			✓	
RS D50+		RS Cure, 3Decure, P Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		Otoflash	✓	✓	✓	✓	✓	✓			✓	
Straumann	P20+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		Otoflash	✓	✓	✓	✓		✓			✓	
	P10+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		Otoflash	✓	✓	✓	✓		✓			✓	
P50+	P Cure, 3Decure, RS Cure	✓	✓	✓	✓	✓	✓	✓	✓	✓		
	DentaMile Cure RS-XL	✓	✓	✓	✓	✓	✓	✓	✓	✓		
	Otoflash	✓	✓	✓	✓	✓	✓			✓		

	Modelo da impressora	Pós-polimerização	Material								
			LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Comfort	LuxaPrint Ortho Flex	LuxaPrint Model BGE	LuxaPrint Model GRY	LuxaPrint Model IVR	LuxaPrint Model TRA	LuxaPrint Tray GRN
Asiga	MAX UV	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
	PRO 4K	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ackuretta	SOL 405	Curie					✓	✓		✓	✓
Organical	3D Print X1 NS	Wash&Cure									✓
Shining	AccuFab L4D	Fab Cure 2	✓				✓	✓	✓	✓	✓
	AccuFab D1s	Fab Cure 2	✓	✓			✓	✓	✓	✓	✓
Huvitz	Lilivis Print	Lilivis Cure						✓			✓