Perfect partners for the best bonding.
PermaCem Universal & LuxaBond Universal

- A reliable combination with very good bonding properties on all common materials
- A wide range of indications
- Simple and straightforward system pairing
- Familiar working procedure; low risk of error



Product overview

PermaCem Universal

Packaging: 1 Smartmix Syringe @ 9 g, 10 Smartmix Tips Short

Shade A1	REF 213401	Shade A2.5	REF 213402	Shade Transparent	REF 213404

Accessories:

50 Smartmix Tips Short	REF 212040

LuxaBond Universal

Packaging:

50 Single-use Brushes White, 25 Endo Brushes, 1 Mixing Palette 1 Bottle @ 5 ml Bond A, 1 Bottle @ 4.5 ml Bond B	REF 213379

More about
PermaCem Universal

2026-09

DMG
Chemisch-Pharmazeutische Fabrik GmbH
Elbgaustraße 248 22547 Hamburg Germany
Fon: +49 (0) 40-84 006-0 Fax: +49 (0) 40-84 006-222
info@dmg-dental.com www.dmg-dental.com
www.facebook.com/dmgdental

DMG

A smile ahead
together



Universal **brilliance.**

PermaCem Universal
permanent cements



Multi-purpose
For ceramics, metals, metal alloys and composites

Strong adhesive force
Clinically proven, scientifically confirmed

DMG



Perfect partners for the best bonding.

PermaCem Universal & LuxaBond Universal.

The sheer diversity of modern restorative treatment options, with the high demands they place on a wide range of restorative materials and processing methods, calls for effective and reliable zero-compromise luting systems.

The combination of PermaCem Universal and LuxaBond Universal delivers excellent adhesion across the board. An in-vitro study conducted by the University of Giessen* on adhesively cemented composite crowns on tooth stumps demonstrates just how well the two products complement one another.

* "Marginal integrity of modern CAD/CAM materials under alternating stress in-vitro." B. Wöstmann, M. Schlenz, T. Niem, 2018

PermaCem Universal

Universal brilliance.

PermaCem Universal is your adhesive all-rounder for secure luting that is suitable for universal use. This dual-curing composite luting cement provides reliably strong adhesion for the versatile cementation of indirect restorations made from ceramics, metals, metal alloys and composites. Its high abrasion resistance ensures reliable marginal integrity and supports the longevity of the restoration. Minimal water absorption also guarantees high aesthetics and effective control of expansion forces.

An easy and comfortable procedure.
A secure fit with just one cement and one bond for the most common indications and restorations – that's what makes the work easy. The fact that excess PermaCem Universal material is especially easy to remove makes the straightforward procedure even simpler. After initial polymerisation, the excess material can be removed with extreme ease in one piece with the scaler.


Good reasons to choose PermaCem Universal

Strong adhesive force
Reliable across the board

Multi-purpose
For ceramics, metals and composites

CAD/CAM-compatible
All common materials

A variety of indications: one cement



Three radiopaque shades

A1 A2,5 TR

The three radiopaque shades (A1, A2.5, Transparent), together with high shade stability and optimum fluorescence, ensure a natural-looking aesthetic result. All shades are compatible with the Try-In pastes from DMG's Vitique System, allowing a reliable preview of the end result.

High abrasion resistance
For a tight marginal seal and long service life

One-piece excess removal
Comes off cleanly with the scaler

Aesthetics in 3 shades
With matching Try-In pastes

LuxaBond Universal

As universal as can be.

The new LuxaBond Universal makes work easy for you. The dual-curing, two-bottle bonding agent for dentine and enamel can also be used as a primer for indirect restorations, giving you more options than ever before. Because, unlike with many other bonding agents, "universal" here really does mean universal. In every respect.


Universal times four:
- Free choice of etching mode on dental hard tissue (self etch, total etch, selective etch)
- Can be used for direct and indirect restorations
- Suitable for light-curing, self-curing and dual-curing materials
- Dual-curing: light-curing and self-curing

One system, two bottles, no misapplication.
Thanks to the system's universal range of applications, mistakes caused by indications not being met are almost entirely ruled out.

"Do I have to etch first?", "Do I have to light-cure my adhesive?", "Can I also use it with self-curing materials?"
Questions you will no longer have to ask thanks to LuxaBond Universal. **The adhesive system gives you every single option.**


Maximum adhesive force, outstanding aesthetics, easy handling: PermaCem Universal is the ideal universal luting cement for every dental clinic.

Maximum adhesive force for universal use on all materials, in MPa



Internal measurements, DMG, Hamburg 2016. C (CAD/CAM composite), Z (zirconium oxide), S (silicate ceramic (IPS e.max)), N (non-precious alloy (Wirobond)), E (precious alloy (Degunorm))

Maximum adhesive force on enamel and dentine, in MPa



Internal measurements, DMG, Hamburg 2016. SED (self etch dentine), TED (total etch dentine), SEE (self etch enamel), TEE (total etch enamel)